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---

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## Introduction

Read the entire manual before using the device.

### Indications for use

[ ] is intended to provide ventilation for non-dependent, spontaneously breathing adult and pediatric patients [ ] with respiratory insufficiency, or respiratory failure, with or without obstructive sleep apnea. The device is for non-invasive use, or invasive use with an uncuffed or deflated tracheostomy.

### Contraindications

[ ] is contraindicated in patients who are unable to endure more than brief interruptions in ventilation. [ ] is not a life support ventilator.

The use of the device may be contraindicated in patients with:

- pneumothorax or pneumomediastinum
- pathologically low blood pressure, particularly if associated with intravascular volume depletion
- cerebrospinal fluid leak, recent cranial surgery or trauma
- severe bullous lung disease
- dehydration.

The use of the [ ] or pulse oximetry [ ] is contraindicated in an MRI environment.

### Adverse effects

Patients should report unusual chest pain, severe headache or increased breathlessness.

The following side effects may arise during the course of noninvasive ventilation with the device:

- drying of the nose, mouth or throat
- nosebleed
- bloating
- ear or sinus discomfort
- eye irritation
- skin rashes.

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##### [Principle of operation for Electrosurgical Generator]

-  Electrosurgical Generator -

##### [Principle of Operation]

The  () is indicated for use in micro, macro, and endoscopic bipolar cutting and coagulation of tissue and sealing of blood vessels in all areas of surgery. It is designed to cut and coagulate live tissue during surgical procedures. The system includes a generator, a bipedal footswitch, and a connecting cable. The generator is equipped with a voice synthesizer to provide audible indication of changes to power settings; and as an option, will announce the operating mode and power setting each time cutting or coagulation is performed. The connecting cable enables the use of a  with the generator.

Current flows between the active electrodes elements and the return electrode element, both elements or poles located on the itself in a bipolar configuration, providing a localized energy field. The result of this arrangement is controlled energy delivery with minimal collateral tissue damage.

The  works by passing RF energy through a conductive solution () in close proximity to or in contact with the tissue to be treated. The conductive solution forms a thin layer between the

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- 81.1 Theory of Operation
- 81.2 Monopolar Mode
- 81.3 Bipolar Mode
- 81.4 ESU Design
- 81.5 Active Electrodes
- 81.6 Dispersive Electrodes
- 81.7 ESU Hazards
- 81.8 Recent Developments

An electrosurgical unit (ESU) passes high-frequency electric currents through biologic tissues to achieve specific surgical effects such as cutting, coagulation, or desiccation. Although it is not completely understood how electrosurgery works, it has been used since the 1920s to cut tissue effectively while at the same time controlling the amount of bleeding. Cutting is achieved primarily with a continuous sinusoidal waveform, whereas coagulation is achieved primarily with a series of sinusoidal wave packets. The surgeon selects either one of these waveforms or a blend of them to suit the surgical needs. An electrosurgical unit can be operated in two modes, the monopolar mode and the bipolar mode. The most noticeable difference between these two modes is the method in which the electric current enters and leaves the tissue. In the monopolar mode, the current flows from a small active electrode into the surgical site, spreads through the body, and returns to a large dispersive electrode on the skin. The high current density in the vicinity of the active electrode achieves tissue cutting or coagulation, whereas the low current density under the dispersive electrode causes no tissue damage. In the bipolar mode, the current flows only through the tissue held between two forceps electrodes. The monopolar mode is used for both cutting and coagulation. The bipolar mode is used primarily for coagulation.

This chapter begins with the theory of operation for electrosurgical units, outlines various modes of operation, and gives basic design details for electronic circuits and electrodes. It then describes how improper application of electrosurgical units can lead to hazardous situations for both the operator and the patient and how such hazardous situations can be avoided or reduced through proper monitoring methods. Finally, the chapter gives an update on current and future developments and applications.

### 81.1 Theory of Operation

---

In principle, electrosurgery is based on the rapid heating of tissue. To better understand the thermodynamic events during electrosurgery, it helps to know the general effects of heat on biologic tissue. Consider a tissue volume that experiences a temperature increase from normal body temperature to 45°C within a few seconds. Although the cells in this tissue volume show neither microscopic nor macroscopic changes, some cytochemical changes do in fact occur. However, these changes are reversible, and the cells return to their normal function when the temperature returns to normal values. Above 45°C, irreversible changes

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TEST REPORT	
IEC 60601-1 Medical electrical equipment Part 1: General requirements for safety	
Report Reference No. ....	659100-2830-0008/58845
	CB-Ref. No.: DE 1 -
Tested by (name + signature) .....	Dirk Ritter 
Reviewed by (name + signature) .....	Michael Mayer 
Date of issue .....	
Contents .....	105 pages File: Testreport-FG12-8-58645_60601-1.doc
Testing Laboratory .....	VDE Testing and Certification Institute
Address .....	, Germany
Testing location/procedure .....	CB TL <input checked="" type="checkbox"/> RMT <input type="checkbox"/> SMT <input type="checkbox"/> WMT <input type="checkbox"/> TMP <input type="checkbox"/> Section FG12 / 110152
Applicant's name .....	
Address .....	
Manufacturer's name .....	
Address .....	
Factory's name .....	
Address .....	
Test specification:	
Standard .....	IEC 601 1:1988 + A1:1991 + A2:1995
Test procedure .....	CB Scheme
Non-standard test method .....	N/A
Test Report Form No. ....	601-1_C/
TRF Originator .....	Underwriters Laboratories Inc.
Master TRF .....	dated
<p><b>Copyright © 2005 IEC System for Conformity Testing and Certification of Electrical Equipment (IECEE), Geneva, Switzerland. All rights reserved.</b></p> <p><small>This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context.</small></p>	
Test item description .....	
Trade Mark .....	
Model/Type reference .....	
Ratings .....	100 – 240 V~ via external power supply unit (Power Brick) type

IEC60601-1 (ed.2)			
Clause	Requirement + Test	Result - Remark	Verdict
<b>6</b>	<b>IDENTIFICATION, MARKING AND DOCUMENTS</b>		—
6.1	Marking on the outside of equipment or equipment parts		P
6.1 c	Markings of the specific power supply affixed		P
6.1 d	If marking is not practicable due to size or nature of enclosure, information is included in accompanying documents		N/A
6.1 e	Name and/or trademark of the manufacturer or supplier .....		P
6.1 f	Model or type reference .....	l the frame ed on	P
6.1 g	Rated supply voltages or voltage range(s)	Ratings for 100 – 240 V~ 1.3 – 0.7 A 50/60 Hz ratings for Lithium Ion Battery, type: : 10,8V, 1Ah	P
	Number of phases .....	Single phase	P
	Type of current .....	Symbol ~	P
6.1 h	Rated frequency or rated frequency range(s) (Hz) ...	50/60 Hz	P
6.1 j	Rated power input (VA, W or A) .....	1.3 – 0.7 A	P
6.1 k	Power output of auxiliary mains socket-outlets		N/A
6.1 l	Class II symbol		P
	Symbol for degree of protection against ingress of water provided .....	IPX1 for hospital use (indoor) IP32 for "out of hospital"  IP31 all specified positions and IP32 special mounting position with supply leads from the bottom	P
	Symbol for protection against electric shock .....		P
	If equipment has more than one applied part with different degrees of protection, the relevant symbols are clearly marked on such applied parts, or on or near relevant outlets		P
	Symbol for protection of defibrillation-proof applied parts		P
	Symbol 14 from Table D1 for defibrillation-proof with protection partly in patient cable		P
6.1 m	Mode of operation (if no marking, suitable for continuous operation)	Continuous mode	P
6.1 n	Types and rating of external accessible fuses .....	Primary fuses of PowerBrick not accessible by operator	N/A
6.1 p	Ratings of external output .....		N/A
6.1 q	Symbol for physiological effect(s):		-
	- attention, consult accompanying documents		N/A

IEC60601-1 (ed.2)					
Clause	Requirement + Test			Result - Remark	Verdict
<b>56.1</b>	<b>TABLE: lists of critical component parts</b>				
Object/part No			Technical data	Standard	Mark(s) of conformity <sup>1)</sup>
1.1 SMPS			Input: 100-240VAC 50/60 Hz 1.3-0.7 A Output: 24VDC,30Whom. 40W/10min/hr 50W 1min/hr Ambient Temp: 0-65°C	UL, CSA, VDE IEC 60950 and requirements according IEC 60601-1	VDE
1.2 Fuse			Current: 2A characteristic: time lag T breaking capacity: high	EN/IEC 60127-2	UL, CSA
<b>Main - Board</b>					
2.1 LAN transformer			1.5 KV prim to sec		VDE
2.2 C197, C198, C199, C211, C210			Class X1/Y1 250V AC 2200pF	IEC60384-14	UL, CSA
2.8 RT1 overload protection			2,2KΩ to 100 KΩ (25°C)	-	-
2.9 outer enclosure			>PC+ABS+FR< UL flame rating: 94V-0	UL	QMFZ2, QMFZ8.
<b>ECG board</b>					
optocoupler			-5000 V rms/ 1 min UL recognized highest over voltage 6000 V Case Temperature Ts = 175 Input Current Is Input = 230 Output Power Ps Output = 600	UL 1577 IEC/EN/DIN EN 60747-5-2 Approval DIN VDE 0884	FPQU2, File E 55361
transformer			2,27 mH Test voltage 5 KVrms	IEC 85	-
Outer enclosure				IEC 85	QMFZ2.E41871
Potting material			94 V 0	IEC 85	QMFZ2.E96722 QMFZ8.E96722

Korea - Differences to IEC 60601-1:1988 + A1:1991 + A2:1995			
6.1j	Insert the following sub-clause between the second and third sub-clauses: Equipment for one or several RATED voltage or frequency ranges, the RATED input for 220 V, 60 Hz or if applicable for 110 V, 60 Hz shall be separately marked.	The input ratings for 220V, 60 Hz and 110V, 60 HZ will be explicitly marked on the device if the product is to be delivered to the Korean market.	Pass
6.1s	HIGH VOLTAGE TERMINAL DEVICES on the outside of EQUIPMENT which are accessible without the use of a TOOL shall be marked with the symbol "dangerous voltage" (see Appendix D, Table DII, Symbol 6) and wit the Korean language, .		N/A
6.2c	Replace the existing subclause wit the following: The presence of HIGH VOLTAGE PARTS shall be marked with the symbol "dangerous voltage" (see Appendix D, Table DII, Symbol 6) and in the Korean language.		N/A
6.8.1	Insert the following sub-clause after the last paragraph: Language of accompanying documents shall be included in Korean.	The documents will be translated into Korean.	Pass
General	National supply voltages are 110 V, 220 V and 380 V.		Pass
General	Only appliances having supply frequency of 60 Hz or a frequency range including 60 Hz are accepted.		Pass
General	Instruction manuals and appliance markings related to safety, including nameplate, shall be in Korean or graphical symbols in accordance with IEC Publication 417.	Confirmation the instruction manual will be translated into the Korean language provided via declaration letter.	Pass
General	Plugs for connection of the equipment to the supply mains shall comply with the Korean Standard (KSC 8305 and 8300). More details are available from KTL (c/o KTL) on request.	Instructions are provided to allow selection of a plug that complies with Korean Standard (KSC 8305 and 8300).	Pass

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Test Report issued under the responsibility of:

**NCB TÜV SÜD Product Service GmbH**  
Ridlerstr. 65  
D - 80339 München  
Germany

Product Service

<b>TEST REPORT</b> <b>IEC 60601-2-47</b> <b>Medical electrical equipment</b> <b>Part 2: Particular requirements for the safety, including essential performance, of ambulatory electrocardiographic systems</b>	
Report Reference No. ....	052-356765-000
Date of issue .....	Refer to CB certificate
Total number of pages .....	33
CB Testing Laboratory .....	TÜV SÜD PRODUCT SERVICE, Ridlerstrasse 65, D-80339 Munich
Address .....	Ridlerstrasse 65, D-80339 Munich
Applicant's name .....	
Address .....	
<b>Test specification:</b>	
Standard .....	IEC 60601-2-47: 2001 (First Edition) for use in conjunction with IEC 60601-1:1988 + A1:1991 + A2:1995
Test procedure .....	CB Scheme
Non-standard test method .....	N/A
Test Report Form No. ....	IEC60601_2_47C
Test Report Form(s) Originator .....	Underwriters Laboratories Inc.
Master TRF .....	Dated 2006-11
<b>Copyright © 2006 IEC System for Conformity Testing and Certification of Electrical Equipment (IECEE), Geneva, Switzerland. All rights reserved.</b>	
<small>This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context.</small>	
<small>If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the reference to the CB Scheme procedure shall be removed.</small>	
<b>This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.</b>	
Test item description .....	
Trade Mark .....	
Manufacturer .....	
Model/Type reference .....	
Ratings .....	

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TEST REPORT		
Page: 01 DOC.NUM.: R60601100 Rev: 4.0 Date: 1/2000	<b>IEC 60601-1,-1-1, -1-3, -1-4,-2-7, -2-28, 2-32</b> <b>Medical electric equipment</b> <b>Part 1: General requirements for safety</b>	
<b>Report</b>		
Reference No.....	12307250 001	
Compiled by (+ signature) .....	<input type="text"/>	
Approved by (+ signature) .....	<input type="text"/>	
Date of issue .....	2008-07-29	
Contents.....	122 pages	
<small>This report is based on a blank test report that was prepared by KEMA using information obtained from the TRF originator (see below).</small>		
<b>Testing laboratory</b>		
Name .....	TÜV Rheinland Japan Ltd., Global Technology Assessment Center	
Address.....	4-25-2 Kita-Yamata, Tsuzuki-ku, Yokohama 224-0021 Japan	
Testing location.....	<input type="text"/>	
<b>Client</b>		
Name .....	<input type="text"/>	
Address.....	<input type="text"/>	
<b>Test specification</b>		
Standard .....	IEC 60601-1:1988+A1:1999+A2:1995 IEC 60601-1-1:2000 IEC 60601-1-3:1994 IEC 60601-1-4:1996 + A1:1999 IEC 60601-2-7:1998 IEC 60601-2-28:1993 IEC 60601-2-32:1994	
Test procedure .....	CB scheme	
Procedure deviation.....	N.A.	
Non-standard test method.....	N.A.	
<b>Test Report Form/blank test report</b>		
Test Report Form No.....	1001-1-01/IEC60601-1-1-D/IEC60601-1-3B/IEC60601-1-4C /IEC60601-2-28B/IEC60601-2-32C/IEC60601-2-7C	
TRF originator.....	UL / UL / UL / UL / UL / UL / UL	
Master TRF.....	Dated 97-04 / 06-09 / 06-09 / 06-09 / 06-09 / 06-10 / 06-11	
<small>Copyright reserved to the bodies participating in the Committee of Certification Bodies (CCB) and/or the bodies participating in the CENELEC Certification Agreement (CCA).</small>		
<b>Test item</b>		
Description.....	X-RAY TV SYSTEM	
Trademark .....	<input type="text"/>	
Model and/or type reference .....	X-RAY TV SYSTEM SONIALVISION safire17	
Manufacturer.....	Same as applicant	

IEC 60601-1


Factory .....	:	<input type="text"/>
Input Rating(s) .....	:	3N 200/380/400/415/440/480VAC, 133kVA, 50/60Hz 1N 200/220/230/240VAC, 9.8kVA, 50/60Hz
<b>Particulars: test item vs. test requirements</b>		
Classification of installation and use .....	:	Permanently installed
Supply connection.....	:	Permanently installed
<b>Test case verdicts</b>		
Test case does not apply to the test object.....	:	N/A
Test item does meet the requirement .....	:	P(ass)
Test item does not meet the requirement .....	:	F(ail)
<b>Testing</b>		
Date of receipt of test item .....	:	N/A
Date(s) of performance of test.....	:	2008-06-27 – 2008-06-30
<b>General remarks</b>		
<p>This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IEC 60601-1-1:2001.</p> <p>"This test report shall not be reproduced except in full without the written approval of the testing laboratory. The test results presented in this report relate only to the item tested.</p> <p>"(see remark #)" refers to a remark appended to the report.</p> <p>"(see appended table)" refers to a table appended to the report (same number of table as number of clause).</p> <p>Throughout this report a comma is used as the decimal separator.</p> <p><b>Summary of compliance with National Differences:</b>  EN 60601-1:1990+A1:1993+A2:1995+A13:1996  EN 60601-1-1:2001  EN 60601-1-3:1994  EN 60601-1-4: 1996 + A1:1999  EN 60601-2-7:1998  EN 60601-2-28:1993  EN 60601-2-32:1994  UL 60601-1:2003 R6.03  CAN/CSA-C22.2 No. 601.1-M90  JIS T 0601-1:1999  JIS T 0601-1-1:2005</p> <p>Summary of compliance with National Differences (for explanation of codes see below): US, CA, JP  CA=Canada, US=United States of America, JP=Japan</p>		

IEC 60601-1-3			
Clause	Requirement Test	Result – Remark	Verdict
4	<b>General Requirements For Tests</b>		P
4.201.1	Statement of compliance is expressed in accordance with this Standard .....	Statement of compliance is stated in Operation Manual. Refer to clause 6.8 of the ICE60601-1 part of this test report.	P
6	<b>IDENTIFICATION, MARKING, AND DOCUMENTS</b>		P
6.1.202	All sub-assemblies, components, and ACCESSORIES of X-RAY EQUIPMENT that can be removed are marked so that they can be identified and correlated with their ACCOMPANYING DOCUMENTS	All sub-assembly components are adequately marked for correlation with their respective accompanying documents	P
	- Interchangeable devices are distinguishable to the OPERATOR	No interchangeable devices	N/A
	Markings permanently affixed and clearly legible per Clause 6 of Part 1	See general standard IEC60601-1 part of this test report	P
6.8	<b>Accompanying documents</b>		P
6.8.201	ACCOMPANYING DOCUMENTS include the required statements per sub-clauses in Table 202	Operation Manual provided	P
6.8.202	ACCOMPANYING DOCUMENTS identify the items they refer to and include replication of information required to be marked on the items		P
	Information on the location and content of required markings on subassemblies is provided	Information provided	P
	a) For accessible markings, the location of the markings and instructions enabling them to be inspected		P
	b) For markings inaccessible in complete assembly		P
	1) Replication of the information marked on the sub-assemblies, or	Replicated in manual	P
	2) A list of the subassemblies with references to their accompanying documents		P
	ASSEMBLING INSTRUCTIONS for items supplied separately from main assemblies are provided	To be installed by the authorized service personnel. Installation manual provided	P
29	<b>X-RADIATION</b>		P
29.201.1	The operating range of X-RAY TUBE VOLTAGE for dental applications restricted per Table 203 (kV).....	Not intended for dental applications	N/A
29.201.2	First HALF-VALUE LAYERS attained were not less than the minimum values in Table 204 (mm Al) .....	2.78mm Al at 70kV (>2.1mm required)	P

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	<h1>EMC</h1>	Ref. Certif. No. <b>FI 6784</b>
<b>IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE)      CB EMC TEST CERTIFICATE</b>		
Product  Name and address of the applicant  Name and address of the manufacturer  Name and address of the factory  Ratings and principal characteristics  Trademark (if any)  Model / Type Ref.  Additional information (if necessary may also be reported on page 2)  A sample of the product was tested and found to be in conformity with  As shown in the Test Report Ref. No. which forms part of this Certificate	<div style="border: 1px solid black; padding: 5px;"> <p><b>Electronic APEX Locator</b></p> <div style="border: 1px solid black; height: 40px; margin-bottom: 5px;"></div> <p>Same as applicant</p> <p>Same as applicant</p> <p>Battery: 4.5Vd.c., (1.5Vd.c. x 3EA, type: LR3, size: AAA)</p> <div style="border: 1px solid black; height: 40px; margin-bottom: 5px;"></div> <p><b>ROOTOR, i-ROOT 100 mini</b></p> <p><b>IEC 60601-1-2:2007 (Ed 3)</b></p> <p><b>F690501/RF-EMG003932</b></p> </div> <p>This CB Test Certificate is issued by the National Certification Body of Finland</p>	
SGS Fimko's ref.no. 264523-1  Date: 21 June 2011	<div style="text-align: center;"> <b>SGS Fimko Ltd</b>            Ari Honkala          Product Line Manager       </div> <div style="text-align: right; font-size: 2em; font-weight: bold; margin-top: 10px;">  </div>	

Test Report issued under the responsibility of:

<b>TEST REPORT</b> <b>IEC 60601-1-2</b> <b>Medical Electrical Equipment</b> <b>PART 1-2: General Requirements for Basic Safety and Essential Performance</b> <b>Collateral Standard: Electromagnetic Compatibility</b>	
Report Reference No. ....	F690501/RF-EMG003932
Date of issue .....	June 16, 2011
Total number of pages .....	
CB Testing Laboratory .....	SGS Korea Co.,Ltd.
Address .....	#18-34, Sanbon-dong, Gunpo, Gyeonggi-do, Korea
Applicant's name .....	
Address .....	
<b>Test specification:</b>	
Standard .....	IEC 60601-1-2 (ed 3)
Test procedure .....	CB
Non-standard test method .....	N/A
Test Report Form No .....	IEC60601_1_2BEMC
Test Report Form(s) Originator .....	Underwriters Laboratories Inc.
Master TRF .....	Dated 2008-06
<b>Copyright © 2008 IEC System for Conformity Testing and Certification of Electrical Equipment (IECEE), Geneva, Switzerland. All rights reserved.</b>	
<small>This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context.</small>	
<small>If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the reference to the CB Scheme procedure shall be removed.</small>	
<b>This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.</b>	
Test item description .....	Electronic APEX Locator
Trade Mark .....	
Manufacturer .....	
Model/Type reference .....	ROOTOR (added Mode : i-ROOT 100 mini)
Ratings .....	Battery : 4,5 Vd.c.(1,5 Vd.c. x 3 EA, type:LR3, size:AAA)

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
**TABLE OF CONTENTS**

Title Page  
Table of Contents  
Study Summary  
Quality Assurance Statement  
Study Director Signature and Verification Dates

1.0 Purpose  
2.0 References  
3.0 Compliance  
4.0 Identification of Test and Control Articles  
5.0 Identification of Test System  
6.0 Justification of Test System and Route of Administration  
7.0 Experimental Design and Dosage  
8.0 Evaluation Criteria  
9.0 Results  
10.0 Conclusion  
11.0 Records  
12.0 Confidentiality Agreement

**STUDY SUMMARY**

No biological reactivity (Grade 0) was observed in the L929 mammalian cells at 48 hours, post exposure to the test article extract. The observed cellular response obtained from the positive control article extract (Grade 4) and negative control article extract (Grade 0) confirmed the suitability of the test system. Based on the criteria of the protocol, the test article,

, is considered non-cytotoxic and meets the requirements of the Elution Test, ISO 10993-5 guidelines.

---


**QUALITY ASSURANCE STATEMENT**

This study was conducted in compliance with U.S. Food and Drug Administration regulations set forth in 21 CFR, Part 58.

The sections of the regulations not performed by or under the direction of Toxikon Corporation, exempt from this Good Laboratory Practice Statement, included characterization and stability of the test article and its mixture with carriers, 21 CFR, Parts 58.105 and 58.113.

The Quality Assurance Unit conducted inspections on the following dates. The findings were reported to the Study Director and to Toxikon's Management.

INSPECTIONS	DATE OF INSPECTION	DATE REPORTED STUDY DIRECTOR	DATE REPORTED MANAGEMENT
DOSE ADMINISTRATION	08/25/06	08/25/06	08/25/06
RAW DATA	09/01/06	09/01/06	09/01/06
FINAL REPORT	09/01/06	09/01/06	09/01/06

  
Allison Lyons-Hook, B.A.  
Quality Assurance

09/01/06  
Date

---

**STUDY DIRECTOR SIGNATURE AND VERIFICATION DATES**

This study meets the technical requirements of the protocol. The study also meets with the requirements of the Good Laboratory Practice Regulations, 21 CFR, Part 58, with the exemptions as stated in the Quality Assurance Statement.

Protocol Number: AMD/VITRO/003-06/000

Study Director: Patricia S. Grutkoski, Ph.D.

Company: Toxikon Corporation

Signature:

Pat S. Grutkoski

Date:

9/1/06

Study Supervisor:

Patricia S. Grutkoski, Ph.D.

**VERIFICATION DATES:**

Protocol Effective Date: 08/11/06

Test Article Receipt: 08/22/06

Project Log Date: 08/22/06

Extraction Dates: 08/24/06 – 08/25/06

Technical Initiation: 08/25/06

Technical Completion: 08/27/06

**1.0 PURPOSE**

The purpose of the study was to determine the biological reactivity of a mammalian cell culture (L929) in response to the test article. The study was designed for the evaluation of test article extracts (e.g., polymeric materials, high-density materials, etc.).

**2.0 REFERENCES**

The study was conducted based upon the following references:

- 2.1 ISO 10993-5, 1999, Biological Evaluation of Medical Devices – Part 5: Tests for *In Vitro* Cytotoxicity.
- 2.2 ISO 10993-12, 2002, Biological Evaluation of Medical Devices – Part 12: Sample Preparation and Reference Materials.
- 2.3 ISO/IEC 17025, 2005, General Requirements for the Competence of Testing and Calibration Laboratories.

**3.0 COMPLIANCE**

The study conformed to the current FDA 21 CFR, Part 58 – Good Laboratory Practice for Non-Clinical Laboratory Studies.

**4.0 IDENTIFICATION OF TEST AND CONTROL ARTICLES**

The Sponsor supplied the following information on a test requisition form or other correspondence, wherever applicable (excluding confidential or trade secret information). The Sponsor was responsible for all test article characterization data as specified in the GLP regulations.

**4.1 Test Article:**

Test Article Name:

CAS/Code #: Not Supplied by Sponsor (N/S)

Lot/Batch #:

Physical State: N/S

Color: N/S

Expiration Date: N/S

Density: N/S

Stability: N/S

Solubility: N/S

pH: N/S

Storage Conditions: Room Temperature

Safety Precautions: Standard Toxikon Laboratory Safety Precautions

#### 4.2 Control Articles (Toxikon Supplied):

##### 4.2.1 Negative Control Article Name: Negative Control High Density Polyethylene (Negative Control Plastic)

Toxikon QC #: CSC-04-05-009-CC

Physical State: Solid

Color: White

Storage Conditions: Room Temperature

Safety Precautions: Standard Laboratory Safety Precaution

##### 4.2.2 Positive Control Article Name: Natural Rubber

Toxikon QC #: CSC-05-03-005-CC

Physical State: Solid

Color: Amber

Stability: Stable at Room Temperature

Storage Conditions: Room Temperature

Safety Precautions: Standard Laboratory Safety Precautions

### 5.0 IDENTIFICATION OF TEST SYSTEM

The test system was mouse fibroblast L929 cells. The cell line was obtained from the American Type Culture Collection, Manassas, Virginia.

### 6.0 JUSTIFICATION OF TEST SYSTEM AND ROUTE OF ADMINISTRATION

6.1 Historically, mouse fibroblast L929 cells had been used for cytotoxicity studies because they demonstrate sensitivity to extractable cytotoxic articles.

6.2 The test article was extracted and administered *in vitro* to mouse fibroblast L929 cells through a solvent compatible with the test system. This was the optimal route of administration available in this test system as recommended in the ISO 10993-5 guidelines.

### 7.0 EXPERIMENTAL DESIGN AND DOSAGE

#### 7.1 Preparation of Test and Control Articles:

7.1.1 The extraction vehicle was Serum-Supplemented (Complete) Minimum Essential Medium (MEM) supplemented with 10% fetal bovine serum. Per Sponsor's request, representative samples of each component of the test article were pooled and extracted. The test article (2.0 g) was extracted in 10 mL of medium (0.2 g per 1.0 mL ratio) at  $37 \pm 1$  °C for  $24 \pm 2$  hours.

7.1.2 Extracts prepared with complete MEM were tested at 100% (neat) concentration.

7.1.3 The positive control article (Natural Rubber) and negative control article (Negative Control Plastic) were similarly prepared and extracted.

**7.2 Pre-Dose Procedure:****7.2.1 Cell Culture Preparation:**

Cultures of L929 cells were initiated in complete MEM not less than 24 hours prior to testing and incubated at  $37 \pm 1$  °C, in a humidified atmosphere containing  $5 \pm 1\%$  carbon dioxide.

7.2.2 Extracts were evaluated for a change in pH. No adjustment was necessary since the pH of the extract was within  $\pm 0.2$  of the control medium. The extracts were filter sterilized by passage through a  $0.22 \mu\text{m}$  pore filter, prior to being applied to the cell monolayer.

**7.3 Dose Administration:**

A volume of 3 mL of the test article or control article extract was used to replace the maintenance medium of the cell culture. All extracts were tested in triplicate.

**7.4 Post-Dose Procedure:**

All cultures were incubated for  $48 \pm 2$  hours, at  $37 \pm 1$  °C, in a humidified atmosphere containing  $5 \pm 1\%$  carbon dioxide.

**8.0 EVALUATION CRITERIA**

8.1 The response of the cell monolayer is evaluated under a microscope. Cytochemical stains are used in the evaluation. The biological reactivity (cellular degeneration and malformation) is rated on a scale of 0 to 4. The test system is considered suitable if no signs of cellular reactivity (Grade 0) are noted for the negative control article and the positive control article showed greater than a Mild reactivity (Grade 2). The test article meets the requirements of the test if none of the cultures treated with the test article show greater than a Mild reactivity (Grade 2). If the suitability of the system is not confirmed, the test is repeated.

<b>Grade</b>	<b>Reactivity</b>	<b>Description of Reactivity Zone</b>
0	None	Discrete intracytoplasmic granules; no cell lysis
1	Slight	Not more than 20% of the cells are round, loosely attached, and without intracytoplasmic granules; occasional lysed cells are present
2	Mild	Not more than 50% of the cells are round and devoid of intracytoplasmic granules; no extensive cell lysis and empty areas between cells
3	Moderate	Not more than 70% of the cell layers contain rounded cells or are lysed
4	Severe	Nearly complete destruction of the cell layers

8.2 The study and its design employ methodology to minimize uncertainty of measurement and control of bias for data collection and analysis.

## 9.0 RESULTS

At the 48 hour observations, no signs of cellular reactivity (Grade 0) were observed for the negative control article and the test article extracts. Severe signs of reactivity (Grade 4) were noted for the positive control article extract at the 48 hour observation.

REACTIVITY GRADES

Time	Date	Test Article			Controls								
					Medium			Negative			Positive		
		A	B	C	A	B	C	A	B	C	A	B	C
0 Hours	08/25/06	0	0	0	0	0	0	0	0	0	0	0	0
24 Hours	08/26/06	0	0	0	0	0	0	0	0	0	3	3	3
48 Hours	08/27/06	0	0	0	0	0	0	0	0	0	4	4	4

## 10.0 CONCLUSION

No biological reactivity (Grade 0) was observed in the L929 mammalian cells at 48 hours, post exposure to the test article extract. The observed cellular response obtained from the positive control article extract (Grade 4) and negative control article extract (Grade 0) confirmed the suitability of the test system. Based on the criteria of the protocol, the test article, \_\_\_\_\_ is considered non-cytotoxic and meets the requirements of the Elution Test, ISO 10993-5 guidelines.

## 11.0 RECORDS

- 11.1 Original raw data are archived at Toxikon Corporation.
- 11.2 A copy of the final report and any report amendments is archived at Toxikon Corporation.
- 11.3 The original final report, and a copy of any protocol amendments or deviations, is forwarded to the Sponsor.
- 11.4 All unused test article shall be disposed of by Toxikon, per Sponsor's request.

## 12.0 CONFIDENTIALITY AGREEMENT

Statements of confidentiality were as agreed upon prior to study.

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Page 1 of 4

English version of the  
“**Certificato di Conformità alla Buona Pratica di Laboratorio**”  
issued by the Italian Ministry of Health (“Ministero della Salute”)  
*This translation has been prepared by Eurofins Biolab S.r.l.  
the original Italian document bear witness to the authenticity of this document.*

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**STATEMENT OF COMPLIANCE WITH THE GOOD LABORATORY  
PRACTICE**

(Legislative Decree 02<sup>nd</sup> March 2007, N. 50 – Directive 2004/9/EC)  
Rome, October 12th 2010 Certification: 121/2010

In conformity with the article 4 paragraph 1 of the Legislative Decree N. 50 of 2<sup>nd</sup> March 2007, with the article 5 of ministerial decree of 4<sup>th</sup> July 1997 and considering positive results of inspection performed on June the 21<sup>st</sup> to 22<sup>nd</sup> 2010

**WE CERTIFY**

that the Test Facility

Is qualified to carry out in compliance with the Good Laboratory Practice as set forth in the Annex II of the aforesaid Legislative Decree:

C.SOC. € 100.000 i.v.  
P. IVA 00762140960  
C.F. 03765750157  
REA MI 966696  
D-U-N-S 429117112  
CIT005 4-385

### 1) Physical-Chemical Tests

- Percent purity
- Impurities characterization, including isomers and secondary products
- Principal impurities percent
- Identification and dosing of stabilizing compounds, including inhibitors or other additives,
- Spectra characterisation
- Elemental analysis
- Melting point
- Boiling point and boiling range
- Relative density
- Vapour pressure
- Water Solubility
- Liposolubility
- Partition coefficient
- Hydrolysis rate as pH function
- Organic solvent solubility
- Dissociation constant in water
- Thermal stability and degradation products identification
- Flash point – solids (EEC methods)
- Flammability
- Oxidising properties – liquids (Federal Register Method)
- Oxidising properties – solids (EEC method)
- Appearance
- Acidity/Alcalinity and pH
- Wettability
- Foaming persistence
- Suspension properties and stability
- Sieving test and indication on granules ponderal distribution
- Active ingredient content in or on particles of baits, granules or treated seeds
- Emulsion, re-emulsion and emulsion stability properties
- Physical and chemical compatibility with other properties
- Analytical methods to determine phytosanitary product composition
- Apparent density with or without compaction
- Loss on drying
- Sulphuric Ashes
- Water content determination
- Granulometry
- Dissolution rate of hydrosoluble bags

### 2) Toxicity Studies

- Acute and Repeated dose (28 days) Toxicity Studies
- *In vitro* and *in vivo* skin and eye irritation and corrosion
- Skin sensitization Studies

C.SOC. € 100.000 i.v.  
P. IVA 00762140960  
C.F. 03765750157  
REA MI 966696  
D-U-N-S 429117112  
CIT005 4-385

### 3) Mutagenesis Studies

- Salmonella Typhimurium reversion assay

### 6) Residual Studies

Analysis methods, included recovery rates and residual limits of detection in:

- Treated plant, vegetal product, food and feed
- Soil
- Water (including potable water)
- Biological fluids as well as animal and human tissue

### 9) Others

#### Biocompatibility Studies

Studies are performed both in vivo and in vitro:

- Citotoxicity studies
- Cutaneous irritation studies
- Intracutaneous reactivity
- Mucosae irritation
- Skin sensitization studies
- Systemic Toxicity studies
- Implant Studies

#### Stability Studies

- These kind of studies are performed at different environmental conditions

#### Efficacy Evaluation Studies

- Bactericidal activity
- Fungicidal activity
- Sporicidal activity
- Evaluation of the efficacy of the preserving system
- Palatability and rodenticide efficacy studies
- Virucidal activity
- Mycobactericidal activity

#### Set up and validation of analytical methods studies

- Analytical methods to determine the pure active substance, its degradation products, isomers, impurities and additives
- Microbiological methods
- Biological methods for in vitro and in vivo dosage
- Methods for endotoxins determination (LAL Assay)

#### Ecotoxicological Studies

- Ready biodegradability (OECD 301)
- Aerobic biodegradability (ISO 14593)
- Inherent biodegradability (OECD 302)
- Final aerobic biodegradability (ISO11734)

C.SOC. € 100.000 i.v.  
P. IVA 00762140960  
C.F. 03765750157  
REA MI 966696  
D-U-N-S 429117112  
CIT005 4-385

- Active Sludge Respirometry (OECD209)
- Inhibition assay on Algae growth (OECD 201)
- Immobilization on Daphnie (OECD 202)
- Fish acute toxicity assay, static conditions (OECD 203)
- Fish prolonged toxicity assay, semistatic conditions (OECD204)
- Daphnie reproduction assay (OECD 211)

**The present certificate is considered valid during 24 months from the issue date.**


This certification will not be renewed, after its expiry, if no request is presented for renewal, according to paragraph 2 of article 4 of D.M. 4th July 1997 (issued on G.U. N. 229 of 1/10/1997), **and at the latest 90 days before the expiry.**

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Assessment of conformity with GLP according to Chemikaliengesetz and Directive 2004/9/EC at:

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The above mentioned test facility/test site is included in the national GLP Compliance Programme and is inspected on a regular basis

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Based on the inspection report it can be confirmed that this test facility/test site is able to conduct the aforementioned studies in compliance with the Principles of GLP

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- ISO 10993-5:2009
- USP, Biological Reactivity Tests, *In Vitro*

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3.4 음성대조물질 : High-density Polyethylene RM

3.5 양성대조물질 : Natural Rubber

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: Mouse Fibroblast Cells (ATCC CCL 1, Clone 929, of Strain L, or equivalent source)

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Grade	Reactivity	Conditions of all Cultures
0	None	Noncytotoxic: Discrete intracytoplasmic granules; no cell lysis
1	Slight	Slightly cytotoxic: Not more than 20% of the cells are round, loosely attached, and without intracytoplasmic granules; occasional lysed cells are present
2	Mild	Mildly cytotoxic: Not more than 50% of the cells are round and devoid of intracytoplasmic granules;no extensive cell lysis and empty areas between cells
3	Moderate	Moderately cytotoxic: Not more than 70% of the cell layers contain rounded cells or are lysed
4	Severe	Severely cytotoxic: Nearly complete destruction of the cell layers

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	Cytotoxicity Score								
	Initial			After 24hrs			After 48hrs		
초음파펄스진단장치 <input type="text"/> ( <input type="text"/> 의 구성품))	0	0	0	0	0	0	0	0	0
Negative control	0	0	0	0	0	0	0	0	0
Positive control	0	0	0	4	4	4	4	4	4

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	Document form: Software Development Version: <input style="width: 50px;" type="text"/>	
<h2 style="margin: 0;">Software Requirements Specification</h2> <div style="text-align: center; margin-top: 10px;"> <input style="width: 200px; height: 20px;" type="text"/>  <input style="width: 100px; height: 20px;" type="text"/> </div>		
Project business name: <input style="width: 100px;" type="text"/>	Internal project name: <input style="width: 50px;" type="text"/>	
<b>Made by:</b> Position: <input style="width: 150px;" type="text"/> Name: <input style="width: 300px;" type="text"/> Date: <input style="width: 80px;" type="text"/> Signature: <input style="width: 300px;" type="text"/>		
<b>Approved by:</b> Position: <input style="width: 80px;" type="text"/> <input style="width: 150px;" type="text"/> Name: <input style="width: 150px;" type="text"/> Date: 3 <input style="width: 80px;" type="text"/> Signature <input style="width: 150px;" type="text"/>	<b>Approved by:</b> Position: <input style="width: 100px;" type="text"/> Name: <input style="width: 150px;" type="text"/> Date: <input style="width: 80px;" type="text"/> Signature <input style="width: 150px;" type="text"/>	<b>Approved by:</b> Position: <input style="width: 100px;" type="text"/> Name: <input style="width: 150px;" type="text"/> Date: 3 <input style="width: 80px;" type="text"/> Signature: <input style="width: 150px;" type="text"/>
Valid from: <input style="width: 100px;" type="text"/>		
Document revision:		

RECORD OF CHANGES						
Revision	Date	Last version			Change Description	Approved by
		Version	Valid from	Valid to		

**1. Table of contents:**

- 1. General context..... 2
  - 1.1. User Interface ..... 3
  - 1.2. Remote Control/..... 5
  - 1.3. Recording..... 6
    - 1.3.1. converter data reading – ..... 6
    - 1.3.2. converter data reading – ..... 7
    - 1.3.3. Runtime checks and warnings..... 8
    - 1.3.4. File system and ..... 8

## 2. General context

The software/[ ] shall be implemented as an [ ] into a specific hardware environment. The environment is specified in Tab 1.1. The software shall have [ ] as described on state diagram in Fig. 1.1.

1. User interface - setting up and starting a recording, [ ]
2. Recording – recording of the ecg signal to the [ ]
3. Remote control/online transmission – communication with [ ] connection


Fig.1.1: State diagram – basic software states and [ ]

Tab 1.1: Hardware specification

Component	Specification
CPU	[ ]
File system	[ ]
RAM	[ ]
Display	[ ]
Storage	[ ]
Peripheral Port	[ ]
Wireless Connection	[ ]

### 1.1. User Interface

All the user interface shall be implemented as a set of screens and [ ]. Because of the limited operation memory size only one common [ ] and one modal [ ] can be allocated at a time. [ ] shall be connected to a message loop from buttons and shall be routed from application [ ] of the particular [ ]. [ ] is responsible for creating a new one before [ ].

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Document form: Software Development  
Version:

## Software Requirements Specification

**Project business name:**  **Internal project name:**

**Made by:**  
**Position:**   
**Name:**   
**Date:**   
**Signature:**

<p><b>Approved by:</b>  <b>Position:</b> <input type="text"/>  <b>Name:</b> <input type="text"/>  <b>Date:</b> <input type="text"/>  <b>Signature:</b> <input type="text"/></p>	<p><b>Approved by:</b>  <b>Position:</b> <input type="text"/>  <b>Name:</b> <input type="text"/>  <b>Date:</b> <input type="text"/>  <b>Signature:</b> <input type="text"/></p>	<p><b>Approved by:</b>  <b>Position:</b> <input type="text"/>  <b>Name:</b> <input type="text"/>  <b>Date:</b> <input type="text"/>  <b>Signature:</b> <input type="text"/></p>
---	---	---

**Valid from:**

**Document revision:**

RECORD OF CHANGES						
Revision	Date	Last version			Change Description	Approved by
		Version	Valid from	Valid to		

**CONTENT**

- 1 FUNCTIONAL REQUIREMENTS ..... 3
  - 1.1 Performance ..... 3
    - 1.1.1 Purpose of the software ..... 3
    - 1.1.2 [ ] requirements ..... 3
  - 1.2 Physical characteristics ..... 3
  - 1.3 Computing environment under which the software is to perform ..... 3
- 2 SOFTWARE SYSTEM INPUTS AND OUTPUTS ..... 4
  - 2.1 Inputs: ..... 4
  - 2.2 Outputs: ..... 4
- 3 INTERFACES ..... 5
  - 3.1 [ ] format ..... 5
    - 3.1.1 [ ] format requirements ..... 5
    - 3.1.2 Main format characteristics ..... 5
    - 3.1.3 [ ] Types ..... 5
  - 3.2 Online signal transmission ..... 7
    - 3.2.1 Online transmission protocol requirements ..... 8
    - 3.2.2 Main format characteristics ..... 8
    - 3.2.3 Recording settings via [ ] ..... 9
- 4 Graphical user interface ..... 9
  - 4.1 Starting a new recording ..... 9
  - 4.2 Continue a [ ] recording ..... 9
  - 4.3 Recording settings ..... 9
  - 4.4 Unit configuration ..... 9
  - 4.5 Information ..... 10
- 5 SOFTWARE-DRIVEN ALERTS ..... 10
  - 5.1 [ ] ..... 10
  - 5.2 [ ] ..... 10

# 1 FUNCTIONAL REQUIREMENTS

## 1.1 Performance

### 1.1.1 Purpose of the software

The [Redacted] software is intended to be used by trained operators under direct supervision of a licensed health care practitioner on adult and pediatric patients. The software shall consist of [Redacted]. The first part drives the hardware components of the unit to acquire and store the [Redacted] measured on its input and the second part provides the operator with a user interface.

The user interface part of the software shall allow the user to setup the recording, start the recording, setup the unit, display information about the unit. Also the software shall signal when [Redacted] get loose or fallen off a patient and allows short voice recordings on user's demand.

Data acquisition functions of the software are sampling and reading [Redacted], [Redacted] detection, detection of [Redacted] and saving all data to a [Redacted].

The software does not perform any medical related data analysis.

### 1.1.2 [Redacted] requirements

Software [Redacted] shall use unit's internal [Redacted] for all actions especially for signal sampling. Standard sampling frequency shall be [Redacted].

## 1.2 Physical characteristics

The software shall be implemented as an embedded system in [Redacted] language. All the functionality shall be implemented specifically for the hardware provided which is specified in the next chapter.

### 1.3 Computing environment under which the software is to perform

Component	Specification
CPU	[Redacted]
File system	[Redacted]
RAM	[Redacted]
Display	[Redacted]
Storage	[Redacted]
Peripheral Port	[Redacted]
Wireless Connection	[Redacted]

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<b>SOFTWARE VALIDATION REPORT</b>			
Project business name: <input style="width: 80%;" type="text"/>		Internal project name: <input style="width: 80%;" type="text"/>	
		Document revision: <input style="width: 40%;" type="text"/>	Date: <input style="width: 40%;" type="text"/>
Document Title: <b>Software Validation Report</b>			
Approval			
Title	Name	Signature	Date
<b>Project manager</b>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>
<b>Head of Development Department</b>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>
<input style="width: 80%;" type="text"/> manager	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>

**Record of Changes**

Rev.	Date	Author	Comments

Software Validation Report

**CONTENT**

1	Validation of <input style="width: 150px;" type="text"/>	3
2	FUNCTIONAL REQUIREMENTS	3
2.1	Performance PASSED!	3
2.1.1	Purpose of the software	3
2.1.2	<input style="width: 50px;" type="text"/> requirements PASSED!	3
2.1.3	Physical characteristics PASSED!	3
2.1.4	Computing environment under which the software is to perform PASSED!	4
2.2	Software system inputs and outputs PASSED!	4
2.2.1	Inputs: PASSED!	4
2.2.2	Outputs: PASSED!	4
2.3	Interfaces PASSED!	4
2.3.1	<input style="width: 100px;" type="text"/> PASSED!	4
2.3.2	<input style="width: 100px;" type="text"/> transmission PASSED!	6
2.3.3	Graphical user interface PASSED!	8
2.3.4	Software-driven alerts PASSED!	9
3	Test report	9
4	Risk analysis	9
5	Verification of <input style="width: 100px;" type="text"/>	10
6	Result of validation	10

## 1 Validation of SW [ ]

Software validation refers to establishing, by objective evidence, that the software conforms to the user needs and intended uses of the device. Software validation is a part of design validation of the finished device. It involves checking for proper operation of the software in its actual or simulated use environment, including integration into the final device where appropriate. Software validation is highly dependent upon comprehensive software testing and other verification tasks previously completed at each stage of the software development life cycle. Planning, verification, configuration management and many other aspects of good software engineering are important activities that together help to support the conclusion that the software is validated.

Validation is done according to Std 1012-1998 IEEE Standard for Software *Verification and Validation*, *IEEE Standards Collection: Software Engineering*. [ ]

Document used for validation of [ ]:

- Software requirements specification [ ] - definition of functionality of system [ ]
- Test report version [ ] according Test protocol [ ]
- Risk analysis [ ] during [ ] development
- Verifications of SW [ ]

Validation is done according to internal directive [ ] version [ ]

Risk analysis is done according to [ ] - Risk Analysis version [ ]

The following paragraphs refer to Software requirements specification. All points described there were checked during validation.

## 2 FUNCTIONAL REQUIREMENTS

### 2.1 Performance

**PASSED!**

#### 2.1.1 Purpose of the software

The [ ] unit software is intended to be used by [ ] operators under direct supervision of a licensed health care practitioner on adult and pediatric patients. The software shall consist of two parts. The first part drives the hardware components of the unit to acquire and store the [ ] measured on its input and the second part provides the operator with a user interface.

The user interface part of the software shall allow the user to setup the recording, start the recording, setup the unit, display information about the unit. Also the software shall signal when [ ]

Data acquisition functions of the software are sampling and reading [ ] detection, detection of artificial cardiac [ ]

The software does not perform any medical related data analysis.

#### 2.1.2 [ ] requirements

**PASSED!**

Software [ ] shall use unit's internal [ ] for all actions especially for signal sampling. Standard sampling frequency shall be [ ]

#### 2.1.3 Physical characteristics

**PASSED!**

The software shall be implemented as an embedded system in [ ] language. All the functionality shall be implemented specifically for the hardware provided which is specified in the next chapter.

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Software Validation Report

## Software Validation Report

Software Name :  (Ver )  
Ref No. :

	Drawing	Check by	Approval
Approval	<input type="text"/>	<input type="text"/>	<input type="text"/>

Corp.

# INDEX

1.	Level of Concern .....	7
2.	Software Description .....	8
2.1.	Software Name : <input style="width: 100px; height: 15px;" type="text"/> .....	9
2.1.1.	Software Function .....	9
2.1.2.	<input style="width: 100px; height: 15px;" type="text"/> .....	10
2.2.	Hardware Requirement .....	10
2.2.1.	System requirement .....	10
2.2.2.	<input style="width: 100px; height: 15px;" type="text"/> .....	10
2.3.	Development Language and Program Size .....	11
2.3.1.	Program Language .....	11
2.3.2.	<input style="width: 100px; height: 15px;" type="text"/> .....	11
2.3.3.	Memory Management .....	11
2.3.4.	<input style="width: 100px; height: 15px;" type="text"/> .....	11
3.	Software Requirements Specification .....	12
3.1.	Interface .....	12
3.1.1.	User Interface .....	12
3.1.1.1.	<input style="width: 100px; height: 15px;" type="text"/> .....	12
3.1.2.	Software Interface .....	12
3.1.3.	System Interface .....	12
3.1.4.	<input style="width: 100px; height: 15px;" type="text"/> User Interface .....	13
3.2.	Functional Requirement .....	20
3.2.1.	Function Essential Factor by Module .....	20
3.3.	Description of Process .....	22
3.3.1.	Scanner Module .....	22
3.3.2.	<input style="width: 100px; height: 15px;" type="text"/> Module .....	23
3.3.3.	Database Module .....	23
3.3.4.	<input style="width: 100px; height: 15px;" type="text"/> Module .....	23
3.3.5.	Image Processing Module .....	24
3.3.6.	Viewer Module .....	24
3.3.7.	<input style="width: 100px; height: 15px;" type="text"/> Module .....	24
3.3.8.	Study List Module .....	25
3.3.9.	Transfer Information List Module .....	25
3.3.10.	<input style="width: 100px; height: 15px;" type="text"/> .....	25
3.3.11.	CD / DVD Write Module .....	25
3.3.12.	Patient Registration and Editing Module .....	26
3.4.	Limitation .....	26
3.5.	Recovery .....	26
3.6.	Inside checking .....	27
3.7.	Error and Interrupt .....	27
3.8.	Software Error .....	27
3.9.	Timing .....	27
4.	Application Algorithm and Technology .....	28
4.1.	Calibration .....	28

4.1.1.	Geometry Correction .....	28
4.1.2.	Offset Correction .....	28
4.1.3.	[REDACTED] .....	28
4.1.4.	Defect Correction.....	28
4.2.	Post-processing.....	28
4.2.1.	[REDACTED] .....	29
4.2.2.	Contrast and Edge Enhancement.....	29
4.2.3.	Dynamic Range Reduction .....	29
4.2.4.	[REDACTED] .....	30
4.2.5.	Tone Scaling.....	30
<b>5.</b>	<b>[REDACTED] Design Chart .....</b>	<b>32</b>
5.1.	Software Diagram and Description .....	32
5.1.1.	Description.....	32
5.1.2.	Diagram.....	33
5.2.	Flow Chart.....	34
5.2.1.	Initialization.....	35
5.2.2.	[REDACTED]Registration .....	35
5.2.3.	Image Grabbing.....	35
5.2.4.	[REDACTED]Search .....	35
5.2.5.	Data Selection and Load .....	35
5.2.6.	Image Display.....	35
5.2.7.	[REDACTED] .....	35
5.2.8.	Save and Transmission .....	36
5.2.9.	Finalization .....	36
5.3.	[REDACTED] .....	36
5.3.1.	Data Class Relationship Diagram.....	36
5.3.2.	Interface Class relationship diagram .....	36
5.4.	Database Description .....	37
5.4.1.	[REDACTED] .....	37
5.4.2.	Parameter Database .....	37
<b>6.</b>	<b>Software Design Specification .....</b>	<b>38</b>
6.1.	Total Process [REDACTED] .....	38
6.1.1.	Create Study from Work-List .....	38
6.1.2.	Patient Registration Process .....	38
6.1.3.	[REDACTED] Registration .....	39
6.1.4.	Image Acquisition from Scanner and [REDACTED] .....	39
6.1.5.	Image Processing.....	39
6.1.6.	[REDACTED] .....	39
6.1.7.	Image Transfer and Print.....	39
6.1.8.	Image Storing in CD / DVD.....	40
6.1.9.	Open Current Study by Image Viewer .....	40
6.1.10.	Close Study .....	40
6.1.11.	Verifying Transfer or Print Status in Transfer Status Window. ...	40
6.2.	Software component.....	40


6.3.1.	Scanner Module .....	42
6.3.2.	Calibration .....	51
6.3.3.	Database Module .....	55
6.3.4.	[REDACTED] .....	62
6.3.5.	Image Processing Module .....	84
6.3.6.	Viewer Module .....	90
6.3.7.	Work List Module .....	97
6.3.8.	[REDACTED] Module .....	99
6.3.9.	Transfer Information List Module .....	101
6.3.10.	Print Module .....	103
6.3.11.	[REDACTED] .....	105
6.3.12.	Patient Registration and Editing Module .....	107
<b>7.</b>	<b>Software Development Environment Description .....</b>	<b>109</b>
7.1.	Introduction .....	109
7.2.	Concepts .....	109
7.2.1.	Responsibilities .....	109
7.2.1.1.	Management .....	109
7.3.	Software Life Cycle .....	110
7.4.	Verification .....	111
7.4.1.	[REDACTED] .....	111
7.4.2.	System Verification .....	111
7.4.3.	[REDACTED] .....	111
7.5.	Validation .....	112
7.6.	Amendments .....	112
7.7.	Documentation / Main Product File / Risk Management File .....	112
<b>8.</b>	<b>Software Verification and Validation Test .....</b>	<b>114</b>
8.1.	Verification Test .....	114
8.1.1.	[REDACTED] .....	115
8.1.2.	Calibration Module .....	121
8.1.3.	Database Module .....	124
8.1.4.	[REDACTED] .....	132
8.1.5.	Image Processing Module .....	170
8.1.6.	Viewer Module .....	176
8.1.7.	[REDACTED] .....	187
8.1.8.	Study List Module .....	190
8.1.9.	Transfer Information List Module .....	193
8.1.10.	Print Module .....	195
8.1.11.	CD / DVD Write Module .....	197
8.1.12.	Patient Registration and Editing Module .....	200
8.2.	Validation Test .....	202
8.2.1.	Scanner Module .....	203
8.2.2.	Calibration Module .....	209
8.2.3.	[REDACTED] .....	211
8.2.4.	[REDACTED] .....	217

8.2.8. Study List Module .....	257
8.2.9. Transfer Information List Module .....	259
8.2.10. Print Module .....	261
8.2.11. CD / DVD Write Module.....	263
8.2.12. Patient Registration and Editing Module.....	265
<b>9. Traceability Analysis .....</b>	<b>266</b>
<b>10. Hazard Analysis.....</b>	<b>300</b>
10.1. Severity .....	300
10.2. <input type="text"/> .....	301
10.3. Risk Acceptability.....	302
10.4. Hazard Analysis.....	303
11. <input type="text"/> .....	317
12. <input type="text"/> .....	318
12.1. Scope .....	318
12.2. <input type="text"/> .....	318

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## 6. Software Design Specification

### 6.1. Total Process Flowchart

#### 6.1.1. Create Study from Work-List

- ✓ Retrieve work list from Work-list server by query.
- ✓ Display retrieved work list in list control.
- ✓
- ✓ Start search by input data in search field and press enter.
- ✓ Selected study is started when the study is double clicked by left mouse button or click
- ✓ In the list control,
- ✓

#### 6.1.2. Patient Registration Process

- ✓ New patient can be added through new patient registration process.
- ✓
- ✓ Date of Birth: Input in a form of
- ✓ Display error message and let user re-input when there are errors in compulsory input field or date of birth.

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## 8. Software Verification and Validation Test

### 8.1. Verification Test

✓ Test Categories

- Scanner Module
- Database Module
- 
- Image Processing Module
- Viewer Module
- 
- Study List Module
- Transfer Information List Module
- 
- 
- Patient Registration and

✓ Procedures

- Check to ensure all components are correctly connected.
- Indicate which components have been used; indicate manufacturer and trade / model name of product, version number of software etc.
- Execute the tests in the table below.
- Record the output for each test in the  column.
- 
- After executing all tests, compare the  with the  and then mark either Pass (P) or Fail (F) in the
- If any additional comments need to be made, note the comments in the
- Attach any related notes, references, test data used, and all printed reports with the test report, and indicate the related Test Number on each reference.
- Upon completion of the tests, the executor should sign under
- 

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<div style="border: 2px solid red; padding: 5px; display: inline-block;"> <input type="text"/> Corporation         </div>		<input type="text"/>	
<b>Document title</b> Specification Verification Test Procedure - <input type="text"/>	<b>Document No.</b> T - <input type="text"/>	<b>Revision</b> 6	<b>Page</b> 1 of 16

AUTHORIZATION			
	NAME	SIGNATURE	DATE
Issued by:	<input type="text"/>	<input type="text"/>	<input type="text"/>
Approved by:	<input type="text"/>	<input type="text"/>	<input type="text"/>
Approved by:	<input type="text"/>	<input type="text"/>	<input type="text"/>

REVISION HISTORY		
REVISION #	DESCRIPTION	DATE
1	New Issue	<input type="text"/>
2	Change of tidal volume and rate based on new specifications – revision 7	12/19/ <input type="text"/>
3	Change of rate, pressure and addition of the disconnect alarm	9/20/ <input type="text"/>
4	Sensitivity test changed to match product specifications. Units of measure added to all tests. All CPAP test criteria are changed to $\pm 2$ cm H <sub>2</sub> O. Pass / fail criteria added for disconnect alarm time delay. Low gas supply alarm restoration test added	12/9/ <input type="text"/>
5	Modify test procedure for calibrated rate control	2/7/ <input type="text"/>
6	Modify test procedure for Model S with calibrated rate control	1/11/ <input type="text"/>

<input type="text"/> Corporation		<input type="text"/>	
Document title Specification Verification Test Procedure – <input type="text"/>	Document No. T - <input type="text"/>	Revision 6	Page 2 of 16

## 1.0 Purpose

To test the  Ventilator with the new enclosure and rate control calibration scale. These tests confirm that the device meets the requirements of the product specification.

## 2.0 Scope

ventilator testing

## 3.0 Definitions

- 3.1 Accuracy – How close the output is to desired measured by a referenced measurement device
- 3.2 Precision – Repeatability of a control
- 3.3 Specificity – Effect of one control on another
- 3.4 Sensitivity – How close a control can be set
- 3.5 Linearity – Equal incremental change in a control

## 4.0 Applicable Documents

- 4.1  Rev 8  Ventilator Specification

## 5.0 Procedure

- 5.1. Test conditions to be room temperature, humidity, barometric pressure
- 5.2. The Pneuton is to be set up as follows:
  - 5.2.1. a standard disposable patient circuit with bacteria filter
  - 5.2.2. compressed oxygen supplied from a cylinder. The ability to vary the input pressure from 41 to 66 psi is required.
  - 5.2.3. the patient circuit is attached to the ventilator and a  test lung. The  is set for a compliance of 50 ml/cm H<sub>2</sub>O and resistance of 20 cm H<sub>2</sub>O/L/sec.
  - 5.2.4. a calibrated ventilator tester (FlowAnalyzer ) connected to a PC) is placed between the ventilator and the test lung. This device is used as an external measuring instrument for inspiratory pressure / flow and the generation of graphic pressure / flow traces.

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## Certificate of Registration

The following organization's quality management system has been assessed and registered by Intertek Testing Services NA Ltd., a **CMDCAS recognized registrar**, as conforming to the requirements of:

### ISO 13485:2003

Organization:

Corporation

Main Site: , USA

Certificate Number

Initial Certification Date  
March 1, 2007

Certificate Effective Date  
March 1, 2010

Certificate Expiry Date  
March 1, 2013

**The Quality Management System is applicable to:**

Design, manufacture, distribution and service of Ventilators, CPAP Systems, breathing circuits and other respiratory support accessories.

In the issuance of this certificate, Intertek assumes no liability to any party other than to the Client, and then only in accordance with the agreed upon Certification Agreement.

  
 David S. Muil – General Manager - Canada  
 Intertek Testing Services NA Ltd. – Lachine, QC, Canada

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**Material Safety Data Sheet**  
acc. to ISO/DIS 11014

Printing date

Reviewed on

**1 Identification of substance**

**Product details**

Trade name:

Article number: 154903

Manufacturer/Supplier:

**Emergency information:**

**2 Composition/Data on components**

**Chemical characterization**

Description: Polyamide 12

Dangerous components: Void

**3 Hazards identification**

Hazard description: Not applicable.

**Information pertaining to particular dangers for man and environment:**

The product does not have to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.

**Classification system:**

NFPA ratings (scale 0 - 4)

HMIS-ratings (scale 0 - 4)

HEALTH	0	Health = 0
FIRE	0	Fire = 0
REACTIVITY	0	Reactivity = 0

**4 First aid measures**

General information: No special measures required.

After inhalation: Seek medical treatment in case of complaints.

**After skin contact:**

Do not pull solidified product away from the skin.

After contact with the molten product, cool rapidly with cold water.

**After eye contact:**

Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

After swallowing: If symptoms persist consult doctor.

**5 Fire fighting measures**

**Suitable extinguishing agents:**

CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

(Contd. on page 2)

USA

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**FORT WAYNE METALS**  
An ISO 9001:2000 Registered Company  
Certified by  TUV Rheinland of North America, Inc.

9609 INDIANAPOLIS ROAD  
FORT WAYNE, IN 46809  
(260) 747-4154  
FAX NO. (260) 747-0398  
www.fwmetals.com

---

MANUFACTURERS OF PRECISION FINE STAINLESS STEEL WIRE

Customer:

Order No:

Part No:

Amount Shipped: 315,655 kg

Product Description: 5 mm x 2500 mm

Dia./Thick Tol. Length  
+ 0.0 mm 50 mm  
- .018 mm 0.0 mm

P/N: 55TE19685A01100

FWMR Prod. No:

Shipment No: 2

Date of Shipment:

FWMR Inv. No:

Ti-6Al-4V ELI  
Anid

**Certification of Compliance and Chemical Composition**

This is to certify that the material in this shipment conforms to Specification:  
ASTM F136 and was produced from our Raw Material Lot No. DA-83488 having the following analysis as furnished by our supplier:

<u>Chemical Composition (percent)</u>											
<u>Subject to tolerance and check analysis per specification submitted.</u>											
C	Mn	P	S	Si	Cr	Ni	Mo	Fe	Ti		
.01								18	Bal		
Al	N	Cu	Co	Se	H	O	VANADIUM				
6.23	.01	<input type="text"/>			*.0023	.11	3.99				

Resistance (Ω/m)

Yield Strength (MPa) 827

Yield Load (Kg) 1657

Tensile Strength (MPa) 972

BreakLoad (Kg) 1948

Elongation (%) 22.0

Gage Length (mm) 51

Heat Number T5137G

Tested and Accepted

\*Average final hydrogen analysis  
Centerless ground & polished.  
Beta Transus Temp 1783°F  
Microstructure is acceptable to ASTM  
F136 standard.  
Reduction of Area 43%

Certification By 

F-QM-28 rev

**FORT WAYNE METALS**

An ISO 9001:2000 Registered Company  
 Certified by TÜV Rheinland of North America, Inc.

9509 ARDMORE AVENUE  
 FORT WAYNE, IN 46809  
 (260) 747-4154  
 FAX NO. (260) 747-0398  
 www.fwmetals.com

MANUFACTURERS OF PRECISION STAINLESS STEEL WIRE

**Certification of Compliance**

Customer:   
 Purch Order No:   
 Customer Part No:   
 FWMR Part No: 85TE

FWMR Cert. No.

Amount Shipped	521.6 Pounds	236.594 Kg
	<input type="text"/>	<input type="text"/>

FWMR Prod. No. 144810-01.01

Condition: Anil

Type: TI-6Al-4V ELI

Lot#

.15748x98.425

4 mm x 2500 mm

	+	-		+	-
Dia/Thick Tolerance (in)	0	.0007	Dia/Thick Tolerance (mm)	0	.018
Length Tolerance (in)	2	0	Length Tolerance (mm)	58	0

The material in this shipment is certified to comply with specification number ASTM F136

Date:

Certified By:

*John Beckman Sr.*

**Chemical Composition (percent)** Subject to tolerance and check analysis per specification submitted.  
 The material in this shipment was produced from the following raw material lot numbers whose analysis was furnished by the supplier

Lot Type	TI-6Al-4V ELI				
Heat Number	T8112G				
Raw Material Lot #	83771				
Description					
Carbon	<.01				
Manganese					
Phosphorus					
Sulfur					
Silicon					
Chromium					
Nickel					
Molybdenum					
Selenium					
Titanium	Bal				
Aluminum	6.23				
Nitrogen	.01				
Copper					
Cobalt					
Iron	.18				
Hydrogen	.003				
Oxygen	.11				
VANADIUM-V	4.10				

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195 Museum Road  
Washington, PA 15307  
Phone 724.228.1000  
Fax 724.228.4156

PRODUCT CERTIFICATION

WORK ORDER  
045272

SALES ORDER / RLS

SOLD TO

, Inc

Quality System is registered  
to ISO 9001:2000 and AS9100:2001

CUSTOMER P.O. <input type="text"/>	CUSTOMER PART	QUANTITY 567 lbs	COILS/TAGS 2	LADING NO 00045963	SHIPMENT DATE 12/15/2010					
PART INFORMATION:121DD4ABD01574D TI 6AL 4V ELI Precision Bar - 0.1574" Dia. Specifications: ASTM F 136 Rev: 08e1 Size: 0.1574" Dia. +/- 0.00015" Length: 8' +/- 6"										
Chemical Report (Wt.%)										
Heat 8-33-06036	Al 6.08	V 3.94	C .02	Fe .19	O .11	N .003	Y <.001	Ot E <.05	TOE <.30	Ti BAL
Hydrogen (As Shipped)										
Hydrogen (%) 0.0035										
As Shipped Properties										
UTS #1 (psi) 149100	YS 0.2% #1 (psi) 131400	EL #1 (%) 19	RA #1 (%) 45							
Hardness (HRC) 31										
Metallography Results										
Microstructure <input type="text"/>										
TEST Alpha <input type="text"/>	RESULT Material free of alpha-case									

(v5.2)

Page 1

Date Printed

CUSTOMER P.O. [ ]	CUSTOMER PART	QUANTITY 567 lbs	COILS/TAGS 2	LADING NO [ ]	SHIPMENT DATE [ ]
----------------------	---------------	---------------------	-----------------	------------------	----------------------

**Material Data**

Melt Origin	Bt Tr °F
Russia	1792
Melt Source	Bt Tr °C
[ ]	978

**CERTIFICATION REQUIREMENTS**

**Material Condition: Annealed 704°C 2 Hours A/C - As Shipped**

End of Certification

This is to certify that the above test results are correct as contained in the records of the company.

Certification shall not be reproduced except in full without the approval of Dynamet Inc. The recording of false, fictitious or fraudulent statements or entries on this document may be punished as a felony under Federal Law.

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August ,

Certificate of Analysis

This is to certify that we have performed the analysis of the product below:

Commodity :  POWDER  
 Grade :   
 Lot No. :   
 Analysis Date : August ,

		Result of Analysis	
Y <sub>2</sub> O <sub>3</sub>	wt%	:	5. <input type="text"/>
Al <sub>2</sub> O <sub>3</sub>	wt%	:	0.050
SiO <sub>2</sub>	wt%	:	≤ 0. <input type="text"/>
Fe <sub>2</sub> O <sub>3</sub>	wt%	:	≤ 0.002
Na <sub>2</sub> O	wt%	:	≤ 0. <input type="text"/>
Ig-loss	wt%	:	<input type="text"/> .60

Since the purchaser's actual use of the products is beyond the control of  Corporation and its related, subsidiary and affiliated companies (collectively, "related companies"),  Corporation and its related companies are neither responsible nor liable for results obtained from the use of the products.  Corporation and its related companies make no representations or warranties as to the appropriateness of the products for any use intended or made by the purchaser. Each purchaser must conduct its own testing, safety and regulatory evaluations.

Since this product is a trial product manufactured in accordance with the purchaser's designation,   and its related companies are neither responsible nor liable for any infringement of third parties' rights.

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**TEST RESULT CERTIFICATE**

<b>Sponsor Address</b>	<div style="border: 1px solid black; width: 150px; height: 20px;"></div>	<b>Technical Initiation</b>	7/31/2006
		<b>Technical Completion</b>	8/9/2006
<b>Contact P.O. Number</b>	60446	<b>Report Date</b>	8/11/2006
		<b>Project Number</b>	<div style="border: 1px solid black; width: 50px; height: 15px;"></div>

<b>Test Article</b>	<div style="border: 1px solid black; width: 150px; height: 20px;"></div>	<b>Ratio</b>	4 g/20 mL
<b>Lot/Batch #</b>	<div style="border: 1px solid black; width: 150px; height: 20px;"></div>	<b>Vehicle</b>	<div style="border: 1px solid black; width: 30px; height: 15px;"></div> Water
<b>Study</b>	Physicochemical Test for Plastics – USP	<b>Extraction Conditions</b>	70 ± 2 °C for 24 ± 2 hours
<b>Comments</b>	SAMPLE: Please use surface area extraction ratios. Individual specimen surface area is approximately 50.27cm <sup>2</sup>		

**REFERENCES:** The study was conducted based upon the following references: USP 29, NF 24, 2006. Monograph <661> Containers, Physicochemical Tests–Plastics.

ISO/IEC 17025, 2005, General Requirements for the Competence of Testing and Calibration Laboratories.

**GENERAL PROCEDURE:** The test article was extracted in purified water after rinsing in  water. The following tests were conducted in order to determine physical and chemical properties of the test article's extracts.

**RESULTS:**

TEST	ACCEPTABLE LEVEL	TEST RESULT
Nonvolatile Residue	≤ <div style="border: 1px solid black; width: 15px; height: 15px;"></div> mg	0.4 mg, Meets Criteria
Residue on Ignition	≤ <div style="border: 1px solid black; width: 15px; height: 15px;"></div> mg	Not Applicable
Heavy Metals	≤ <div style="border: 1px solid black; width: 15px; height: 15px;"></div> ppm	< 1 ppm, Meets Criteria
Buffering Capacity	≤ <div style="border: 1px solid black; width: 15px; height: 15px;"></div> mL	0.00 mL, Meets Criteria

**CONCLUSION:** The test article meets the USP Physicochemical Criteria for Plastics based upon the methods employed.

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## 2. 참고문헌 (Reference)

- 1) ISO 11607-1:2006      Packaging for terminally sterilized medical devices - Part 1. Requirements for materials, sterile barrier systems and packaging systems
- 2) ISO 11607-2:2006      Packaging for terminally sterilized medical devices - Part 2. Validation requirements for forming, sealing and assembly processes
- 3) ISO 11737-2:1998      Sterilization of medical devices - Microbiological methods - Part 2. Tests of sterility performed in the validation of a sterilization process
- 4) ASTM F1980              Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
- 5) ASTM D882-09          Tensile Properties of Thin Plastic Sheeting
- 6) ASTM F1929-98         Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
- 7) [ ]                          공정밸리데이션 절차서
- 8) [ ]                          [ ]

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
Monopolar is the most commonly used electrosurgical modality. This is due to its versatility and clinical effectiveness. In monopolar electro-surgery, the active electrode is in the surgical site. The patient return electrode is somewhere else on the patient's body. The current passes through the patient as it completes the circuit from the active electrode to the patient return electrode.

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##### Monopolar Circuit

This picture represents a common monopolar circuit. There are four components to the monopolar circuit: Generator, Active Electrode, Patient, and Patient Return Electrode.

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
[redacted] Coagulators ([redacted])

The proprietary design of [redacted] the shaft temperature as compared to existing [redacted] coagulators and may reduce the potential for oral burns.

Inadvertent thermal injuries during routine otolaryngology procedures can result from inadequate device insulation.<sup>1</sup>

The shaft temperature of [redacted] is lower than that of the competitive devices evaluated (see chart below).<sup>2,3</sup>

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1. Richter [redacted], 2008;72(7):1013-1021.
2. Internal testing of [redacted] coagulator, [redacted] devices. Data on file at [redacted]
3. [redacted] and the [redacted] are used in procedures that require a smaller diameter and do not have the additional insulation.

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---

[ ]

**Performed by:** The performer of the test report is listed below:

Function, Title	Name	Signature	Date
Clinical Liaison	[ ]	[ ]	[ ]

**Document Approvers:** The document approvers are listed below:

Function, Title	Name	Signature	Date
NBP Program Manager	[ ]	[ ]	[ ]

# ANSI/AAMI SP10 Clinical Evaluation Report

## 1 Summary

The [ ] in combination with Mobile [ ] cuffs ([ ]), [ ], Mobile CL Reusable cuffs ([ ]), [ ] and the cradle [ ] was evaluated according ANSI/AAMI SP10:2002/(R)2008+A1:2003/(R)2008 using the auscultatory reference method.

485 valid paired measurements were obtained from 125 subjects including 19 pediatric subjects and compared to indirect cuff/stethoscope auscultation. Two trained observers determined systolic and diastolic blood pressure using the Korotkoff method. The 5<sup>th</sup> Korotkoff sound was used to determine diastolic pressure. The sequential single-arm method was used because the stepwise cuff deflation of the instrument didn't allow simultaneous instrument / observer measurements. Differences between the two observers were verified to meet the criteria of ANSI/AAMI SP10:2002/(R)2008+A1:2003/(R)2008. The measurement results of the device under test had to yield a mean difference versus the averaged observer measurements of no greater than  $\pm 5$  mmHg, with a maximum standard deviation of 8 mmHg.

### Summary data for data evaluation according to ANSI/AAMI SP10:2002/(R)2008 +A1:2003/(R)2008:

The limb circumferences ranged from [ ] cm to 41 cm, with a distribution of 20.8% below 25 cm and 2.8% above 35 cm. The auscultatory systolic pressure ranged from 71.5 mmHg to 153.5 mmHg, with an average of 114,6 mmHg and with a distribution of 19,2% below 100 mmHg and 0.0% above 160 mmHg. The auscultatory diastolic pressure ranged from 40.5 mmHg to 102.0 mmHg, with an average of 72.0 mmHg and with a distribution of 17.6% below 60 mmHg and 0.8% above 100 mmHg. 16.0% of the subjects were measured supine, 54.4% were measure seated and 29.6% were measured standing. The mean errors and standard deviations in the table below demonstrate that the mean error and the standard deviation of all device measurements compared to the respective observer reference measurements are within the limits of the ANSI/AAMI SP10:2002/(R)2008 +A1:2003/(R)2008:

	Method 1 (individual measurements)		Method 2 (subject averages)	
	Systolic [mmHg]	Diastolic [mmHg]	Systolic [mmHg]	Diastolic [mmHg]
Mean Difference	4.10 ( $\pm 5$ )	-1.63 ( $\pm 5$ )	4.11 ( $\pm 5$ )	-1.65 ( $\pm 5$ )
Standard Deviation	6.81 (<8)	6.87 (<8)	4.49 (<5.24)	5.03 (<6.74)
# of measurements	485	485	125	125

The values in ( ) are the limits according to ANSI/AAMI SP10:2002/(R)2008+A1:2003/(R)2008. For method 2 the maximum standard deviation depends on the actual mean difference.

The obtained clinical data together with the laboratory validation data demonstrate that the [ ] in combination with Mobile [ ] cuffs [ ], [ ] cuffs ([ ]), [ ] and the cradle [ ] meets the requirements (mean error of less than  $\pm 5$

# ANSI/AAMI SP10 Clinical Evaluation Report

## 2 Data Collection and Clinical Evaluation Plan

### 2.1 Non-Clinical Data Collection

Site	
Department	

### 2.2 Locations for Clinical Data Collection

Hospital	
Department	

### 2.3 Objectives

The objective of this clinical evaluation was to determine that the overall accuracy of the [redacted] [redacted] meets the requirements as expressed in the ANSI/AAMI SP10:2002/(R)2008 +A1:2003/(R)2008 for adult and pediatric subjects. The auscultatory method is used as the reference method.

For systolic and diastolic pressures, treated separately, the mean difference of paired measurements of the device under test and the auscultatory method shall be  $\pm 5$  mmHg or less, with a standard deviation of 8 mmHg or less.

### 2.4 Legal and Ethical Considerations

This study conforms to the terms of the Declaration of Helsinki (latest revision), [redacted] ~~Device Directive, FDA, applicable IIR requirements, and each site's country, state, and local IRB~~ [redacted] The patient safety of equipment [redacted] [redacted]

Laboratory performance evaluation of the device under test has been conducted using an advanced pressure simulator as a source. The overall accuracy remains to be evaluated and compared under clinical conditions.

There were no additional risks for the volunteer or patient because:

- No invasive procedures, or drawing blood or tissue samples were carried out due to this study
- No additional, or different medical treatment to the subjects, of any kind, was required in conjunction with this study

# ANSI/AAMI SP10 Clinical Evaluation Report

## 3 Results

### Subjects


Overall 125 subjects including 19 pediatric subjects entered the clinical study. All collected data were evaluated.

Number of subjects:	125
Number of measurements:	485

*Table 1 - Number of subjects and paired measurements*

### Range and distribution of limb circumferences

Figure 1 summarizes the distribution of limb circumferences as seen in this clinical evaluation. The limb circumferences range from  cm to 41 cm, with a distribution of 20.8% below 25 cm and 2.8% above 35 cm.

*Figure 1 - Distribution of limb circumference*

### Range and distribution of systolic blood pressure

Figures 2 summarizes the distribution of systolic blood pressures as seen in this clinical evaluation. The auscultatory systolic pressure ranges from 71.5 mmHg to 153.5 mmHg, with an average of 114,6 mmHg and with a distribution of 19,2% below 100 mmHg and 0.0% above 160 mmHg.

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# ANSI/AAMI SP10 Clinical Evaluation Report

## Bland-Altman graphical plots of systolic differences against the systolic reference values


Figure 4 -  graphical display of systolic differences against systolic values determined with the Korotkoff method

## Bland-Altman graphical plots of diastolic differences against the diastolic reference values


Figure 5 -  graphical display of diastolic differences against diastolic values determined with the Korotkoff method

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# ANSI/AAMI SP10 Clinical Evaluation Report

---

## Percentages of readings with differences within 5, 10 and 15 mmHg

	Systolic	Diastolic
NBP readings within $\pm 5$ mmHg of reference	54.6%	64.1%
NBP readings within $\pm 10$ mmHg of reference	81.2%	87.0%
NBP readings within $\pm 15$ mmHg of reference	94.0%	95.3%

**Table 2** - Distribution of paired measurement differences in percentages

## Summary discussion:

The obtained clinical data together with the laboratory validation data demonstrate that the [redacted] [redacted] in combination with [redacted] cuffs [redacted], [redacted], [redacted] cuffs [redacted], [redacted] and the cradle [redacted] meets the specified accuracy (mean error of less than  $\pm 5$  mmHg and a standard deviation of less than 8 mmHg) over the specified blood pressure and pulse rate range.

## 4 Adverse Events

As expected no adverse events caused by the study itself have been reported.

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Artificial Organs  
20(12):1278-1281, Blackwell Science, Inc., Boston  
© 1996 International Society for Artificial Organs

## Clinical Use of Omniflow Vascular Graft as Arteriovenous Bridging Graft for Hemodialysis

Shoei-Shen Wang and Shu-Hsun Chu

Department of Surgery, National Taiwan University, Taipei, Taiwan

**Abstract:** From April 1992 to November 1995, 109 bridging arteriovenous fistulas were performed in 89 patients for hemodialysis at National Taiwan University Hospital. These included 61 Omniflow vascular grafts (OVG) in 47 patients and 48 GoreTex polytetrafluoroethylene (PTFE) grafts in 42 patients. There were 36 male and 53 female patients, whose ages ranged from 24 to 84 years with a mean of  $59 \pm 2$  years. The operation was performed under lidocaine local anesthesia with antibiotic prophylaxis. The anatomic configurations included 50 (45.9%) forearm loop grafts, 19 (17.4%) forearm straight grafts, and 40 (36.7%) thigh loop grafts. Hemodialysis via the bridging graft was started 2 weeks after the operation. The patients received hemodialysis 3 times a week or more when indicated. No patients suffered from limb loss. Pseudoan-

eurysm and graft infection requiring operation occurred 0.96% and 1.92% per patient year, respectively, in the OVG group, and pseudoaneurysm and graft infection both occurred 2.38% per patient in the PTFE group. In the OVG group, the thigh loop graft had better patency than the forearm straight graft. The forearm loop graft had a better patency rate than the forearm straight graft. In the OVG group, the 1 to 4 year patency rate was  $77 \pm 6\%$ ,  $58 \pm 8\%$ ,  $48 \pm 9\%$ , and  $34 \pm 13\%$ , respectively, whereas in the PTFE group, the 1 to 4 year patency rate was  $69 \pm 7\%$ ,  $50 \pm 9\%$ ,  $43 \pm 10\%$ , and  $26 \pm 15\%$ , respectively. Therefore, the Omniflow is a good alternative bridging graft for hemodialysis. **Key Words:** Blood vessel prosthesis—Vascular patency—Omniflow—Arteriovenous fistula—GoreTex—Polytetrafluoroethylene.

In chronic renal failure patients, the autogenous Brescia-Cimino fistula is usually the most popular and satisfactory vascular access for hemodialysis (1,2). When autogenous arteriovenous fistula (AVF) is no longer feasible, the use of a prosthetic conduit to form a bridging AVF is the best alternative. The expanded polytetrafluoroethylene graft was the most commonly used vascular prosthesis (3-5), yet it had a high rate of thrombosis, infection, or puncture complications (6,7). Because the patency rate of biosynthetic grafts was reported to be higher than the patency rate of expanded polytetrafluoroethylene graft (8,9), we introduced the biosynthetic graft as a bridging graft for hemodialysis.

The Omniflow II vascular graft (OVG) is a biosynthetic vascular prosthesis composed of glutaraldehyde-tanned ovine collagen which is grown around a polyester mesh template. The polyester

mesh is an integral part of the graft wall and is encapsulated by Type 3 collagen whereas the outer layers are formed by longitudinally oriented Type 1 collagen (10-12). The collagen may provide the biocompatibility and hemocompatibility whereas the polyester mesh may enhance the durability. Therefore, we started to use the OVG as an arteriovenous bridging graft for hemodialysis in 1992. In this article, we report our 3 year experiences with this biosynthetic composite vascular graft.

### MATERIALS AND METHODS

From April 1992 to November 1995, at National Taiwan University Hospital, 591 operations in 482 patients were performed to create hemodialysis vascular access, including 109 (18.4%) bridging AVF and 482 (81.6%) autogenous AVF. This included all patients with new-onset renal failure and many other patients with complicated access histories who were referred for final care after many operations had been tried elsewhere and had failed. When our operations failed, many patients shifted from hemodialysis to continuous ambulatory peri-

Received August 1996.  
Address correspondence and reprint requests to Dr. Shoei-Shen Wang, National Taiwan University Hospital, No. 7 Chung-Shan South Road, Taipei, Taiwan.

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■ Articles

## ➤ Percutaneous closure of the left atrial appendage versus warfarin therapy for prevention of stroke in patients with atrial fibrillation: a randomised non-inferiority trial

David R Holmes, Vivek Y Reddy, Zoltan G Turi, Shephal K Doshi, Horst Sievert, Maurice Buchbinder, Christopher M Mullin, Peter Sick, for the PROTECT AF Investigators\*

Lancet 2009; 374: 534-42

See Editorial page 501

See Comment page 504

\*Members listed at end of paper

Mayo Clinic College of Medicine, Rochester, MN, USA (Prof D R Holmes MD); Mount Sinai School of Medicine, New York, NY, USA (V Y Reddy MD); Cooper Hospital, Camden, NJ, USA (Z G Turi MD); Pacific Heart Institute/St Johns Hospital, Santa Monica, CA, USA (S K Doshi MD); CardioVascular Center Frankfurt, Sankt Katharinen, Frankfurt, Germany (Prof H Sievert MD); Foundation for Cardiovascular Medicine, La Jolla, CA, USA (M Buchbinder MD); Integra Group, Brooklyn Park, MN, USA (C M Mullin MS); and Krankenhaus der Barmherzigen Brüder, Regensburg, Germany (P Sick MD)

### Summary

**Background** In patients with non-valvular atrial fibrillation, embolic stroke is thought to be associated with left atrial appendage (LAA) thrombi. We assessed the efficacy and safety of percutaneous closure of the LAA for prevention of stroke compared with warfarin treatment in patients with atrial fibrillation.

**Methods** Adult patients with non-valvular atrial fibrillation were eligible for inclusion in this multicentre, randomised non-inferiority trial if they had at least one of the following: previous stroke or transient ischaemic attack, congestive heart failure, diabetes, hypertension, or were 75 years or older. 707 eligible patients were randomly assigned in a 2:1 ratio by computer-generated randomisation sequence to percutaneous closure of the LAA and subsequent discontinuation of warfarin (intervention; n=463) or to warfarin treatment with a target international normalised ratio between 2.0 and 3.0 (control; n=244). Efficacy was assessed by a primary composite endpoint of stroke, cardiovascular death, and systemic embolism. We selected a one-sided probability criterion of non-inferiority for the intervention of at least 97.5%, by use of a two-fold non-inferiority margin. Serious adverse events that constituted the primary endpoint for safety included major bleeding, pericardial effusion, and device embolisation. Analysis was by intention to treat. This study is registered with Clinicaltrials.gov, number NCT00129545.

**Findings** At 1065 patient-years of follow-up, the primary efficacy event rate was 3.0 per 100 patient-years (95% credible interval [CrI] 1.9-4.5) in the intervention group and 4.9 per 100 patient-years (2.8-7.1) in the control group (rate ratio [RR] 0.62, 95% CrI 0.35-1.25). The probability of non-inferiority of the intervention was more than 99.9%. Primary safety events were more frequent in the intervention group than in the control group (7.4 per 100 patient-years, 95% CrI 5.5-9.7, vs 4.4 per 100 patient-years, 95% CrI 2.5-6.7; RR 1.69, 1.01-3.19).

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# Intracorneal inlay to correct presbyopia: Long-term results

Ömer F. Yılmaz, MD, Neşe Alagöz, MD, Gökhan Pekel, MD, Engin Azman, MD,  
Ebru F. Aksoy, MD, Hanefi Çakır, MD, Ercüment Bozkurt, MD, Ahmet Demirok, MD

**PURPOSE:** To evaluate the long-term visual results of Acufocus ACI-7000 (now Kamra) intracorneal inlay implantation in presbyopic phakic patients.

**SETTING:** Beyoğlu Eye Training and Research Hospital, Istanbul, Turkey.

**DESIGN:** Clinical trial.

**METHODS:** This study comprised patients with emmetropic or post-laser in situ keratomileusis (LASIK) presbyopia. Patients had an uncorrected near visual acuity (UNVA) of 20/40 or worse, correctable to 20/25 or better at distance. The inlay was implanted on the stromal bed after the LASIK flap was relifted or a flap created. The inlay was centered on the presurgical position of the first Purkinje reflex. The main outcome measures were distance and near vision and the complication rate.

**RESULTS:** The study enrolled 39 patients aged 45 to 60 years. At the 4-year follow-up, all patients (N = 22) had 2 or more lines of improvement in UNVA with no significant loss in distance vision. The mean final UNVA was 20/20 (Jaeger [J1]); 96% of patients could read J3 or better. The uncorrected distance acuity was 20/40 or better in all eyes. Five patients had cataract progression, and 2 had a change in refractive status. No eye with an intracorneal inlay had intraoperative complications during cataract extraction. Four inlays were explanted during the study. There were no severe corneal complications that affected final vision.

**CONCLUSION:** Intracorneal inlay implantation was an effective, safe, and reversible procedure for the long-term surgical treatment of presbyopia.

**Financial Disclosure:** No author has a financial or proprietary interest in any material or method mentioned. Additional disclosure is found in the footnotes.

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**DEPARTMENT OF HEALTH & HUMAN SERVICES** Public Health Service

Food and Drug Administration  
10903 New Hampshire Ave.  
Silver Spring, MD 20993

**Certificate No. 4317-8-2010**


CERTIFICATE TO FOREIGN GOVERNMENT

In order to allow the importation of United States products into foreign countries, the U.S. Food and Drug Administration (FDA) certifies the following information concerning the product(s) to be exported listed below:

<u>Name of Product(s)</u>	<u>Name of Manufacturer/Distributor Address</u>
See Attached List (1 page)	


The product(s) described above (and the manufacturing/distribution site(s) which produces/distributes it) is subject to the jurisdiction of the FDA under the Federal Food, Drug, and Cosmetic Act.

It is certified that the above product(s) may be marketed in, and legally exported from, the United States of America at this time. The manufacturing plant(s) in which the product(s) is produced is subject to periodic inspections. The last such inspection showed that the plant(s), at that time, appeared to be in substantial compliance with current good manufacturing practice requirements for the products(s) listed above.

  
Chester T. Reynolds  
Acting Deputy Director,  
Division of Risk Management Operations  
Office of Compliance  
Center for Devices and Radiological Health

**This certificate expires 24 months from the date notarized.**

COUNTY OF MONTGOMERY  
STATE OF MARYLAND

Subscribed and sworn to before me this 21 day of sept month 2010 year. '   

Regulatory Status

The  has been approved for commercial distribution in the following countries or regulatory jurisdictions:

<b>Country or Regulatory Jurisdiction</b>	<b>Commercial Approval</b>
European Union (CE Mark)	November 2005
Bahamas	<input type="checkbox"/> 2006
Israel	October 2006
South Africa	March 2007
Saudi Arabia	October 2007
Bulgaria	February 2008 <sup>1</sup>
United States	April 2008 <input type="checkbox"/> January 2010 <input type="checkbox"/>
Canada	May 2009 <input type="checkbox"/> October 2010 <input type="checkbox"/>
Thailand	July 2009
Taiwan	May 2010
Australia	<input type="checkbox"/>
Turkey	March 2011
Russia	March 2011

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### Section 8: Regulatory Information

#### 1.0 US Regulatory Approvals

[Redacted] was approved by FDA under premarket approval (PMA) application P080032 on April 27, 2010. A copy of the approval order can be found under [Redacted]

#### 2.0 Canada Regulatory Approvals

[Redacted] was issued a Class III Medical Device Licence by Health Canada (Licence No. [Redacted]) by 26, 2010, with an update to include the Accessory Kit on September 19, 2011. A copy of the current device licence can be found under [Redacted]

#### 3.0 EU Regulatory Approvals

[Redacted] System is supported by Technical File TF001.

##### 3.1 Declaration of Conformity

Reference PDM # 90677941 for the Declaration of Conformity.

##### 3.2 EC-Certificate

Reference PDM # 90677944 for the Annex II.3 EC-Certificate.

#### 4.0 Quality System Certificate

Reference PDM # 90677944 for the EN ISO 13485:2003 certificate.

Reference PDM # 90677944 for the CMDCAS ISO 13485:2003 certificate.

#### 5.0 List of Standards to which the Product Complies

A list of standards to which the Alair System complies can be found in Attachment 1.

#### 6.0 CFG/CFS

Reference PDM # 90743333 for the Certificate to Foreign Government for the "severe asthma" indication.

Reference PDF # 90743331 for the Certificate of Free Sale for the "asthma" indication.

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DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

Ms. [Redacted]  
Vice [Redacted]ory Affairs  
Asth [Redacted]  
888 [Redacted]oor  
Sun [Redacted]4089

APR 27 2010

Re:

[Redacted]

Dear Ms. Brown:

The Center for Devices and Radiological Health (CDRH) of the Food and Drug Administration (FDA) has completed its review of your premarket approval application (PMA) for the [Redacted] treatment of se [Redacted] controlled with inhaled corticosteroids and long acting beta agonists. We are pleased to inform you that the PMA is approved. You may begin commercial distribution of the device in accordance with the conditions of approval described below.

The sale and distribution of this device are restricted to prescription use in accordance with 21 CFR 801.109 and under section 515(d)(1)(B)(ii) of the Federal Food, Drug, and Cosmetic Act (the act). The device is further restricted under section 515(d)(1)(B)(ii) of the act insofar as the labeling must specify the specific training or experience practitioners need in order to use the device. FDA has determined that these restrictions on sale and distribution are necessary to provide reasonable assurance of the safety and effectiveness of the device. Your device is therefore a restricted device subject to the requirements in sections 502(q) and (r) of the act, in addition to the many other FDA requirements governing the manufacture, distribution, and marketing of devices.

Expiration dating for this device has been established and approved at one year.

Continued approval of this PMA is contingent upon the submission of periodic reports, required under 21 CFR 814.84, at intervals of one year from the date of approval of the original PMA. Two copies of this report, identified as "Annual Report" and bearing the applicable PMA reference number, should be submitted to the address below. The Annual Report should indicate the beginning and ending date of the period covered by the report and should include the information required by 21 CFR 814.84.

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IFIKAT ◆ CERTIFICATE ◆ 認証証書 ◆ СЕРТИФИКАТ ◆ CERTIFICADO ◆ CERTIFICAT

Product Service

## EC Certificate

### Full Quality Assurance System

Directive 93/42/EEC on Medical Devices (MDD), Annex II excluding (4)  
(Devices in Class IIa, IIb or III)

No. G1 11 06 77336 001

Manufacturer:

EC-Representative:

S.A.

Product  
Category(ies):

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for design, manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex II. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class III devices an additional Annex II (4) certificate is mandatory. See also notes overleaf.

Report No.:

Valid from:

2011-08-01

Valid until:

2014-04-30

Hans-Heiner Junker

Date, 2011-07-21

TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123

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**Date de modification:**

**Device Class/Classe de l'instrument: 3**

This Licence is issued in accordance with the  
Medical Devices Regulations, Section 36,  
for the following medical device:

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**Licence Type/Type d'homologation:**

System / Système

**Intended Purpose/Usage/Condition(s):**

LICENCE WITH CONDITIONS

**Reason for Amendment/Raison de la modification**

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순번	인정 번호	기관명	공통규격 (전기안전)	보조규격	전자파규격	개별규격
8	KT123	한국에스지에스 (주) 군포시험소	CAN/CSA-C22.2No.6 0601.1:2008 UL60601-1:2003(R2006) EN60601-1:1990+A1:1993+A2:1995 EN60601-1:2006 IEC60601-1:1998+A1:1991+A2:1995 IEC60601-1:2005(ThirdEdition) IEC60601-1:2005(ThirdEdition) KSPIEC60601-1:2002 KSCIEC60601-1:2008	EN60601-1-1:2001 EN60601-1-3:1994 EN60601-1-8:2004+A1:2006 IEC60601-1-1:2000 IEC60601-1-8:2003A1:2006 (Firstedition) KSCIEC60601-1-1:2002 KSCIEC60601-1-3:2002 KSCIEC60601-1-8:2005	EN60601-1-2:2007 IEC60601-1-2:2007	EN60601-2-10:2000+Am1:2001 EN60601-2-10:2000+Am1:2001 EN60601-2-2:2001 EN60601-2-2:2007 EN60601-2-22:1996 EN60601-2-27:2006 EN60601-2-28:1998 EN60601-2-32:1994 EN60601-2-35:1996 EN60601-2-38:1996+A1:2000 EN60601-2-42:2003 EN60601-2-44:2001+A1:2003 EN60601-2-49:2001 EN60601-2-7:1998 EN60601-2-9:1996 IEC60601-2-10:1987(FirstEd)+Am1:2001 IEC60601-2-2:1998(Thirdedition) IEC60601-2-2:2003(Fourthedition) IEC60601-2-22:1996(Secondedition) IEC60601-2-2-28:1998(Firstedition) IEC60601-2-27:2003(Secondedition) IEC60601-2-32:1994(Firstedition) IEC60601-2-35:1996(Firstedition) IEC60601-2-38:1996(Firstedition)+A1:1999 IEC60601-2-42:2002(Secondedition) IEC60601-2-44:2001+A1:2002(Ediior21) IEC60601-2-49:2001(Firstedition) IEC60601-2-7:1998(Secondedition) IEC60601-2-9:1996(Secondedition) KSCIEC60601-2-10:2002 KSCIEC60601-2-22:2002 KSCIEC60601-2-27:2005 KSCIEC60601-2-28:2002 KSCIEC60601-2-32:2002 KSCIEC60601-2-35:2002 KSCIEC60601-2-38:2002 KSCIEC60601-2-42:2003 KSCIEC60601-2-44:2005 KSCIEC60601-2-49:2005 KSCIEC60601-2-7:2002 KSCIEC60601-2-9:2002 KSPIEC60601-2-2:2002

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