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Introduction

Read the entire manual before using the device.

Indications for use

[] is intended to provide ventilation for non-dependent, spontaneously breathing adult and pediatric patients [] with respiratory insufficiency, or respiratory failure, with or without obstructive sleep apnea. The device is for non-invasive use, or invasive use with an uncuffed or deflated tracheostomy.

Contraindications

[] is contraindicated in patients who are unable to endure more than brief interruptions in ventilation. [] is not a life support ventilator.

The use of the device may be contraindicated in patients with:

- pneumothorax or pneumomediastinum
- pathologically low blood pressure, particularly if associated with intravascular volume depletion
- cerebrospinal fluid leak, recent cranial surgery or trauma
- severe bullous lung disease
- dehydration.

The use of the [] or pulse oximetry [] is contraindicated in an MRI environment.

Adverse effects

Patients should report unusual chest pain, severe headache or increased breathlessness.

The following side effects may arise during the course of noninvasive ventilation with the device:

- drying of the nose, mouth or throat
- nosebleed
- bloating
- ear or sinus discomfort
- eye irritation
- skin rashes.

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Soft Tissue Fixation

HA Suture Anchor

Available in a variety of size options with a choice of suture configurations, Suture Anchors can be used for primary or double row repairs and for the medial row of a suture bridging technique.

Absorbable

Our absorbable HA suture anchor combine with hydroxyapatite (HA) to the mineral compound of bone. Once absorbed, it has the fracture resistance of human cortical bone.* HA is a natural bone mineral shown to promote bone replacement.

REF	Description
	4.5 mm Anchor with two #2 sutures (white, COBRAID-blue), sterile
	4.5 mm Anchor with two #2 sutures (COBRAID-blue, COBRAID-black), sterile
	5.5 mm Anchor with two #2 sutures (white, COBRAID-blue), sterile
	5.5 mm Anchor with two #2 sutures (COBRAID-blue, COBRAID-black), sterile
	5.5 mm Anchor with three #2 sutures (white, COBRAID-blue, COBRAID-black)
	6.5 mm Anchor with two #2 sutures (white, COBRAID-blue), sterile
	6.5 mm Anchor with two #2 sutures (COBRAID-blue, COBRAID-black), sterile
	6.5 mm Anchor with three #2 sutures (white, COBRAID-blue, COBRAID-black)

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System. Skeletal implants for the orthodontic movement of teeth.

Overview

Custom is intended to be implanted intraorally and used as an anchor for orthodontic procedures. It includes anchor screws, anchor plates, screws for plate fixation, instruments, and a module case for storage and sterilization.

- Provides fixed anchorage for improved control of tooth movement
- Compatible with a variety of orthodontic devices including archwires, elastics and springs
- Allows immediate loading
- Eliminates the need for extraoral anchorage (headgear)
- Anchor plates can be adapted to the patient's bony anatomy
- Manufactured from commercially pure (CP) titanium and titanium alloy*

*Ti-6Al-7Nb

Anchor plate

Profile of an archwire in an anchor screw

Anchor screw

Anchor plates

- Allow placement away from tooth roots
- Can be contoured and/or trimmed to conform to patient anatomy
- Implanted using up to five self-drilling titanium screws for stability
- Three designs

Mesh design

Orthodontists can mount preferred dental bracket to anchors using standard adhesive.

Bracket design

Simulates orthodontic bracket and allows attachment of various orthodontic appliances

Domed design

Allows attachment of archwires from three different vectors, or attachment of elastics or springs

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INDICATIONS FOR USE

is indicated for use in the following procedures for reattachment of soft tissue to bone:

Shoulder

- Rotator Cuff
- Biceps Tenodesis

CONTRAINDICATIONS

1. Procedures other than those listed in the indications section.
2. Pathologic conditions of bone such as cystic changes or severe osteopenia that would impair its ability to securely fix the r.
3. Pathologic changes in the soft tissues being fixated to bone that would prevent their secure fixation by the r.
4. Comminuted bone surface that would militate against secure fixation of the r.
5. Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitation, previous infection, etc.
6. Conditions which tend to limit the patient's ability to restrict activities or follow directions during the healing period.
7. The r is not designed for and should never be used to attach artificial ligaments.

GENERAL DEVICE DESCRIPTION

The proposed (See Figure 1) is a one piece implantable cannulated, threaded anchor designed to secure soft tissue to bone in the shoulder. The anchor is provided loaded on a disposable Inserter Driver, with a Threader Tab (a Wire Kite is attached to the device), and a Ramp. The proposed are offered in one size, namely 4.75 mm.

The proposed is molded from absorbable material similar to the (Reference Number:).

The other proposed r (Provided with Absorbable Orthocord Sutures) is molded from non-absorbable PEEK (Polyetheretherketone) material similar to the (Reference Number:).

The proposed l : will be provided packaged with or without sutures. The packaging may include up to one strand of partially absorbable Orthocord () suture.

The proposed : will be provided in one size with major outer diameter of 4.75 x 17.5 mm.

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Clinical Research J Korean Neurosurg Soc 27 : 1370-1378, 1998

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=Abstract=

for the Treatment of Degenerative Lumbar Disease

Degenerative involution of the spine causes destruction of spinal stabilizer which consists of bone, ligament, joint capsule, and disc, which substantially leads to hypermobility and instability of the spine.

Generally the hard fixation system has been used for the treatment of lumbar instability. However, it has many complications, including screw loosening, screw fracture, and instability on, above, and below the fusion segment. These complications of the hard fixation system has brought to the invention of a more physiologic fixation device, the soft fixation system.

We have used the Graf soft fixation system as an instrument for degenerative lumbar disease. 106 cases were operated between August, 1993 and March, 1996. The clinical assessments, radiologic findings, and operative results were analyzed.

The height of disc space significantly increased from 8.84mm to 9.84mm on L3/4, 9.28mm to 10.13mm on L4/5, and 9.44mm to 10.47mm on L5/S1. Flexion instability changed from -6.9° to 5.5° on L3/4, -7.45° to 5.04° on L4/5, -2.09° to 10.81° on L5/S1, translation instability was corrected from 16.8% to 14.9% on L3/4, 19.9% to 12.4% on L4/5, 27.1% to 20.1% on L5/S1 after Graf soft fixation.

The clinical results were as follows : excellent in 56%, 27% good, 9% fair, and 8% poor.

These results suggest that Graf soft fixation system for degenerative lumbar disease would not only be useful and effective, but also be safer in terms of unwanted complications of the hard fixation system.

KEY WORDS : Degenerative lumbar disease • Hard fixation system • Soft fixation system • Graf system • Instability.

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A Principle of Laser Therapy and it's Application

- ABSTRACT -

We discuss the simple concepts of laser physics, absorption of laser radiation by tissue and current application in medical and physical therapy field. The continuing development in the field of laser technology has encouraged collaboration between clinician and physicist. This can be reached through more exactly controlled and regulated laser systems for human. In medicine, surgery and physical therapy a laser should be used when it can be shown clearly established techniques.

The purpose of this article is to provide the connection between pure laser physics and the other clinical application.

Key Words : Laser therapy, Physical therapy

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Preface

Effects of Federal Policies on Extracorporeal Shock Wave Lithotripsy Case Study 36 in the Health Technology Case Study Series. This case study has been prepared in connection with OTA's project on *Payment for Physician Services: Strategies for Medicare*, which was mandated by the Deficit Reduction Act of 1984. The House Committee on Ways and Means and Energy and Commerce and the Senate Committee on Finance have jurisdiction over that part of the law. The Senate Special Committee on Aging also requested the study of physician payment.

Case studies are designed to fulfill two functions. The primary purpose is to provide specific information that can be used in forming general conclusions regarding broader policy issues. The first 19 cases in the Health Technology Case Study Series, for example, were conducted in conjunction with the overall project on *The Implications of*

By examining the 19 cases as a group and looking for common problems or strengths in the techniques of cost-effectiveness or cost-benefit analysis, the study was able to better analyze the potential contribution that those techniques might make to the management of medical technology and health care costs and quality.

The second function of the case studies is to provide useful information on the specific technologies covered. The design and the funding levels of most of the case studies are such that they should be read primarily in the context of the associated overall projects. Nevertheless, in many instances, the case studies do represent extensive reviews of the literature on the efficacy, safety, and costs of the specific technologies and as such can stand on their own as a useful contribution to the field.

Case studies are prepared in some instances because they have been specifically requested by congressional committees and in others because they have been selected through an extensive review process involving staff and consultations with the congressional staffs, advisory panel to the associated overall project, the Health Program Advisory Committee, and other experts in various fields. Selection criteria were developed to ensure that case studies provide the following:

- examples of types of technologies by func-

tion (preventive, diagnostic, therapeutic, and rehabilitative);

- examples of types of technologies by physical nature (drugs, devices, and procedures);
- examples of technologies in different stages of development and diffusion (new, emerging, and established);
- examples from different areas of medicine (e.g., general medical practice, pediatrics, radiology, and surgery);
- examples addressing medical problems that are important because of their high frequency or significant impacts (e.g., cost);
- examples of technologies with associated high costs either because of high volume (for low-cost technologies) or high individual costs;
- examples that could provide information material relating to the broader policy and methodological issues being examined in the particular overall project; and
- examples with sufficient scientific literature.

Case studies were either prepared by staff, commissioned by staff and performed under contract by experts (generally in academia), or written by staff on the basis of contractors' papers.

Each subject of each case study is reviewed extensively by staff and by members of the advisory panel to the associated project. For commissioned cases, comments are provided to authors, along with suggestions for revisions. Subsequent drafts are sent to numerous experts for review and comment. Each case is seen by at least 30 reviewers, and sometimes by 80 or more outside reviewers. These individuals may be from relevant Government agencies, professional societies, consumer and public interest groups, medical practice, and academic medicine. Academicians such as economists, sociologists, decision analysts, biologists, and so forth, as appropriate, also review the cases.

Although cases are not statements of official position, the review process is designed to satisfy the concern with each case study's scientific quality and objectivity. During the various stages of the review and revision process, therefore, the process encourages, and to the extent possible requires, authors to present balanced information and recognize divergent points of view.

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2. Product description

2.1 System components

Designation	Art. no.
Sagittal saw	
Cover	
Rechargeable battery	
Sterile funnel	
Instructions for use	

2.2 Components necessary for use

Designation	Art. no.
Sagittal saw	
Cover	
Sterile funnel	
Saw blade according to indication	

2.3 Intended use

The is used for intraoperative separation of bones and cartilage in orthopedic surgery, traumatology and cardiothoracic surgery.

2.4 Operating principle

The incorporates an electric motor supplied by a rechargeable battery, and electronic circuitry for controlling the rotational speed of the motor. The oscillation of the saw can be adjusted continuously by pressing oscillation frequency control trigger 5.

The converts the motor rotation into oscillation of the saw blade. The bone is separated by the oscillating movement of the saw blade.

No keys or wrenches are required for coupling the saw blades. The pivotable saw head 2 of can be set to eight different angular positions, in steps of 45°.

3. Working with the Ac

3.1 System set-up

Connecting the accessories

Combinations of accessories that are not mentioned in the ed if the: and safety if th safety cha your omer

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INTENDED USE

CAUTION:

Federal law restricts this device to sale by or on the order of a physician.

The _____ supplies _____ by applying sequential compression to the lower extremity through a variety of sleeves. The _____ is both small and mobile. During prevention treatment the patient can easily move about without medical staff intervention. The system supplies several treatment options and sleeve combinations.

One or two sleeves (any combination) can be connected to the device at the same time.

All the components of the system are _____-Free.

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behind ear thermometer instruction manual

The behind ear thermometer has been carefully developed for accurate, safe and fast temperature measurements.

However, as with any thermometer, proper technique is critical in obtaining accurate temperatures. Therefore, read the use instructions carefully and thoroughly.

Warnings and precautions

1. Do not use thermometer for 20 minutes on skin if exposed to hot or cold temperatures such as direct sunlight, fireplace heat, cold (compress), direct air conditioning air flow, swimming, bathing, etc.
2. Do not place the thermometer on scarred tissue or tissue compromised by skin disorders.
3. Do not use on patients in trauma or patients treated with certain drug therapies.
4. Do not take temperature from the side on which the patient was laying on - use the opposite side.
5. Do not expose the thermometer to temperature extremes below -20°C or over 50°C nor excessive humidity (> 85% non-condensing relative humidity).
6. Do not use this thermometer if the operating ambient temperature range is below 10C° or over 40C° .
7. Do not mount the battery with the positive sign facing down. The thermometer will not work.
8. This thermometer is intended for household use only.
9. Use of this thermometer is not intended as a substitute for consultation with your physician.
10. Readings from the right side of the patient's head may differ from the left - therefore it is preferable to take all temperatures on the same side of the head.

How does behind ear thermometer work?

behind ear thermometer is a hand-held, battery powered device designed to measure human body temperature with a special sensor that accurately and gently detects heat flow through the skin directly behind the ear lobe using the heat conduction effect and prediction techniques. The thermometer calculates heat flow through the skin and by use of a mathematical algorithm estimates the internal body temperature. The temperature reading range of the thermometer is from 34°C to 42.2 °C and the reading time is about 1 second.

Why measure behind the ear?

The skin area behind the earlobe is over the posterior auricular artery which is a branch of the carotid artery. The carotid artery carries blood to the brain and therefore the site is the best external place to measure temperature.

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Description of the device's principle of operation

is an adjunct hemostatic agent for use in surgical procedures, when control of bleeding from capillary, venous, or arteriolar vessels by pressure, ligature and other conventional procedures is either ineffective or impractical.

consists of hydrophilic particles that rapidly absorb water from the blood and initiate a dehydration process. Blood solids (platelets, blood cells) and blood proteins are concentrated so that natural hemostasis is accelerated. In contact to blood, the hemostatic powder forms a gel matrix that acts as an instant barrier against further bleeding, independent of the patients clotting status. Absorption normally requires several days and depends amongst others on the amount applied and the site of use.

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FOREWORD

GENERAL INFORMATION

The _____ unit and the accompanying components and individual elements fulfill, as a unit, the currently valid safety standards and comply with the stipulations of _____ and the medical products regulations.

The unit and its external components (accessory elements) are safe if used properly and in compliance with the explanations and instructions provided in this documentation. Nevertheless, the unit or its external components can pose dangers. Therefore, we urgently recommend that anyone operating the shortwave diathermy unit become aware of the potential dangers of the unit and its external components before beginning work.

Please read and observe all safety instructions in this operating manual.

The following features are available on the _____ unit:

Clear Touch Screen LCD

Allows the operator to select an option on the screen by simply touching it. Guides the operator through the setup process providing continuous feedback about treatment settings. Gives the operator optimal visibility during attended procedures. Allows the operator to turn on the unit simply by pressing the LCD.

Thermal Dosimetry

Allows the operator to set up a treatment using the four treatment dosages according to Schliephake.

User Protocols

User protocols allow you to set, save, and change the parameters of each program (protocol) in order to tailor it to meet your patients' specific needs. You may create up to 200 user protocols.

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



DESCRIPTION OF SYMBOLS

-  Consult Instructions for Use
-  Product Reference Number
-  Prescription Use Only
-  Serial Number
-  Lot Number
-  Type BF Applied Part
-  Waste electrical and electronic equipment (WEEE) should not be disposed as unsorted municipal waste; WEEE should be collected separately.
-  Direct Current
-  9V Alkaline Battery
-  Classified by Underwriters Laboratories Medical Equipment With Respect to Electric Shock, Fire, and Mechanical Hazards Only. In accordance with UL 60601-1, CAN/CSA C22.2 No. 601.1, EN 60601-1, IEC 60601-1, and IEC 60601-2-10.
-  Manufacturer
-  Authorized Representative in European Community

INDICATIONS FOR USE

The Urgent® PC Neuromodulation System is intended to treat patients suffering from urinary urgency, urinary frequency, and urge incontinence. The Urgent PC is also indicated for the treatment of fecal incontinence.

PRODUCT DESCRIPTION

The Urgent PC Neuromodulation System is a minimally invasive neuromodulation system designed to deliver retrograde access to the sacral nerve through percutaneous electrical stimulation of the tibial nerve. The method of treatment is referred to as Percutaneous Tibial Nerve Stimulation (PTNS).

4

The Urgent PC Neuromodulation System is a combination of the Urgent PC Stimulator (Stimulator) and the Urgent PC Stimulation Lead Set (Lead Set). The Stimulator and Lead Set are sold separately.

The Stimulator is a battery powered, external pulse generator and is designed, constructed, and manufactured for multiple use. The Stimulator is to be used only in conjunction with the single-use Lead Set. The Lead Set (comprised of the Lead Wire, Needle Electrode, and Alcohol Pad) transfers the electrical current from the Stimulator to the tibial nerve via the Needle Electrode. The only components of the Urgent PC Neuromodulation System provided sterile are the Needle Electrodes.

CONTRAINDICATIONS FOR USE

1. In order for treatment to be effective and to avoid any possible problems or complications, the device is contraindicated for use on patients who have the following history or conditions:
 - Patients with pacemakers or implantable defibrillators
 - Patients prone to excessive bleeding
 - Patients with nerve damage that could impact either percutaneous tibial nerve or pelvic floor function
2. The Stimulator is not intended for intra-cardiac or trans-thoracic use.
3. Do not use the Stimulator on patients who are pregnant or planning to become pregnant while using this product.
4. Concurrent use of medical monitoring equipment during stimulation is not recommended.
5. This device is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

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Cages for Vertebral Body Replacement in the Thoracolumbar Spine

Study Design. An *in vitro* biomechanical study of expandable cages for vertebral body replacement in the human thoracolumbar spine.

Objectives. The purpose of this study was to compare the *in vitro* biomechanical properties of 3 different expandable cages with a nonexpandable cage.

Summary and Background Data. Recently, there has been a rapid increase in the use and the commercial availability of expandable cages for vertebral body replacement in the thoracolumbar spine. Although all 3 expandable cages, evaluated in this study, are approved for clinical use in Europe, little information is available concerning the biomechanical properties of these implants.

Material and Methods. Thirty-two human thoracolumbar spines (T11 to L3) were tested in flexion, extension, axial rotation, and lateral bending with a nondestructive loading technique using an unconstrained testing apparatus. Three-dimensional displacement was measured using an optical measurement system. First, all motion segments were tested intact. After complete corpectomy of L1, cages were implanted according to producer's information. The following implants (n = 8/group) were tested: 1) meshed titanium cage (nonexpandable cage, DePuy AcroMed); 2) X-tenz (expandable cage, DePuy AcroMed); 3) Synex (expandable Cage, Synthes); and 4) VBR (expandable cage, Ulrich). Finally, posterior stabilization using the Universal Spine System (Synthes), posterior-anterior stabilization using the Universal Spine System (Synthes), and anterior plating (Locking Compression Plate, Synthes) was applied to each test specimen. The mean apparent stiffness values, range of motion, and neutral and elastic zone were calculated from the corresponding load-displacement curves.

Results. No significant differences could be determined between the *in vitro* biomechanical properties of expandable and nonexpandable cages. In comparison to the intact motion segment, isolated anterior stabilization

using cages and anterior plating significantly decreased stiffness and increased range of motion in all directions. In contrast, additional posterior stabilization significantly increased stiffness and decreased range of motion in all directions compared to the intact motion segment. The combined anterior-posterior stabilization demonstrated greatest stiffness results.

Conclusion. Biomechanical results indicate that design variations of expandable cages for vertebral body replacement are of little importance. Additionally, no significant difference could be determined between the biomechanical properties of expandable and nonexpandable cages. After corpectomy, isolated implantation of expandable cages plus anterior plating was not able to restore normal stability of the motion segment. Therefore, isolated anterior stabilization using cages plus Locking Compression Plate should not be used for vertebral body replacement in the thoracolumbar spine.

Key words: vertebral body replacement, biomechanics, expandable cages. *Spine* 2004;29:1413-1419

From Unfall- und Wiederherstellungschirurgie, Universitätsklinikum Charité der Humboldt Universität Berlin, Campus Virchow-Klinikum, Berlin, Germany.
 Acknowledgment date: August 29, 2002. First revision date: January 3, 2003. Second revision date: June 13, 2003. Acceptance date: August 29, 2003.
 The legal regulatory status of the device(s)/drug(s) that is/are the subject of this manuscript is not applicable in my country. The device(s)/drug(s) that is/are the subject of this manuscript is/are not FDA-approved for this indication and is/are not commercially available in the United States. The device(s)/drug(s) is/are FDA-approved or approved by corresponding national agency for this indication. No funds were received in support of this work. No benefits in any form have been or will be received from a commercial party related directly or indirectly to the subject of this manuscript.
 Address correspondence and reprint requests to Robert Pflugmacher, MD, Unfall- und Wiederherstellungschirurgie, Universitätsklinikum Charité der Humboldt Universität Berlin, Campus Virchow Klinikum, Augustenburgerplatz 1, Berlin 13353, Germany; E-mail: robert.pflugmacher@charite.de

Various incidents such as fractures, tumors, or infections affect the vertebral body and lead to instability of the spine. A posterior stabilization with an angular stable internal fixation device is a widely accepted surgical treatment for instability of the thoracic and lumbar spine.¹⁻⁴ However, isolated posterior stabilization is often associated with a loss of correction and implant failure. Therefore, the anterior approach for restoration of the anterior spine has become popular. Additionally, surgical restoration of the anterior column is important for normal spinal biomechanics.⁵⁻⁸

Up to now, tricortical iliac crest bone graft has been the gold standard for vertebral body replacement, but it is associated with high donor site morbidity.^{2,9-14} Additional problems such as pseudarthrosis, graft collapse with kyphotic deformity, and graft extrusion have led to a rapid increase in the use of cages for vertebral body replacement in the thoracolumbar spine.^{2,9,10,15,16} The meshed titanium cage, according to Lowery and Harms (DePuy AcroMed, Sulzbach, Germany), is one of the well-known cages used for vertebral body replacement.¹⁷ It is a nondistractable designed cage.

Recently, new expandable cages have been developed. In comparison to nonexpandable cages, these implants offer several surgical advantages. However, the biomechanical *in vitro* properties of these implants have not been characterized until now. Therefore, the aim of this study was to evaluate biomechanical properties of expandable cages for vertebral body replacement in the human thoracolumbar spine. Additionally, different in-

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Summary

This summary of 510(k) with the requirements of Incubator. being submitted in accordance

Company making the submission:

Name Address Telephone Fax Toll Free Contact E-mail		or	
---	--	----	--

1. **Device Name:**
 Trade/Proprietary Name:
 Common/Usual Name:
 Classification Name:
2. **Predicate Devices:**
 The _____ Incubator is substantially equivalent to other neonatal incubators on the market such as: _____ and Isolette Infant Incubator.
3. **Intended Use Statement:**
 Indications for Use: The _____ Intensive Care Incubator provides a controlled thermal environment for neonates who are unable to provide their own thermoregulation. The _____ Care Incubator can be used in two operating modes: Air Control Mode and Skin Control Mode with two probes that can be used with twins or for a follow up of central/peripheral temperature variations in a patient as a tool for the early detection of thermal stress-effect on the newborn.

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Includes: 4 syringes (1.25g/1.21ml ea)

RAVES & RANTS

- + Relatively nonrunny
- + Nice blue color
- Gets hot fast
- Only meant for lab use

Consistency and Handling

4.5 Only slightly runny.

Color

Blue.

Packaging

Cardboard box housing a clear plastic tray securing the syringes. The black plastic, push-type syringes have a small white cap on the end of tip that also pulls off. The expiration date is printed on a sticker on the back of the cardboard box.

Directions

Simple instructions printed on the bottom of the cardboard box. >

For Bleaching and Soft Tissue

1

Ultradent

★★★★
(4.3)

Cost: \$37.99 (\$7.91/ml)

Shelf life: 2 years refrigerated

Includes:

- 4 syringes (1.2ml ea)
- \$32.99/4 (\$6.87/ml)
- 10 Micro 20 Tips
- 10 Black Mini Tips

MSDS: Included

RAVES & RANTS

- + Light color good for bleaching
- + Viscosity makes applying on tissue easy
- Runny consistency
- Expensive

Consistency and Handling

3.5 Medium viscosity, but tends to be runny.

Color

Off-white with sparkles.

Exothermic Polymerization Heat

Produced a painful response in one test subject after three seconds of curing. The other subject endured the entire 20 seconds of the test, noting the temperature seemed to peak at four seconds and then regressed to a more or less tolerable level.

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TEST REPORT IEC 60601-1 Medical electrical equipment Part 1: General requirements for safety	
Report Reference No.	659100-2830-0008/58845 CB-Ref. No.: DE 1 -
Tested by (name + signature)	Dirk Ritter
Reviewed by (name + signature).....	Michael Mayer
Date of issue	
Contents	105 pages File: Testreport-FG12-8-58645_60601-1.doc
Testing Laboratory	VDE Testing and Certification Institute
Address	, Germany
Testing location/procedure	CB TL <input checked="" type="checkbox"/> RMT <input type="checkbox"/> SMT <input type="checkbox"/> WMT <input type="checkbox"/> TMP <input type="checkbox"/> Section FG12 / 110152
Applicant's name	
Address	
Manufacturer's name	
Address	
Factory's name	
Address	
Test specification:	
Standard	IEC 601 1:1988 + A1:1991 + A2:1995
Test procedure	CB Scheme
Non-standard test method	N/A
Test Report Form No.	601-1_C/
TRF Originator	Underwriters Laboratories Inc.
Master TRF	dated
<p>Copyright © 2005 IEC System for Conformity Testing and Certification of Electrical Equipment (IECEE), Geneva, Switzerland. All rights reserved.</p> <p><small>This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context.</small></p>	
Test item description	
Trade Mark	
Model/Type reference	
Ratings	100 – 240 V~ via external power supply unit (Power Brick) type

IEC60601-1 (ed.2)			
Clause	Requirement + Test	Result - Remark	Verdict
6	IDENTIFICATION, MARKING AND DOCUMENTS		—
6.1	Marking on the outside of equipment or equipment parts		P
6.1 c	Markings of the specific power supply affixed		P
6.1 d	If marking is not practicable due to size or nature of enclosure, information is included in accompanying documents		N/A
6.1 e	Name and/or trademark of the manufacturer or supplier		P
6.1 f	Model or type reference	l the frame ed on	P
6.1 g	Rated supply voltages or voltage range(s)	Ratings for 100 – 240 V~ 1.3 – 0.7 A 50/60 Hz ratings for Lithium Ion Battery, type: : 10,8V, 1Ah	P
	Number of phases	Single phase	P
	Type of current	Symbol ~	P
6.1 h	Rated frequency or rated frequency range(s) (Hz) ...	50/60 Hz	P
6.1 j	Rated power input (VA, W or A)	1.3 – 0.7 A	P
6.1 k	Power output of auxiliary mains socket-outlets		N/A
6.1 l	Class II symbol		P
	Symbol for degree of protection against ingress of water provided	IPX1 for hospital use (indoor) IP32 for "out of hospital" IP31 all specified positions and IP32 special mounting position with supply leads from the bottom	P
	Symbol for protection against electric shock		P
	If equipment has more than one applied part with different degrees of protection, the relevant symbols are clearly marked on such applied parts, or on or near relevant outlets		P
	Symbol for protection of defibrillation-proof applied parts		P
	Symbol 14 from Table D1 for defibrillation-proof with protection partly in patient cable		P
6.1 m	Mode of operation (if no marking, suitable for continuous operation)	Continuous mode	P
6.1 n	Types and rating of external accessible fuses	Primary fuses of PowerBrick not accessible by operator	N/A
6.1 p	Ratings of external output		N/A
6.1 q	Symbol for physiological effect(s):		-
	- attention, consult accompanying documents		N/A

IEC60601-1 (ed.2)					
Clause	Requirement + Test			Result - Remark	Verdict
56.1	TABLE: lists of critical component parts				
Object/part No			Technical data	Standard	Mark(s) of conformity ¹⁾
1.1 SMPS			Input: 100-240VAC 50/60 Hz 1.3-0.7 A Output: 24VDC,30Whom. 40W/10min/hr 50W 1min/hr Ambient Temp: 0-65°C	UL, CSA, VDE IEC 60950 and requirements according IEC 60601-1	VDE
1.2 Fuse			Current: 2A characteristic: time lag T breaking capacity: high	EN/IEC 60127-2	UL, CSA
Main - Board					
2.1 LAN transformer			1.5 KV prim to sec		VDE
2.2 C197, C198, C199, C211, C210			Class X1/Y1 250V AC 2200pF	IEC60384-14	UL, CSA
2.8 RT1 overload protection			2,2KΩ to 100 KΩ (25°C)	-	-
2.9 outer enclosure			>PC+ABS+FR< UL flame rating: 94V-0	UL	QMFZ2, QMFZ8.
ECG board					
optocoupler			-5000 V rms/ 1 min UL recognized highest over voltage 6000 V Case Temperature Ts = 175 Input Current Is Input = 230 Output Power Ps Output = 600	UL 1577 IEC/EN/DIN EN 60747-5-2 Approval DIN VDE 0884	FPQU2, File E 55361
transformer			2,27 mH Test voltage 5 KVrms	IEC 85	-
Outer enclosure				IEC 85	QMFZ2.E41871
Potting material			94 V 0	IEC 85	QMFZ2.E96722 QMFZ8.E96722

Korea - Differences to IEC 60601-1:1988 + A1:1991 + A2:1995			
6.1j	Insert the following sub-clause between the second and third sub-clauses: Equipment for one or several RATED voltage or frequency ranges, the RATED input for 220 V, 60 Hz or if applicable for 110 V, 60 Hz shall be separately marked.	The input ratings for 220V, 60 Hz and 110V, 60 HZ will be explicitly marked on the device if the product is to be delivered to the Korean market.	Pass
6.1s	HIGH VOLTAGE TERMINAL DEVICES on the outside of EQUIPMENT which are accessible without the use of a TOOL shall be marked with the symbol "dangerous voltage" (see Appendix D, Table DII, Symbol 6) and wit the Korean language, .		N/A
6.2c	Replace the existing subclause wit the following: The presence of HIGH VOLTAGE PARTS shall be marked with the symbol "dangerous voltage" (see Appendix D, Table DII, Symbol 6) and in the Korean language.		N/A
6.8.1	Insert the following sub-clause after the last paragraph: Language of accompanying documents shall be included in Korean.	The documents will be translated into Korean.	Pass
General	National supply voltages are 110 V, 220 V and 380 V.		Pass
General	Only appliances having supply frequency of 60 Hz or a frequency range including 60 Hz are accepted.		Pass
General	Instruction manuals and appliance markings related to safety, including nameplate, shall be in Korean or graphical symbols in accordance with IEC Publication 417.	Confirmation the instruction manual will be translated into the Korean language provided via declaration letter.	Pass
General	Plugs for connection of the equipment to the supply mains shall comply with the Korean Standard (KSC 8305 and 8300). More details are available from KTL (c/o KTL) on request.	Instructions are provided to allow selection of a plug that complies with Korean Standard (KSC 8305 and 8300).	Pass

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TEST REPORT		
Page: 01 DOC.NUM.: R60601100 Rev: 4.0 Date: 1/2000	IEC 60601-1,-1-1, -1-3, -1-4,-2-7, -2-28, 2-32 Medical electric equipment Part 1: General requirements for safety	
Report		
Reference No.....	12307250 001	
Compiled by (+ signature)	<input type="text"/>	
Approved by (+ signature)	<input type="text"/>	
Date of issue	2008-07-29	
Contents.....	122 pages	
<small>This report is based on a blank test report that was prepared by KEMA using information obtained from the TRF originator (see below).</small>		
Testing laboratory		
Name	TÜV Rheinland Japan Ltd., Global Technology Assessment Center	
Address	4-25-2 Kita-Yamata, Tsuzuki-ku, Yokohama 224-0021 Japan	
Testing location.....	<input type="text"/>	
Client		
Name	<input type="text"/>	
Address.....	<input type="text"/>	
Test specification		
Standard	IEC 60601-1:1988+A1:1999+A2:1995 IEC 60601-1-1:2000 IEC 60601-1-3:1994 IEC 60601-1-4:1996 + A1:1999 IEC 60601-2-7:1998 IEC 60601-2-28:1993 IEC 60601-2-32:1994	
Test procedure	CB scheme	
Procedure deviation.....	N.A.	
Non-standard test method.....	N.A.	
Test Report Form/blank test report		
Test Report Form No.....	1001-1-01/IEC60601-1-1-D/IEC60601-1-3B/IEC60601-1-4C /IEC60601-2-28B/IEC60601-2-32C/IEC60601-2-7C	
TRF originator.....	UL / UL / UL / UL / UL / UL / UL	
Master TRF.....	Dated 97-04 / 06-09 / 06-09 / 06-09 / 06-09 / 06-10 / 06-11	
<small>Copyright reserved to the bodies participating in the Committee of Certification Bodies (CCB) and/or the bodies participating in the CENELEC Certification Agreement (CCA).</small>		
Test item		
Description	X-RAY TV SYSTEM	
Trademark	<input type="text"/>	
Model and/or type reference	X-RAY TV SYSTEM SONIALVISION safire17	
Manufacturer.....	Same as applicant	

IEC 60601-1

Factory	:	<input type="text"/>
Input Rating(s)	:	3N 200/380/400/415/440/480VAC, 133kVA, 50/60Hz 1N 200/220/230/240VAC, 9.8kVA, 50/60Hz
Particulars: test item vs. test requirements		
Classification of installation and use	:	Permanently installed
Supply connection.....	:	Permanently installed
Test case verdicts		
Test case does not apply to the test object.....	:	N/A
Test item does meet the requirement	:	P(ass)
Test item does not meet the requirement	:	F(ail)
Testing		
Date of receipt of test item	:	N/A
Date(s) of performance of test.....	:	2008-06-27 – 2008-06-30
General remarks		
<p>This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.</p> <p>"This test report shall not be reproduced except in full without the written approval of the testing laboratory. The test results presented in this report relate only to the item tested. "(see remark #)" refers to a remark appended to the report. "(see appended table)" refers to a table appended to the report (same number of table as number of clause). Throughout this report a comma is used as the decimal separator.</p>		
Summary of compliance with National Differences:		
EN 60601-1:1990+A1:1993+A2:1995+A13:1996		
EN 60601-1-1:2001		
EN 60601-1-3:1994		
EN 60601-1-4: 1996 + A1:1999		
EN 60601-2-7:1998		
EN 60601-2-28:1993		
EN 60601-2-32:1994		
UL 60601-1:2003 R6.03		
CAN/CSA-C22.2 No. 601.1-M90		
JIS T 0601-1:1999		
JIS T 0601-1-1:2005		
Summary of compliance with National Differences (for explanation of codes see below): US, CA, JP CA=Canada, US=United States of America, JP=Japan		

IEC 60601-1-3				
Clause	Requirement	Test	Result – Remark	Verdict

4	General Requirements For Tests			P
4.201.1	Statement of compliance is expressed in accordance with this Standard		Statement of compliance is stated in Operation Manual. Refer to clause 6.8 of the ICE60601-1 part of this test report.	P


6	IDENTIFICATION, MARKING, AND DOCUMENTS			P
6.1.202	All sub-assemblies, components, and ACCESSORIES of X-RAY EQUIPMENT that can be removed are marked so that they can be identified and correlated with their ACCOMPANYING DOCUMENTS		All sub-assembly components are adequately marked for correlation with their respective accompanying documents	P
	- Interchangeable devices are distinguishable to the OPERATOR		No interchangeable devices	N/A
	Markings permanently affixed and clearly legible per Clause 6 of Part 1		See general standard IEC60601-1 part of this test report	P
6.8	Accompanying documents			P
6.8.201	ACCOMPANYING DOCUMENTS include the required statements per sub-clauses in Table 202		Operation Manual provided	P
6.8.202	ACCOMPANYING DOCUMENTS identify the items they refer to and include replication of information required to be marked on the items			P
	Information on the location and content of required markings on subassemblies is provided		Information provided	P
	a) For accessible markings, the location of the markings and instructions enabling them to be inspected			P
	b) For markings inaccessible in complete assembly			P
	1) Replication of the information marked on the sub-assemblies, or		Replicated in manual	P
	2) A list of the subassemblies with references to their accompanying documents			P
	ASSEMBLING INSTRUCTIONS for items supplied separately from main assemblies are provided		To be installed by the authorized service personnel. Installation manual provided	P

29	X-RADIATION			P
29.201.1	The operating range of X-RAY TUBE VOLTAGE for dental applications restricted per Table 203 (kV).....		Not intended for dental applications	N/A
29.201.2	First HALF-VALUE LAYERS attained were not less than the minimum values in Table 204 (mm Al)		2.78mm Al at 70kV (>2.1mm required)	P

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	<h1>EMC</h1>	<p>Ref. Certif. No. FI 6784</p>
<p>IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT (IECEE) CB EMC TEST CERTIFICATE</p>		
<p>Product</p> <p>Name and address of the applicant</p> <p>Name and address of the manufacturer</p> <p>Name and address of the factory</p> <p>Ratings and principal characteristics</p> <p>Trademark (if any)</p> <p>Model / Type Ref.</p> <p>Additional information (if necessary may also be reported on page 2)</p> <p>A sample of the product was tested and found to be in conformity with</p> <p>As shown in the Test Report Ref. No. which forms part of this Certificate</p> <p>This CB Test Certificate is issued by the National Certification Body of Finland</p>	<p>Electronic APEX Locator</p> <div style="border: 1px solid black; height: 40px; width: 100%;"></div> <p>Same as applicant</p> <p>Same as applicant</p> <p>Battery: 4.5Vd.c., (1.5Vd.c. x 3EA, type: LR3, size: AAA)</p> <div style="border: 1px solid black; height: 40px; width: 100%;"></div> <p>ROOTOR, i-ROOT 100 mini</p> <p>IEC 60601-1-2:2007 (Ed 3)</p> <p>F690501/RF-EMG003932</p>	
<p>SGS Fimko's ref.no. 264523-1</p> <p>Date: 21 June 2011</p>	<p>Signature: </p> <p>Ari Honkala Product Line Manager</p>	

Test Report issued under the responsibility of:

TEST REPORT IEC 60601-1-2 Medical Electrical Equipment PART 1-2: General Requirements for Basic Safety and Essential Performance Collateral Standard: Electromagnetic Compatibility	
Report Reference No.	F690501/RF-EMG003932
Date of issue	June 16, 2011
Total number of pages	
CB Testing Laboratory	SGS Korea Co.,Ltd.
Address	#18-34, Sanbon-dong, Gunpo, Gyeonggi-do, Korea
Applicant's name	
Address	
Test specification:	
Standard	IEC 60601-1-2 (ed 3)
Test procedure	CB
Non-standard test method	N/A
Test Report Form No	IEC60601_1_2BEMC
Test Report Form(s) Originator	Underwriters Laboratories Inc.
Master TRF	Dated 2008-06
Copyright © 2008 IEC System for Conformity Testing and Certification of Electrical Equipment (IECEE), Geneva, Switzerland. All rights reserved.	
This publication may be reproduced in whole or in part for non-commercial purposes as long as the IECEE is acknowledged as copyright owner and source of the material. IECEE takes no responsibility for and will not assume liability for damages resulting from the reader's interpretation of the reproduced material due to its placement and context.	
If this Test Report Form is used by non-IECEE members, the IECEE/IEC logo and the reference to the CB Scheme procedure shall be removed.	
This report is not valid as a CB Test Report unless signed by an approved CB Testing Laboratory and appended to a CB Test Certificate issued by an NCB in accordance with IECEE 02.	
Test item description	Electronic APEX Locator
Trade Mark	
Manufacturer	
Model/Type reference	ROOTOR (added Mode : i-ROOT 100 mini)
Ratings	Battery : 4,5 Vd.c.(1,5 Vd.c. x 3 EA, type:LR3, size:AAA)

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1. PURPOSE / OBJECTIVE

The purpose of this report is to explain why the results from [redacted] can be applied to [redacted] [redacted] / stent [redacted] and therefore demonstrate that the [redacted] is non-cytotoxic.

2. DESCRIPTION OF THE ORIGINAL TEST

Test article: [redacted] stent material (silicon carbide coated CoCr stent material, fully processed)

A single extract of the test article was prepared using MEM supplemented with 5% serum and 2% antibiotics based on the ration of 120 cm² : 20 mL. This test extract was placed onto three separate monolayers of L-929 mouse fibroblast cells. Three separate monolayers were prepared for the reagent control, negative control and for the positive control. All monolayers were incubated at 37°C for 48 hours. The monolayer in the test article, reagent control, negative control and positive control wells was examined microscopically at 48 hours to determine any change in cell morphology. The biological reactivity (cytotoxic effect) was graded on a scale of 0 to 4 (none, slight, mild, moderate or severe). Please refer to attachment 1 for details.

3. JUSTIFICATION FOR TRANSFER OF RESULTS

The cobalt chromium stent body raw material for [redacted] and [redacted] stents is the same material sourced from the same suppliers to the same specifications (except dimensions). The stents undergo the same manufacturing processes from raw material to finished product. Therefore the test article (silicon carbide coated CoCr stent material, fully processed) is representative of both stents in their finished state.

Therefore, the results of the biocompatibility testing of the [redacted] stent are also valid for the [redacted] stent.

4. CONCLUSION

The MEM test article extract showed no evidence of causing cell lysis or toxicity (grade 0). The MEM test article extract met the requirements of the test since the grade was less than a grade 2 (mild reactivity). The reagent control, negative control and the positive control performed as anticipated.

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Author

Final Report Date
September 1, 2006

COMPLIANCE
21 CFR, Part 58
Good Laboratory Practice for Non-Clinical Laboratory Studies

MANAGEMENT OF THE STUDY

Performing Laboratory
Toxikon Corporation
15 Wiggins Avenue
Bedford, MA 01730

Sponsor

TABLE OF CONTENTS

Title Page
Table of Contents
Study Summary
Quality Assurance Statement
Study Director Signature and Verification Dates

1.0 Purpose
2.0 References
3.0 Compliance
4.0 Identification of Test and Control Articles
5.0 Identification of Test System
6.0 Justification of Test System and Route of Administration
7.0 Experimental Design and Dosage
8.0 Evaluation Criteria
9.0 Results
10.0 Conclusion
11.0 Records
12.0 Confidentiality Agreement

STUDY SUMMARY

No biological reactivity (Grade 0) was observed in the L929 mammalian cells at 48 hours, post exposure to the test article extract. The observed cellular response obtained from the positive control article extract (Grade 4) and negative control article extract (Grade 0) confirmed the suitability of the test system. Based on the criteria of the protocol, the test article, _____, is considered non-cytotoxic and meets the requirements of the Elution Test, ISO 10993-5 guidelines.


QUALITY ASSURANCE STATEMENT

This study was conducted in compliance with U.S. Food and Drug Administration regulations set forth in 21 CFR, Part 58.

The sections of the regulations not performed by or under the direction of Toxikon Corporation, exempt from this Good Laboratory Practice Statement, included characterization and stability of the test article and its mixture with carriers, 21 CFR, Parts 58.105 and 58.113.

The Quality Assurance Unit conducted inspections on the following dates. The findings were reported to the Study Director and to Toxikon's Management.

INSPECTIONS	DATE OF INSPECTION	DATE REPORTED STUDY DIRECTOR	DATE REPORTED MANAGEMENT
DOSE ADMINISTRATION	08/25/06	08/25/06	08/25/06
RAW DATA	09/01/06	09/01/06	09/01/06
FINAL REPORT	09/01/06	09/01/06	09/01/06

Allison Lyons-Hook, B.A.
Quality Assurance

Date

STUDY DIRECTOR SIGNATURE AND VERIFICATION DATES

This study meets the technical requirements of the protocol. The study also meets with the requirements of the Good Laboratory Practice Regulations, 21 CFR, Part 58, with the exemptions as stated in the Quality Assurance Statement.

Protocol Number: AMD/VITRO/003-06/000

Study Director: Patricia S. Grutkoski, Ph.D.

Company: Toxikon Corporation

Signature:

Pat S. Grutkoski

Date:

9/1/06

Study Supervisor:

Patricia S. Grutkoski, Ph.D.

VERIFICATION DATES:

Protocol Effective Date: 08/11/06
Test Article Receipt: 08/22/06
Project Log Date: 08/22/06
Extraction Dates: 08/24/06 – 08/25/06
Technical Initiation: 08/25/06
Technical Completion: 08/27/06

1.0 PURPOSE

The purpose of the study was to determine the biological reactivity of a mammalian cell culture (L929) in response to the test article. The study was designed for the evaluation of test article extracts (e.g., polymeric materials, high-density materials, etc.).

2.0 REFERENCES

The study was conducted based upon the following references:

- 2.1 ISO 10993-5, 1999, Biological Evaluation of Medical Devices – Part 5: Tests for *In Vitro* Cytotoxicity.
- 2.2 ISO 10993-12, 2002, Biological Evaluation of Medical Devices – Part 12: Sample Preparation and Reference Materials.
- 2.3 ISO/IEC 17025, 2005, General Requirements for the Competence of Testing and Calibration Laboratories.

3.0 COMPLIANCE

The study conformed to the current FDA 21 CFR, Part 58 – Good Laboratory Practice for Non-Clinical Laboratory Studies.

4.0 IDENTIFICATION OF TEST AND CONTROL ARTICLES

The Sponsor supplied the following information on a test requisition form or other correspondence, wherever applicable (excluding confidential or trade secret information). The Sponsor was responsible for all test article characterization data as specified in the GLP regulations.

4.1 Test Article:

Test Article Name:

CAS/Code #: Not Supplied by Sponsor (N/S)

Lot/Batch #:

Physical State: N/S

Color: N/S

Expiration Date: N/S

Density: N/S

Stability: N/S

Solubility: N/S

pH: N/S

Storage Conditions: Room Temperature

Safety Precautions: Standard Toxikon Laboratory Safety Precautions

4.2 Control Articles (Toxikon Supplied):

4.2.1 Negative Control Article Name: Negative Control High Density Polyethylene (Negative Control Plastic)

Toxikon QC #: CSC-04-05-009-CC

Physical State: Solid

Color: White

Storage Conditions: Room Temperature

Safety Precautions: Standard Laboratory Safety Precaution

4.2.2 Positive Control Article Name: Natural Rubber

Toxikon QC #: CSC-05-03-005-CC

Physical State: Solid

Color: Amber

Stability: Stable at Room Temperature

Storage Conditions: Room Temperature

Safety Precautions: Standard Laboratory Safety Precautions

5.0 IDENTIFICATION OF TEST SYSTEM

The test system was mouse fibroblast L929 cells. The cell line was obtained from the American Type Culture Collection, Manassas, Virginia.

6.0 JUSTIFICATION OF TEST SYSTEM AND ROUTE OF ADMINISTRATION

6.1 Historically, mouse fibroblast L929 cells had been used for cytotoxicity studies because they demonstrate sensitivity to extractable cytotoxic articles.

6.2 The test article was extracted and administered *in vitro* to mouse fibroblast L929 cells through a solvent compatible with the test system. This was the optimal route of administration available in this test system as recommended in the ISO 10993-5 guidelines.

7.0 EXPERIMENTAL DESIGN AND DOSAGE

7.1 Preparation of Test and Control Articles:

7.1.1 The extraction vehicle was Serum-Supplemented (Complete) Minimum Essential Medium (MEM) supplemented with 10% fetal bovine serum. Per Sponsor's request, representative samples of each component of the test article were pooled and extracted. The test article (2.0 g) was extracted in 10 mL of medium (0.2 g per 1.0 mL ratio) at 37 ± 1 °C for 24 ± 2 hours.

7.1.2 Extracts prepared with complete MEM were tested at 100% (neat) concentration.

7.1.3 The positive control article (Natural Rubber) and negative control article (Negative Control Plastic) were similarly prepared and extracted.

7.2 Pre-Dose Procedure:**7.2.1 Cell Culture Preparation:**

Cultures of L929 cells were initiated in complete MEM not less than 24 hours prior to testing and incubated at 37 ± 1 °C, in a humidified atmosphere containing $5 \pm 1\%$ carbon dioxide.

7.2.2 Extracts were evaluated for a change in pH. No adjustment was necessary since the pH of the extract was within ± 0.2 of the control medium. The extracts were filter sterilized by passage through a $0.22 \mu\text{m}$ pore filter, prior to being applied to the cell monolayer.

7.3 Dose Administration:

A volume of 3 mL of the test article or control article extract was used to replace the maintenance medium of the cell culture. All extracts were tested in triplicate.

7.4 Post-Dose Procedure:

All cultures were incubated for 48 ± 2 hours, at 37 ± 1 °C, in a humidified atmosphere containing $5 \pm 1\%$ carbon dioxide.

8.0 EVALUATION CRITERIA

8.1 The response of the cell monolayer is evaluated under a microscope. Cytochemical stains are used in the evaluation. The biological reactivity (cellular degeneration and malformation) is rated on a scale of 0 to 4. The test system is considered suitable if no signs of cellular reactivity (Grade 0) are noted for the negative control article and the positive control article showed greater than a Mild reactivity (Grade 2). The test article meets the requirements of the test if none of the cultures treated with the test article show greater than a Mild reactivity (Grade 2). If the suitability of the system is not confirmed, the test is repeated.

Grade	Reactivity	Description of Reactivity Zone
0	None	Discrete intracytoplasmic granules; no cell lysis
1	Slight	Not more than 20% of the cells are round, loosely attached, and without intracytoplasmic granules; occasional lysed cells are present
2	Mild	Not more than 50% of the cells are round and devoid of intracytoplasmic granules; no extensive cell lysis and empty areas between cells
3	Moderate	Not more than 70% of the cell layers contain rounded cells or are lysed
4	Severe	Nearly complete destruction of the cell layers

8.2 The study and its design employ methodology to minimize uncertainty of measurement and control of bias for data collection and analysis.

9.0 RESULTS

At the 48 hour observations, no signs of cellular reactivity (Grade 0) were observed for the negative control article and the test article extracts. Severe signs of reactivity (Grade 4) were noted for the positive control article extract at the 48 hour observation.

REACTIVITY GRADES

Time	Date	Test Article			Controls								
					Medium			Negative			Positive		
		A	B	C	A	B	C	A	B	C	A	B	C
0 Hours	08/25/06	0	0	0	0	0	0	0	0	0	0	0	0
24 Hours	08/26/06	0	0	0	0	0	0	0	0	0	3	3	3
48 Hours	08/27/06	0	0	0	0	0	0	0	0	0	4	4	4

10.0 CONCLUSION

No biological reactivity (Grade 0) was observed in the L929 mammalian cells at 48 hours, post exposure to the test article extract. The observed cellular response obtained from the positive control article extract (Grade 4) and negative control article extract (Grade 0) confirmed the suitability of the test system. Based on the criteria of the protocol, the test article, _____ is considered non-cytotoxic and meets the requirements of the Elution Test, ISO 10993-5 guidelines.

11.0 RECORDS

- 11.1 Original raw data are archived at Toxikon Corporation.
- 11.2 A copy of the final report and any report amendments is archived at Toxikon Corporation.
- 11.3 The original final report, and a copy of any protocol amendments or deviations, is forwarded to the Sponsor.
- 11.4 All unused test article shall be disposed of by Toxikon, per Sponsor's request.

12.0 CONFIDENTIALITY AGREEMENT

Statements of confidentiality were as agreed upon prior to study.

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	Document form: Software Development Version: <input style="width: 50px;" type="text"/>	
<h2 style="margin: 0;">Software Requirements Specification</h2> <div style="text-align: center; margin-top: 5px;"> <input style="width: 200px; height: 15px;" type="text"/> <input style="width: 100px; height: 15px;" type="text"/> </div>		
Project business name: <input style="width: 100%; height: 20px;" type="text"/>	Internal project name: <input style="width: 50%; height: 20px;" type="text"/>	
Made by: Position: <input style="width: 150px;" type="text"/> Name: <input style="width: 300px;" type="text"/> Date: <input style="width: 80px;" type="text"/> Signature: <input style="width: 300px; height: 20px;" type="text"/>		
Approved by: Position: <input style="width: 60px;" type="text"/> Name: <input style="width: 150px;" type="text"/> Date: <input style="width: 60px;" type="text"/> Signature: <input style="width: 150px; height: 20px;" type="text"/>	Approved by: Position: <input style="width: 60px;" type="text"/> Name: <input style="width: 150px;" type="text"/> Date: <input style="width: 60px;" type="text"/> Signature: <input style="width: 150px; height: 20px;" type="text"/>	Approved by: Position: <input style="width: 60px;" type="text"/> Name: <input style="width: 150px;" type="text"/> Date: <input style="width: 60px;" type="text"/> Signature: <input style="width: 150px; height: 20px;" type="text"/>
Valid from: <input style="width: 100px;" type="text"/>		
Document revision:		

RECORD OF CHANGES						
Revision	Date	Last version			Change Description	Approved by
		Version	Valid from	Valid to		

1. Table of contents:

- 1. General context..... 2
 - 1.1. User Interface 3
 - 1.2. Remote Control/..... 5
 - 1.3. Recording..... 6
 - 1.3.1. converter data reading – 6
 - 1.3.2. converter data reading – 7
 - 1.3.3. Runtime checks and warnings..... 8
 - 1.3.4. File system and 8

2. General context

The software/[] shall be implemented as an [] into a specific hardware environment. The environment is specified in Tab 1.1. The software shall have [] as described on state diagram in Fig. 1.1.

1. User interface - setting up and starting a recording, []
2. Recording – recording of the ecg signal to the []
3. Remote control/online transmission – communication with [] connection


Fig.1.1: State diagram – basic software states and []

Tab 1.1: Hardware specification

Component	Specification
CPU	[]
File system	[]
RAM	[]
Display	[]
Storage	[]
Peripheral Port	[]
Wireless Connection	[]

1.1. User Interface

All the user interface shall be implemented as a set of screens and []. Because of the limited operation memory size only one common [] and one modal [] can be allocated at a time. [] shall be connected to a message loop from buttons and shall be routed from application [] of the particular []. [] is responsible for creating a new one before [].

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Document form: Software Development
Version:

Software Requirements Specification

Project business name: **Internal project name:**

Made by:
Position:
Name:
Date:
Signature:

<p>Approved by: Position: <input type="text"/> Name: <input type="text"/> Date: <input type="text"/> Signature: <input type="text"/></p>	<p>Approved by: Position: <input type="text"/> Name: <input type="text"/> Date: <input type="text"/> Signature: <input type="text"/></p>	<p>Approved by: Position: <input type="text"/> Name: <input type="text"/> Date: <input type="text"/> Signature: <input type="text"/></p>
---	---	---

Valid from:

Document revision:

RECORD OF CHANGES						
Revision	Date	Last version			Change Description	Approved by
		Version	Valid from	Valid to		

CONTENT

- 1 FUNCTIONAL REQUIREMENTS 3
 - 1.1 Performance 3
 - 1.1.1 Purpose of the software 3
 - 1.1.2 [] requirements 3
 - 1.2 Physical characteristics 3
 - 1.3 Computing environment under which the software is to perform 3
- 2 SOFTWARE SYSTEM INPUTS AND OUTPUTS 4
 - 2.1 Inputs: 4
 - 2.2 Outputs: 4
- 3 INTERFACES 5
 - 3.1 [] format 5
 - 3.1.1 [] format requirements 5
 - 3.1.2 Main format characteristics 5
 - 3.1.3 [] Types 5
 - 3.2 Online signal transmission 7
 - 3.2.1 Online transmission protocol requirements 8
 - 3.2.2 Main format characteristics 8
 - 3.2.3 Recording settings via [] 9
- 4 Graphical user interface 9
 - 4.1 Starting a new recording 9
 - 4.2 Continue a [] recording 9
 - 4.3 Recording settings 9
 - 4.4 Unit configuration 9
 - 4.5 Information 10
- 5 SOFTWARE-DRIVEN ALERTS 10
 - 5.1 [] 10
 - 5.2 [] 10

1 FUNCTIONAL REQUIREMENTS

1.1 Performance

1.1.1 Purpose of the software

The [Redacted] software is intended to be used by trained operators under direct supervision of a licensed health care practitioner on adult and pediatric patients. The software shall consist of [Redacted]. The first part drives the hardware components of the unit to acquire and store the [Redacted] measured on its input and the second part provides the operator with a user interface.

The user interface part of the software shall allow the user to setup the recording, start the recording, setup the unit, display information about the unit. Also the software shall signal when [Redacted] get loose or fallen off a patient and allows short voice recordings on user's demand.

Data acquisition functions of the software are sampling and reading [Redacted], [Redacted] detection, detection of [Redacted] and saving all data to a [Redacted].

The software does not perform any medical related data analysis.

1.1.2 [Redacted] requirements

Software [Redacted] shall use unit's internal [Redacted] for all actions especially for signal sampling. Standard sampling frequency shall be [Redacted].

1.2 Physical characteristics

The software shall be implemented as an embedded system in [Redacted] language. All the functionality shall be implemented specifically for the hardware provided which is specified in the next chapter.

1.3 Computing environment under which the software is to perform

Component	Specification
CPU	[Redacted]
File system	[Redacted]
RAM	[Redacted]
Display	[Redacted]
Storage	[Redacted]
Peripheral Port	[Redacted]
Wireless Connection	[Redacted]

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SOFTWARE VALIDATION REPORT			
Project business name: <input style="width: 80%;" type="text"/>		Internal project name: <input style="width: 80%;" type="text"/>	
Document Title: Software Validation Report		Document revision: <input style="width: 50px;" type="text"/>	Date: <input style="width: 100px;" type="text"/>
Approval			
Title	Name	Signature	Date
Project manager	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>
Head of Development Department	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>
<input style="width: 50px;" type="text"/> manager	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>

Record of Changes

Rev.	Date	Author	Comments

Software Validation Report

CONTENT

1	Validation of <input style="width: 100px;" type="text"/>	3
2	FUNCTIONAL REQUIREMENTS	3
2.1	Performance PASSED!	3
2.1.1	Purpose of the software	3
2.1.2	<input style="width: 50px;" type="text"/> requirements PASSED!	3
2.1.3	Physical characteristics PASSED!	3
2.1.4	Computing environment under which the software is to perform PASSED!	4
2.2	Software system inputs and outputs PASSED!	4
2.2.1	Inputs: PASSED!	4
2.2.2	Outputs: PASSED!	4
2.3	Interfaces PASSED!	4
2.3.1	<input style="width: 50px;" type="text"/> PASSED!	4
2.3.2	<input style="width: 50px;" type="text"/> transmission PASSED!	6
2.3.3	Graphical user interface PASSED!	8
2.3.4	Software-driven alerts PASSED!	9
3	Test report	9
4	Risk analysis	9
5	Verification of <input style="width: 100px;" type="text"/>	10
6	Result of validation	10

1 Validation of SW []

Software validation refers to establishing, by objective evidence, that the software conforms to the user needs and intended uses of the device. Software validation is a part of design validation of the finished device. It involves checking for proper operation of the software in its actual or simulated use environment, including integration into the final device where appropriate. Software validation is highly dependent upon comprehensive software testing and other verification tasks previously completed at each stage of the software development life cycle. Planning, verification, configuration management and many other aspects of good software engineering are important activities that together help to support the conclusion that the software is validated.

Validation is done according to Std 1012-1998 IEEE Standard for Software *Verification and Validation*, *IEEE Standards Collection: Software Engineering*. [];

Document used for validation of []:

- Software requirements specification [] - definition of functionality of system []
- Test report version [] according Test protocol []
- Risk analysis [] during [] development
- Verifications of SW []

Validation is done according to internal directive [] version []

Risk analysis is done according to [] - Risk Analysis version []

The following paragraphs refer to Software requirements specification. All points described there were checked during validation.

2 FUNCTIONAL REQUIREMENTS

2.1 Performance

PASSED!

2.1.1 Purpose of the software

The [] unit software is intended to be used by [] operators under direct supervision of a licensed health care practitioner on adult and pediatric patients. The software shall consist of two parts. The first part drives the hardware components of the unit to acquire and store the [] measured on its input and the second part provides the operator with a user interface.

The user interface part of the software shall allow the user to setup the recording, start the recording, setup the unit, display information about the unit. Also the software shall signal when []

Data acquisition functions of the software are sampling and reading [] detection, detection of artificial cardiac []

The software does not perform any medical related data analysis.

2.1.2 [] requirements

PASSED!

Software [] shall use unit's internal [] for all actions especially for signal sampling. Standard sampling frequency shall be []

2.1.3 Physical characteristics

PASSED!

The software shall be implemented as an embedded system in [] language. All the functionality shall be implemented specifically for the hardware provided which is specified in the next chapter.

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Software Validation Report

Software Validation Report

Software Name : (Ver)
Ref No. :

	Drawing	Check by	Approval
Approval	<input type="text"/>	<input type="text"/>	<input type="text"/>

Corp.

INDEX

1.	Level of Concern	7
2.	Software Description	8
2.1.	Software Name : <input style="width: 100px; height: 15px;" type="text"/>	9
2.1.1.	Software Function	9
2.1.2.	<input style="width: 100px; height: 15px;" type="text"/>	10
2.2.	Hardware Requirement	10
2.2.1.	System requirement	10
2.2.2.	<input style="width: 100px; height: 15px;" type="text"/>	10
2.3.	Development Language and Program Size	11
2.3.1.	Program Language	11
2.3.2.	<input style="width: 100px; height: 15px;" type="text"/>	11
2.3.3.	Memory Management	11
2.3.4.	<input style="width: 100px; height: 15px;" type="text"/>	11
3.	Software Requirements Specification	12
3.1.	Interface	12
3.1.1.	User Interface	12
3.1.1.1.	<input style="width: 100px; height: 15px;" type="text"/>	12
3.1.2.	Software Interface	12
3.1.3.	System Interface	12
3.1.4.	<input style="width: 100px; height: 15px;" type="text"/> User Interface	13
3.2.	Functional Requirement	20
3.2.1.	Function Essential Factor by Module	20
3.3.	Description of Process	22
3.3.1.	Scanner Module	22
3.3.2.	<input style="width: 100px; height: 15px;" type="text"/> Module	23
3.3.3.	Database Module	23
3.3.4.	<input style="width: 100px; height: 15px;" type="text"/> Module	23
3.3.5.	Image Processing Module	24
3.3.6.	Viewer Module	24
3.3.7.	<input style="width: 100px; height: 15px;" type="text"/> Module	24
3.3.8.	Study List Module	25
3.3.9.	Transfer Information List Module	25
3.3.10.	<input style="width: 100px; height: 15px;" type="text"/>	25
3.3.11.	CD / DVD Write Module	25
3.3.12.	Patient Registration and Editing Module	26
3.4.	Limitation	26
3.5.	Recovery	26
3.6.	Inside checking	27
3.7.	Error and Interrupt	27
3.8.	Software Error	27
3.9.	Timing	27
4.	Application Algorithm and Technology	28
4.1.	Calibration	28

4.1.1.	Geometry Correction	28
4.1.2.	Offset Correction	28
4.1.3.	[REDACTED]	28
4.1.4.	Defect Correction.....	28
4.2.	Post-processing.....	28
4.2.1.	[REDACTED]	29
4.2.2.	Contrast and Edge Enhancement.....	29
4.2.3.	Dynamic Range Reduction	29
4.2.4.	[REDACTED]	30
4.2.5.	Tone Scaling	30
5.	[REDACTED] Design Chart	32
5.1.	Software Diagram and Description	32
5.1.1.	Description.....	32
5.1.2.	Diagram	33
5.2.	Flow Chart.....	34
5.2.1.	Initialization.....	35
5.2.2.	[REDACTED]Registration	35
5.2.3.	Image Grabbing.....	35
5.2.4.	[REDACTED]Search	35
5.2.5.	Data Selection and Load	35
5.2.6.	Image Display.....	35
5.2.7.	[REDACTED]	35
5.2.8.	Save and Transmission	36
5.2.9.	Finalization	36
5.3.	[REDACTED]	36
5.3.1.	Data Class Relationship Diagram.....	36
5.3.2.	Interface Class relationship diagram	36
5.4.	Database Description	37
5.4.1.	[REDACTED]	37
5.4.2.	Parameter Database	37
6.	Software Design Specification	38
6.1.	Total Process [REDACTED]	38
6.1.1.	Create Study from Work-List	38
6.1.2.	Patient Registration Process	38
6.1.3.	[REDACTED] Registration	39
6.1.4.	Image Acquisition from Scanner and [REDACTED]	39
6.1.5.	Image Processing.....	39
6.1.6.	[REDACTED]	39
6.1.7.	Image Transfer and Print.....	39
6.1.8.	Image Storing in CD / DVD.....	40
6.1.9.	Open Current Study by Image Viewer	40
6.1.10.	Close Study	40
6.1.11.	Verifying Transfer or Print Status in Transfer Status Window. ...	40
6.2.	Software component.....	40


6.3.1.	Scanner Module	42
6.3.2.	Calibration	51
6.3.3.	Database Module	55
6.3.4.	[REDACTED]	62
6.3.5.	Image Processing Module	84
6.3.6.	Viewer Module	90
6.3.7.	Work List Module	97
6.3.8.	[REDACTED] Module	99
6.3.9.	Transfer Information List Module	101
6.3.10.	Print Module	103
6.3.11.	[REDACTED]	105
6.3.12.	Patient Registration and Editing Module	107
7.	Software Development Environment Description	109
7.1.	Introduction	109
7.2.	Concepts	109
7.2.1.	Responsibilities	109
7.2.1.1.	Management	109
7.3.	Software Life Cycle	110
7.4.	Verification	111
7.4.1.	[REDACTED]	111
7.4.2.	System Verification	111
7.4.3.	[REDACTED]	111
7.5.	Validation	112
7.6.	Amendments	112
7.7.	Documentation / Main Product File / Risk Management File	112
8.	Software Verification and Validation Test	114
8.1.	Verification Test	114
8.1.1.	[REDACTED]	115
8.1.2.	Calibration Module	121
8.1.3.	Database Module	124
8.1.4.	[REDACTED]	132
8.1.5.	Image Processing Module	170
8.1.6.	Viewer Module	176
8.1.7.	[REDACTED]	187
8.1.8.	Study List Module	190
8.1.9.	Transfer Information List Module	193
8.1.10.	Print Module	195
8.1.11.	CD / DVD Write Module	197
8.1.12.	Patient Registration and Editing Module	200
8.2.	Validation Test	202
8.2.1.	Scanner Module	203
8.2.2.	Calibration Module	209
8.2.3.	[REDACTED]	211
8.2.4.	[REDACTED]	217

8.2.8. Study List Module	257
8.2.9. Transfer Information List Module	259
8.2.10. Print Module	261
8.2.11. CD / DVD Write Module.....	263
8.2.12. Patient Registration and Editing Module.....	265
9. Traceability Analysis	266
10. Hazard Analysis.....	300
10.1. Severity	300
10.2. <input type="text"/>	301
10.3. Risk Acceptability.....	302
10.4. Hazard Analysis.....	303
11. <input type="text"/>	317
12. <input type="text"/>	318
12.1. Scope	318
12.2. <input type="text"/>	318

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6. Software Design Specification

6.1. Total Process Flowchart

6.1.1. Create Study from Work-List

- ✓ Retrieve work list from Work-list server by query.
- ✓ Display retrieved work list in list control.
- ✓
- ✓ Start search by input data in search field and press enter.
- ✓ Selected study is started when the study is double clicked by left mouse button or click
- ✓ In the list control,
- ✓

6.1.2. Patient Registration Process

- ✓ New patient can be added through new patient registration process.
- ✓
- ✓ Date of Birth: Input in a form of
- ✓ Display error message and let user re-input when there are errors in compulsory input field or date of birth.

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8. Software Verification and Validation Test

8.1. Verification Test

✓ Test Categories

- Scanner Module
- Database Module
-
- Image Processing Module
- Viewer Module
-
- Study List Module
- Transfer Information List Module
-
-
- Patient Registration and

✓ Procedures

- Check to ensure all components are correctly connected.
- Indicate which components have been used; indicate manufacturer and trade / model name of product, version number of software etc.
- Execute the tests in the table below.
- Record the output for each test in the column.
-
- After executing all tests, compare the with the and then mark either Pass (P) or Fail (F) in the
- If any additional comments need to be made, note the comments in the
- Attach any related notes, references, test data used, and all printed reports with the test report, and indicate the related Test Number on each reference.
- Upon completion of the tests, the executor should sign under
-

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<input type="text"/> Corporation		<input type="text"/>	
Document title Specification Verification Test Procedure - <input type="text"/>	Document No. T - <input type="text"/>	Revision 6	Page 1 of 16

AUTHORIZATION			
	NAME	SIGNATURE	DATE
Issued by:	<input type="text"/>	<input type="text"/>	<input type="text"/>
Approved by:	<input type="text"/>	<input type="text"/>	<input type="text"/>
Approved by:	<input type="text"/>	<input type="text"/>	<input type="text"/>

REVISION HISTORY		
REVISION #	DESCRIPTION	DATE
1	New Issue	<input type="text"/>
2	Change of tidal volume and rate based on new specifications – revision 7	12/19/ <input type="text"/>
3	Change of rate, pressure and addition of the disconnect alarm	9/20/ <input type="text"/>
4	Sensitivity test changed to match product specifications. Units of measure added to all tests. All CPAP test criteria are changed to ± 2 cm H ₂ O. Pass / fail criteria added for disconnect alarm time delay. Low gas supply alarm restoration test added	12/9/ <input type="text"/>
5	Modify test procedure for calibrated rate control	2/7/ <input type="text"/>
6	Modify test procedure for Model S with calibrated rate control	1/11/ <input type="text"/>

<input type="text"/> Corporation		<input type="text"/>	
Document title Specification Verification Test Procedure — <input type="text"/>	Document No. T - <input type="text"/>	Revision 6	Page 2 of 16

1.0 Purpose

To test the Ventilator with the new enclosure and rate control calibration scale. These tests confirm that the device meets the requirements of the product specification.

2.0 Scope

ventilator testing

3.0 Definitions

- 3.1 Accuracy – How close the output is to desired measured by a referenced measurement device
- 3.2 Precision – Repeatability of a control
- 3.3 Specificity – Effect of one control on another
- 3.4 Sensitivity – How close a control can be set
- 3.5 Linearity – Equal incremental change in a control

4.0 Applicable Documents

- 4.1 Rev 8 Ventilator Specification

5.0 Procedure

- 5.1. Test conditions to be room temperature, humidity, barometric pressure
- 5.2. The Pneuton is to be set up as follows:
 - 5.2.1. a standard disposable patient circuit with bacteria filter
 - 5.2.2. compressed oxygen supplied from a cylinder. The ability to vary the input pressure from 41 to 66 psi is required.
 - 5.2.3. the patient circuit is attached to the ventilator and a test lung. The is set for a compliance of 50 ml/cm H₂O and resistance of 20 cm H₂O/L/sec.
 - 5.2.4. a calibrated ventilator tester (FlowAnalyzer) connected to a PC) is placed between the ventilator and the test lung. This device is used as an external measuring instrument for inspiratory pressure / flow and the generation of graphic pressure / flow traces.

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Certificate of Registration

The following organization's quality management system has been assessed and registered by Intertek Testing Services NA Ltd., a **CMDCAS recognized registrar**, as conforming to the requirements of:

ISO 13485:2003

Organization:

Corporation

Main Site: , USA

Certificate Number

Initial Certification Date
March 1, 2007

Certificate Effective Date
March 1, 2010

Certificate Expiry Date
March 1, 2013

The Quality Management System is applicable to:

Design, manufacture, distribution and service of Ventilators, CPAP Systems, breathing circuits and other respiratory support accessories.

In the issuance of this certificate, Intertek assumes no liability to any party other than to the Client, and then only in accordance with the agreed upon Certification Agreement.

 David S. Muil – General Manager - Canada
 Intertek Testing Services NA Ltd. – Lachine, QC, Canada

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(Manual 3)

Instruction for Use

Manufacturer has been granted certificates of ISO9001,ISO13485

Please read the Instruction for Use carefully before using the product.

Product Name: Hollow Fiber Dialyzer

Model: 

Product reference No.: 

Indications:  Dialyzer can be used for the hemodialysis treatment of acute and chronic renal failure.

Contraindications: No absolute contraindication for hemodialysis treatment. Strict monitoring shall be achieved for the patients who have tendencies of blood bleeding or cruor during the treatment. In case of any complication that would affect the stable condition, the treatment shall be terminated.

Warning:

- Air shall be removed from dialysate.
- The dialysate and the blood should flow in counter-current.
- The tolerance of blood compartment volume is $\pm 10\%$.
- The maximum fixed flow rate of dialysate is 800 mL/min.
- The effective membrane area shall not be less than 95% of the regulated value.
- During the treatment, the transmembrane pressure shall not exceed 500 mmHg.
- The blood flow rate shall not be less than 150 mL/min, but not more than 500 mL/min.
- The preparation for hemodialysis and treatment should be carried out under aseptic condition to avoid renal infection.
-  Dialyzer shall be used under certain medical supervision. To avoid bacterial and pyrogenic contamination, it is suggested to use together with hemodialysis machine and dialysate which are in accordance with the national and international standards, and most importantly, with the water, concentrated solution and dialysate which are in accordance with the national and international standards.
- This dialyzer is sterile, nonpyrogenic and nontoxic, sterilized by Gamma ray, with a shelf life of 3 years, which shall be checked before use. To avoid infection, do not use the expired product. Do not use the product if the package is damaged.
- This dialyzer is for single use and reuse is strictly prohibited. Reprocessing of this product may

immediately corrected by clamping off the tubing and replacing the dialyzer. Appropriate amount of blood shall be provided for blood loss if necessary.

- In case of blood leakage from the inner membrane of dialyzer, proper measures shall be immediately taken to stop the treatment and replace the dialyzer.

Note: The blood which might have been polluted by the dialysate is prohibited to return to the patient. Appropriate amount of blood shall be provided for blood loss if necessary.

3. Air Infiltration

- During the treatment, the dialyzer shall be replaced if there is air inside.
- To prevent any solution with air from entering the patient, it is necessary to use the dialysis equipment with air monitoring function according to the current regulations.

4. Hypersensitive Reaction

- If hypersensitive reaction occurs to the patient, proper measures shall be immediately taken, stop the treatment as well as provide proper medications. Change the dialyzer and the blood line

Product Performance: This dialyzer has reliable performance, which can be used for hemodialysis. The basic parameters of product performance and the laboratory data of this series will be provided as follows for reference.

Note: The laboratory data of this dialyzer was measured according to the standards ISO8637

Table 1 Basic Parameters of Product Performance

Model					
Product ref.					
Sterilization Way	Gamma ray	Gamma ray	Gamma ray	Gamma ray	Gamma ray
Effective membrane area (m ²)	1.4	1.6	1.8	2.0	2.2
Maximum TMP (mmHg)	500	500	500	500	500
Pressure drop of blood compartment (mmHg) Q _B =200mL/min	<75	<60	<50	<45	<45
Pressure drop of blood compartment (mmHg) Q _B =300mL/min	<105	<95	<85	<75	<70
Pressure drop of blood compartment (mmHg) Q _B =400mL/min	<110	<100	<90	<80	<70
Pressure drop of dialysate compartment (mmHg) Q _D =500mL/min	<45	<40	<40	<40	<40
Priming volume (mL)	76	88	98	110	120

Table 2 Series H (Clearance)

Model					
Product ref.					
Test Conditions: Q _D = 500mL/min, temperature: 37°C ±1°C, Q _F = 10 mL/ min					

Clearance (mL/min) $Q_B = 200 \text{ mL/min}$	Urea	193	195	197	198	199
	Creatinine	181	185	190	195	199
	Phosphate	176	181	185	190	195
	Vitamin B ₁₂	147	155	165	175	180
	Inulin	115	118	125	141	156
Clearance (mL/min) $Q_B = 300 \text{ mL/min}$	Urea	260	267	275	280	285
	Creatinine	235	250	260	270	280
	Phosphate	220	235	250	262	267
	Vitamin B ₁₂	170	175	195	208	220
	Inulin	120	130	140	160	180
Clearance (mL/min) $Q_B = 400 \text{ mL/min}$	Urea	310	320	330	340	350
	Creatinine	280	295	305	317	330
	Phosphate	260	265	285	297	310
	Vitamin B ₁₂	180	190	200	212	223
	Inulin	130	142	155	167	183

Remark: The tolerance of the clearance data is $\pm 10\%$.

Table 3 Series H (Ultrafiltration Coefficient)

Model					
Product ref.					
$Q_B = 300 \text{ mL/min}$	62 mL/hr/ mmHg	71 mL/hr/ mmHg	78 mL/hr/ mmHg	84 mL/hr/ mmHg	89 mL/hr/ mmHg


The Ultrafiltration Coefficients were measured in the condition of $Q_B = 300 \text{ mL/min}$, $\text{TMP} = 50 \text{ mmHg}$.

Table 4 Series H (Sieving Coefficient)

Model					
Product ref.					
Sieving coefficients (S.C.) $Q_B = 200 \text{ mL/min}$ $Q_F = 30 \text{ mL/min}$	Inulin: $1.0 \pm 10\%$ β_2 microglobulin: ≥ 0.7 Myoglobin: ≥ 0.55 Albumin: ≤ 0.01				

Main Structure:

Component	Housing	Dialysis membrane	Cover	O-ring	Sealant	Protective cap
Material						

Symbol clarification:

	Single use only		Sterilized using irradiation
	Manufacture date		Lot Number
	Use by		Temperature limitation
	Caution		Do not use if package is damaged
	CE MARK		Authorised Representative in European Community

Quality Assurance and Responsibility Limitation

- The manufacturer hereby guarantees that this dialyzer has been manufactured completely according to the manufacturer's enterprise standard, it complies with the manufacturing process specified in the standards ISO 8637.
- If the manufacturer receives any notice within 3 years since the date of sterilization, the defective dialyzer can be replaced according to the lot number printed on the package of the product.
- The manufacturer is free of responsibilities on condition that: a) The dialyzer is lost in the place where the product is used b) The dialyzer is directly or indirectly damaged artificially. c) The product failure results from the operation error of the user, carelessly reading and comprehending the correspond instruction and indication for use before using.
- The manufacturer is responsible for the quality of the product when it is firstly used only if the preparation and operation described in this instruction for use are implemented.

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[Redacted]

Mechanical Test Report

Manufacturer: [Redacted]
Product name: [Redacted]
Brand name: [Redacted]
Reference name: [Redacted]

1. Introduction

The target of this test is to verify the compliance of the performance of ready-to-sale, i.e. sterilized hemodialyzers with the normative requirements regarding their mechanical properties. The following mechanical properties are tested:

- Effective membrane area
- Particle pollution
- Positive pressure test
- Negative pressure test
- The occlusive of the blood chamber
- Maximum transmembrane pressure
- Temperature resistance
- Clearance
- Seiving coefficient
- Ultrafiltration coefficient
- Blood chamber volume
- Pressure drop of the blood chamber
- Pressure of the dialysate compartment

2. Test Period

The start time of the test: March. 10th, 2011
The end time of the test: May. 10th, 2011

3. Reference standard:

3.1 ISO 8637:2010 Cardiovascular implants and artificial organs-Haemodialysers, haemodiafilters, haemofilters and haemoconcentrators
3.2 In-house test standards

4. Mechanical Test

4.1 Effective membrane area
4.1.1 Test Standard
Test according to in-house test standard
4.1.2 Test method
Calculate the effective membrane area according to the below formula

$$S = \pi D L n \cdot 10^{-9}$$

[Redacted]

stops increasing and the blood compartment has some blood leaking, record the TMP value, and this value is the Maximum Transmembrane Pressure value.

4.6.3 Test results

Product ref	Test data	Criteria	Conclusion
	528 mmHg	≥500 mmHg	Complies with
	525 mmHg	≥500 mmHg	Complies with
	538 mmHg	≥500 mmHg	Complies with
	522 mmHg	≥500 mmHg	Complies with
	525 mmHg	≥500 mmHg	Complies with

4.7 Temperature resistance

4.7.1 Test Standard

Test according to in-house test standard

4.7.2 Test method

Place the hollow fiber dialyzer into a 0 °C refrigerator for 3min, then place it into a 60 °C constant temperature incubator for 3h, let it return to room temperature after removal to observe, repeat the test of 4.1.

4.7.3 Test results

Product ref	Test data	Criteria	Conclusion
	Deformation and fracture is not observed	Deformation and fracture phenomena should not be observed when the product in the temperature range of 0 – 60 °C	Complies with
	Deformation and fracture is not observed		Complies with
	Deformation and fracture is not observed		Complies with
	Deformation and fracture is not observed		Complies with
	Deformation and fracture is not observed		Complies with

4.8 Clearance


4.8.1 Test Standard


Test according to ISO 8637:2010 as the below method.

4.8.2 Test method

4.8.2.1 Test solution

Fill in the blood compartment with imitation solution containing the substances listed in the following Table 1. Fill in the dialysate compartment of hemodialyzer with dialysate.

4.8.2.3 Calculation Formula of Clearance

Application formula of hemodialysis clearance K

$$K = \left(\frac{C_{in} - C_{out}}{C_{in}} \right) Q_{in} + \frac{C_{out}}{C_{in}} Q_f$$

Where

- C_{in} the solution concentration at blood inlet of hemodialyzer
- C_{out} the solution concentration at blood outlet of hemodialyzer
- Q_{in} the blood flow rate at the product inlet
- Q_f the flow rate of filtration (ultrafiltration rate)
- C_{in} and C_{out} use the same concentration unit.

4.8.3 Test results

Product ref		[Redacted]										Conclusion
Test data & Criteria		Test date	Crite res	Test date	Crite res	Test date	Crite res	Test date	Crite res	Test date	Crite res	
$Q_D = 200 \text{ mL/min}$	Urea											Complies with
	Creatinine											Complies with
	Phosphate											Complies with
	Vitamin B ₁₂											Complies with
	Inulin											Complies with
$Q_D = 300 \text{ mL/min}$	Urea											Complies with
	Creatinine											Complies with
	Phosphate											Complies with
	Vitamin B ₁₂											Complies with
	Inulin											Complies with
$Q_D = 400 \text{ mL/min}$	Urea											Complies with
	Creatinine											Complies with
	Phosphate											Complies with
	Vitamin B ₁₂											Complies with
	Inulin											Complies with


☞ Test Conditions: $Q_D = 500 \text{ mL/min}$, temperature: $37 \text{ }^\circ\text{C} \pm 1 \text{ }^\circ\text{C}$, $Q_f = 10 \text{ mL/min}$

☞ Remark: The tolerance of the clearance data is $\pm 10\%$.

4.9 Seiving coefficient

4.9.1 Test Standard

Test according to ISO 8637:2010 as the below method.

4.9.3 Test results

Product ref.	Test data	Criteria	Conclusion
	Inulin : 0.91 β2- microglobulin : ≥0.80 Myoglobin : ≥0.62 Albumin : ≤0.001	Inulin : 0.9 ± 10% β2- microglobulin : ≥0.7 Myoglobin : ≥0.55 Albumin : ≤0.01	Complies with
	Inulin : 0.91 β2- microglobulin : ≥0.80 Myoglobin : ≥0.64 Albumin : ≤0.001		Complies with
	Inulin : 0.92 β2- microglobulin : ≥0.80 Myoglobin : ≥0.60 Albumin : ≤0.002		Complies with
	Inulin : 0.92 β2- microglobulin : ≥0.80 Myoglobin : ≥0.60 Albumin : ≤0.002		Complies with
	Inulin : 0.92 β2- microglobulin : ≥0.80 Myoglobin : ≥0.65 Albumin : ≤0.002		Complies with

⇨ Q_b = 200mL/min , Q_d = 30mL/min

4.10 Ultrafiltration coefficient

4.10.1 Test Standard

Test according to ISO 8637:2010 as the below method.

4.10.2 Test method

4.10.2.1 Test solution


The test fluid is anticoagulated bovine plasma or artificial blood, with a protein concentration of 60 g/L ± 5g/L. The dialysate compartment or filtration compartment shall not be filled with solution.

4.10.2.2 Test steps

Set up the test circuit as shown in figure 2. Establish stable blood and dialysate flow rate (including temperature, flow rate and pressure). Measure the ultrafiltration rate. Measure the value of ultrafiltration rate according to the transmembrane pressure from low to high.

4.10.3 Test results

Product ref.	Test data	Criteria	Conclusion
	61 mL/hr/mmHg	62mL/hr/mmHg ± 20%	Complies with
	71 mL/hr/mmHg	71 mL/hr/mmHg ± 20%	Complies with
	78 mL/hr/mmHg	78mL/hr/mmHg ± 20%	Complies with
	83 mL/hr/mmHg	84mL/hr/mmHg ± 20%	Complies with
	88 mL/hr/mmHg	89mL/hr/mmHg ± 20%	Complies with

☞ The Ultrafiltration Coefficients were measured in the condition of $Q_b = 300\text{mL/min}$, $\text{TMP} = 50\text{mmHg}$.

4.11 Blood chamber volume

4.11.1 Test Standard

Test according to the method regulated in ISO 8637:2010.

4.11.2 Test method

4.11.3 Test results

Product ref	Test data	Criteria	Conclusion
	74- 78mL	76 mL $\pm 10\%$	Complies with
	86- 90mL	88 mL $\pm 10\%$	Complies with
	96- 100mL	98 mL $\pm 10\%$	Complies with
	108- 112mL	110 mL $\pm 10\%$	Complies with
	118- 122mL	120 mL $\pm 10\%$	Complies with

4.12 Pressure drop of the blood chamber

4.12.1 Test Standard

Test according to the ISO 8637:2010 as below method.

4.12.2 Test method

4.12.2.1 Test solution

Fill the blood chamber with anticoagulant bovine plasma or artificial blood with protein concentration of $60\text{g/L} \pm 5\text{g/L}$.

fill the dialysate chamber with the usual dialysate.

4.12.2.2 Test steps

When $Q_b = 200\text{ml/min}$, 300 ml/min , 400 ml/min at the temperature $37\text{ }^\circ\text{C} \pm 1\text{ }^\circ\text{C}$, read the pressure values at the entrance and exit site of the blood compartment, then calculate the pressure drop.

4.12.3 Test results

Product ref	Test data	Criteria	Conclusion
	$Q_b = 200\text{ml/min} : 60\text{ mmHg}$ $Q_b = 300\text{ ml/min} : 90\text{ mmHg}$ $Q_b = 400\text{ ml/min} : 95\text{mmHg}$	$Q_b = 200\text{ml/min} : <75\text{mmHg}$ $Q_b = 300\text{ ml/min} : <105\text{ mmHg}$ $Q_b = 400\text{ ml/min} : <110\text{ mmHg}$	Complies with
	$Q_b = 200\text{ml/min} : 50\text{ mmHg}$ $Q_b = 300\text{ ml/min} : 85\text{ mmHg}$ $Q_b = 400\text{ ml/min} : 90\text{mmHg}$	$Q_b = 200\text{ml/min} : <60\text{mmHg}$ $Q_b = 300\text{ ml/min} : <95\text{ mmHg}$ $Q_b = 400\text{ ml/min} : <100\text{ mmHg}$	Complies with
	$Q_b = 200\text{ml/min} : 40\text{ mmHg}$ $Q_b = 300\text{ ml/min} : 75\text{mmHg}$ $Q_b = 400\text{ ml/min} : 80\text{mmHg}$	$Q_b = 200\text{ml/min} : <50\text{mmHg}$ $Q_b = 300\text{ ml/min} : <85\text{ mmHg}$ $Q_b = 400\text{ ml/min} : <90\text{ mmHg}$	Complies with
	$Q_b = 200\text{ml/min} : 40\text{ mmHg}$ $Q_b = 300\text{ ml/min} : 65\text{ mmHg}$	$Q_b = 200\text{ml/min} : <45\text{mmHg}$ $Q_b = 300\text{ ml/min} : <75\text{ mmHg}$	Complies with

[Redacted]

	$Q_D = 400 \text{ ml/min} : 70 \text{ mmHg}$	$Q_D = 400 \text{ ml/min} : < 80 \text{ mmHg}$	
[Redacted]	$Q_D = 200 \text{ ml/min} : 34 \text{ mmHg}$ $Q_D = 300 \text{ ml/min} : 60 \text{ mmHg}$ $Q_D = 400 \text{ ml/min} : 62 \text{ mmHg}$	$Q_D = 200 \text{ ml/min} : < 45 \text{ mmHg}$ $Q_D = 300 \text{ ml/min} : < 70 \text{ mmHg}$ $Q_D = 400 \text{ ml/min} : < 70 \text{ mmHg}$	Complies with

4.13 Pressure of the dialysate compartment

4.13.1 Test Standard

Test according to the ISO 8637:2010 as below method

4.13.2 Test method

4.13.2.1 Test solution

Fill in the dialysate compartment with the usual dialysate.

Fill the blood compartment with anticoagulant bovine plasma or artificial blood and seal it.

4.13.2.2 Test steps

When $Q_D = 500 \text{ ml/min}$, at the temperature $37^\circ\text{C} \pm 1^\circ\text{C}$, read the pressure values at the entrance and exit site of the dialysate compartment, then calculate the pressure drop.

4.13.3 Test results

Product ref.	Test data	Criteria	Conclusion
[Redacted]	30 mmHg	<45 mmHg	Complies with
[Redacted]	30 mmHg	<40 mmHg	Complies with
[Redacted]	30 mmHg	<40 mmHg	Complies with
[Redacted]	30 mmHg	<40 mmHg	Complies with
[Redacted]	32 mmHg	<40 mmHg	Complies with

5. Conclusion

The mechanical tests carried out according to ISO 8637:2010 and in-house test standards on the hemodialyzers manufactured by [Redacted] show, that with regard to their mechanical properties of the hemodialyzers, these products clearly comply with the normative requirements.

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Testing & Development Center for Dental Materials, [redacted]
 Use only for submission to National FDA or Quality Assurance. Page 1

Technical Report [redacted] 

Title : Torque Force Test of Implant (스크류회전전단토크)

I. File Information
 Document Number : [redacted] Tester : [redacted]
 Date of test : [redacted] Date of print : [redacted]

II. Product Information
 Product Name : [redacted]
 Materials : [redacted]
 Model Number : [redacted]
 Lot Number : n/s

III. Test Methods
 Room Temperature : 23 ℃
 Normative Reference : GGP 86, Section IB1b.
 Purpose : Measurement of torque force between abutment and fixture by screw.
 Procedure : Torque Gauge
 Center Protocol Code : n/a

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VI. LAB note
 Lab. Meeting : 2007-02-28, OK

VII. Attachments
 1. Nothing.

Respectfully submitted by [redacted]
 Testing and Development Center for Dental Materials
 Testing Laboratory Registered at KPDA

Technical Report Form #: [redacted]

Center for Dental Materials, Kyung Hee University, 21, Haeji-Dong, Dong-Deam-Noon, Seoul, Korea, 130-701, TEL. 02-961-6330, Fax 02-963-2427, E-MAIL: kdm@khu.ac.kr. This report is issued on the test articles submitted by client, and can not used for the purpose of advertisement and judgement etc., and never be copied without the published permission of a head of center.

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Splints and stress transmission to teeth: an in vitro experiment

Received 12 October 1999; accepted 10 March 2000

Abstract

Objective. To determine the influence of hard and soft splints with two thicknesses on the stress transmission to the tooth supporting the splint and the opposite tooth.

Methods. Continuous vertical forces up to 500 N were applied to two opposite first molar phantom teeth using a universal loading machine. Deformation was detected by strain gauges attached to the cervical area of the buccal and lingual aspects of the lower tooth. Strain, as a function of force, was collected and the slope, defined as the compliance (in $\mu\text{S/N}$) of the system, was calculated.

Results. The highest compliance was found with hard splints. When splints were constructed on the upper molar, the highest compressive compliance was registered on the buccal side ($2.8 \mu\text{S/N}$) and tension compliance on the lingual side ($-0.35 \mu\text{S/N}$). When constructed on the lower tooth, the opposite was found. Soft splints resulted in compression on both the buccal and lingual sides when adjusted to the upper or lower tooth. A higher compliance was found on the buccal side ($1.26 \mu\text{S/N}$), while on the lingual side, the values varied ($0.48-0.78 \mu\text{S/N}$).

Conclusions. Soft splints are more efficient in protecting teeth against the damage of bending forces; although there is an increase of compression forces. The tooth opposing a hard splint is exposed to a higher risk of bending forces. © 2000 Elsevier Science Ltd. All rights reserved.

Keywords: Bruxism; Hard occlusal splints; Soft occlusal splints; Stress transmission; Strain gauge

1. Introduction

One of the most problematic parafunctional behaviors is bruxism, which is defined as clenching or grinding of the teeth during sleep, or during daytime when the individual is not chewing or swallowing [1]. Fractures of dental restorations and tooth structures can be the result of bruxism [2]. This can be very frustrating for the patient, as well as for the dentist, to see their work compromised within a few months. Often a dentist deliberates whether to propose complex and expensive rehabilitation, especially implants, to a heavy bruxer patient.

In the past 25 years, epidemiological studies have underscored the high incidence of bruxism. In a study conducted on groups of adults and children, it was concluded that the phenomenon was widespread enough to warrant including it as a public health problem [3].

During nocturnal bruxism activity, bite forces are much greater than those exerted during mastication [4]. The mean maximal bite force registrations among Europeans and Americans is in the range of 600–750 N [5], while

masticatory forces are much lower (about 60–100 N) [6]. During clenching, forces of about 1000–1500 N, or higher, have been reported [7]. In some individuals, the forces generated and the duration of bruxing episodes can produce considerable loads within the masticatory system [8].

The most common damage to the dental unit is the wearing out of tooth surfaces (facets) [9], idiopathic cervical erosions [10], cusp or restoration fractures (fillings, crowns or porcelain bridges), and tooth mobility. Under excessive pressure, as part of the self-protective features of the dentition, the periodontal and/or pulpal receptors induce negative feedback on the activity of the jaw elevators [11]. Apparently, there is greater strain transmission to restorations on osseointegrated implants that have no proprioceptive receptors [12].

Direction of force resulting from a given occlusal load is important. Forces transmitted in an axial direction are best withstood. Horizontal forces are not tolerated as well as vertical forces [13] and can cause bending of the tooth, thus creating compressive and tensile stresses. Both dentin and enamel have high compressive strength, but their ability to withstand tension is limited [14]. The tensile force, which is highest at the tooth fulcrum (cervical area), could cause disruption of the chemical bonds between hydroxyapatite

* Corresponding author. Tel: +972-3-6409112; fax: +972-3-6409250.


Fig. 1. Two phantom first molars mounted on two acrylic blocks in an intercuspal position. Strain gauges bonded to the buccal cervical area of the lower tooth and connected to a strain indicator.


Fig. 2. Hard thick splint fitted to the upper or lower molar and adjusted to the opposite tooth according to the principles of ideal occlusion. Buccal view.

crystals, and could also be responsible for cervical erosion [15].

Occlusal splints have been advocated to reduce the harmful sequel to bruxism [16], especially when the goal is to protect large restorations. Splints may reduce the total amount of nocturnal activity in bruxers by decreasing the number of bruxing events per night. A splint may also reduce daytime postural muscle activity [6]. However, it is not actually known whether splints can absorb some of the stresses transmitted to the teeth. There is no universal agreement as to the effect of hard and soft splints on bruxism. It has been reported that hard stabilization splints significantly decrease muscle activity [17], while the vacuum-formed soft splints increase muscle activity [18,19].

The purpose of this study was to test the effect of splints—made of hard acrylic and soft vacuum formed—on the strain developed at the tooth cervix. This effect was analyzed on the tooth bearing the splint and the opposite tooth, with reference to the material used and the thickness of the splint.

In several studies, the comparative mechanical properties

of different materials have been examined for use as protective mouth devices, during sports activity [20,21,23]. To the best of our knowledge, this is the only *in vitro* study where tooth behavior using splints under simulation of occlusal forces, has been investigated.

2. Materials and methods

Two opposite phantom first molars were mounted onto two acrylic blocks in an intercuspal contact position and in an angulation similar to the natural dentition. Teeth were adjusted by occlusal equilibration to fit the occlusal contact according to the principles of ideal occlusion [13]. Two linear miniature strain gauges (EA-06-031EC-350, Measurements Group, Raleigh, NC, USA) were bonded to the cervical area of the buccal and lingual aspects of the first lower molar and connected to a strain indicator (Vishay 2100, Measurements Group, Raleigh, NC, USA) (Fig. 1). Splints were constructed on either the upper (U) or the lower (L) molar tooth. Two hard (H) splints covered the occlusal surface of the molar in 6 mm (T) and 2 mm (t) thicknesses


Fig. 3. Hard thick splint fitted to the upper or lower molar and adjusted to the opposite tooth according to the principles of ideal occlusion. Buccal view.

from self-curing dental acrylic (Orthocryl 2000-Dentaurum, Ispringen, Germany), and two soft (S) vacuum-formed splints, "150 Mouthgard" for the thick (T) and "040 soft EVA" for the thin (t) splint (DR Dental Resources, Delano, MN). For each combination of material, thickness and location, three initials were used for identification. Ten splints were prepared for HTL and HtL, and five splints for each of the following combinations: STL, STU, StL, StU, HTU and HtU, for a total of 50 splints. Each splint was adjusted according to the principles of ideal occlusion, i.e. where the supporting cusps (upper lingual or lower buccal) made contact with the flat occlusal surface of the splint on the opposite tooth (Figs. 2 and 3).

The teeth were mounted on a special jig connected to a loading machine (Instron 4502, High Wycombe Buckinghamshire, England). Compressive forces up to 500 N were applied using a 5 mm/min cross head speed. Since the teeth were mounted in an angulation similar to the natural dentition, compressive forces are translated at the tooth level as vertical and horizontal components that simulate the clinical condition. Strains, as a function of force, were automatically acquired upon loading by a Data Acquisition

System (ViewDac, Kiethley/Assist, Taunton, MA, USA). Each trial consisted of three consecutive loading episodes. At the end of each loading cycle, the splint was removed and resealed. The mean slope, defined as the compliance in $\mu\text{S}/\text{N}$ of the system, was calculated using the Data Acquisition System (ViewDac program).

Two-way analysis of variance (ANOVA) with data screening was used to determine the mean differences in compliance between thicknesses (thick and thin) and location (upper or lower) of the splints.

3. Results

Fig. 4 presents a characteristic diagram of strain versus force obtained from the HTU splint. All the diagrams presented a linear characteristic and therefore the slope (compliance) of each one could be obtained. The mean ($\pm\text{SD}$) of the buccal and lingual compliance of all the experiments are presented in Fig. 5. Positive values indicate the compression developed at the strain gauge area, while negative values represent tension.

Compressive strains were registered from the buccal and lingual aspects of the lower molar when there was no splint between the teeth. Compliance on the buccal and lingual aspects was 1.27 and 0.15 $\mu\text{S}/\text{N}$, respectively.

When the hard, thick and thin splints were constructed on the upper molar (HTU and HtU), the highest compressive compliance was registered on the buccal side (2.8 $\mu\text{S}/\text{N}$), and tension compliance on the lingual side. Values from the buccal side were about 5-fold higher than those from the lingual side. Furthermore, values on the buccal side were more than twice of those that were obtained without a splint, while the compliance on the lingual side was more than 3-fold higher and in the opposite direction (tension). When the hard splints were adjusted to the lower tooth, the opposite effect was found. Compression was registered on the lingual side and tension on the buccal side. Values on the lingual side were 6-fold higher than tension on the buccal side. Two-way ANOVA revealed a statistically significant main effect for the location ($p < 0.001$) of the hard splint and no effect for the thickness on the buccal and lingual compliance. No interaction was found between the thickness and the location.

Contrary to the hard splints, soft splints showed compression on both buccal and lingual sides, regardless of the location. Higher compliance (between 1.5- and 2.5-fold) was found on the buccal side compared to the lingual. Regarding the buccal compliance, the two-way ANOVA test revealed a statistically significant main effect for the thickness ($p = 0.023$) and no effect for the location of the splint. However, regarding the lingual compliance, the two-way ANOVA test revealed a statistically significant main effect for both the thickness ($p = 0.0001$) and the location ($p = 0.0001$). For both buccal and lingual compliance, no


Fig. 4. Linear characteristic diagram of strain (y-axis) versus force (x-axis) obtained from the hard thick upper splint.


Fig. 5. Mean and SD of the buccal and lingual compliance for each type of splint: N, without splint; H, hard acrylic splint; S, soft vacuum proform splint; T, thick; t, thin; U, splint fitted on the upper molar; L, splint fitted on the lower molar.

interaction was found between the thickness and the location.

4. Discussion

This preliminary *in vitro* study, uses two opposite phantom head teeth to evaluate the differences in stress transmission to the tooth, as a function of the different kinds of customized splints. Phantom head teeth simulate the shape and size of natural teeth and present an advantage to using a matched pair of opposite first molars, as in this experiment. Although these teeth have very different stiffness properties when compared with natural teeth, they only serve as an anchor to the strain gauge and the splint in this model, and enable an analysis of the influence of various splints on stress transmission while all other parameters are steady. The same model could use natural teeth, but it should be remembered that biological homeostasis is impossible to retain in an *in vitro* model. In addition, phantom head teeth and natural extracted unmatched teeth cannot simulate the PDL as the biological shock absorber.

The strain-gauge model is a very sensitive technique used to analyze tooth behavior under physiologic loads [23,24]. A higher compliance expresses greater deformation at the site of the strain gauge for the same applied force. When there is compression on one side of the tooth and tension on the other, it should be implied that bending occurs toward the side of the compression proportional to the compliance difference.

In the present system, the location of the point of force application (occlusal contacts) relative to the strain gauges, and the relative angulation of the two opposing teeth could influence the amount of strain registered by the strain gauges. When there is no splint, the occlusal relationship of the opposite tooth is established following the principles of ideal occlusion (four occlusal contacts on the lower tooth: two on the buccal cusps and two on the central fossa line) providing compressive strains developing at both sides. The closeness of the occlusal contacts to the strain gauge attached to the buccal aspect of the molar explains the greater deformation detected on that side, when compared to the lingual.

A hard splint has a flat occlusal surface that is adjusted to provide precise contact with the supporting opposite cusps (two occlusal contacts). When the hard splint fits over the upper tooth, bending of the lower tooth occurs toward the buccal side, since two contacts are established between the central portion of the occlusal surface of the splint and the buccal cusp tips of the lower tooth. In this case, displacement of the point of force application occurs toward the buccal side of the lower tooth. When the hard splint fits over the lower tooth, occlusal contacts occur between the lingual upper cusps and the central portion of the surface of the splint, causing a shift in force application toward the lingual side.

In the presence of the splint, the same amount of force has to be transferred to the lower tooth through two contact points instead of four, without a splint, causing higher contact stresses. It results in higher compressive forces being developed at the buccal or lingual aspects of the lower tooth.

This *in vitro* single tooth model study has some clinical application. In cases where a hard splint is located on the upper arch, our study model provides a relatively good simulation of the clinical behavior of the teeth that oppose the splint. The tooth reacts as an independent entity with minimal force transmission to the adjacent teeth. When the splint was located on the lower tooth (a simulation of a tooth-bearing splint), it acted independently in this *in vitro* study. However, in the clinic, a full coverage arch splint could have a significant impact on force distribution and could act in a different manner.

With soft splints, there are contact areas with the opposing tooth because of the resilient properties of the material. Therefore, there is a better distribution of the occlusal force, better energy absorption [20,22], and almost no change in compliance, which leads to a decreased bending effect on the teeth. The thick soft splints relative to thin ones decreased buccal compliance and increased the lingual causing the additional reduction of bending forces. Hypothetically, this finding is caused due to the better "spreading" capability of the material over the tooth surface.

From the results of this *in vitro* study, soft splints provided the best protective capacities against destructive bending forces, although the total compressive forces increased.

The use of splints in the treatment of bruxism ensures a reduction in the wearing out of the tooth and leads to a better force distribution in the dental arch supporting the splint. Moreover, the splint increases cognitive awareness and leads to a behavioral modification via the reduction in muscular activity.

The effect of soft splints on muscular activity have been investigated in TMD patients. Soft appliances have been shown to increase nocturnal masseteric EMG activity in 50% of cases, while 80% of cases had a significant reduction of nocturnal muscular activity with hard splints [17]. In other studies, a higher therapeutic effect of hard splints when compared to soft ones for reducing TMD symptoms has been shown [18,19,25]. This effect may be based on behavioral modification. In daytime wearing out of soft splints, subjects are exposed to a considerable amount of conscious or unconscious nervous chewing [20,26]. Our model does not enable an estimation of these behaviors. We strongly believe that clinically controlled studies should focus on the stress transmission to the tooth bearing the splint and the one opposing the splint, to better understand the value of this widely used appliance, at the tooth level. In the light of the surprising results of this *in vitro* preliminary study, one should ask whether the use of splints is justified.

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To Whom it May Concern:

We, [REDACTED] hereby certify that the scope of the performance test report, "[REDACTED] Technical Data Verification Specification and Report" (2012. 04.25) is the examination of the technical data of [REDACTED], [REDACTED], and [REDACTED]

In addition, each test unit used for each performance test listed in below table is basically equivalent to [REDACTED] in terms of software and performance specification.

Therefore, we can conclude that the results of tests are also applicable to [REDACTED] and there have been no significant change of product since the date of each test report that affects the safety and essential performance.

Table

Report no	Title item	Device under Test	Issue date
06-1318-PR07 Rev.00	Climatic Compatibility	[REDACTED]	May 19,2006
QR025003	Environmental condition		September 15,2003
VER2005-MT-1021	Specification Test		November 1,2005
QR025102	Gas delivery requirement		April 30, 2001
VER/2006/VV/105	Weight		July 21, 2006
QR025002E	Electrical and power requirement		April 21,2001
QR025002S, RevY	Operator Interface, Oxygen monitor, Pressure monitor, Volume monitor, Respiratory monitor, Gas Delivery		December 2,2003

Report no	Title item	Device under Test	Issue date
	subsystem, Ventilator, Alarm, Miscellaneous functionality		
QR025002S, RevA	Operator Interface, Oxygen monitor, Pressure monitor, Volume monitor, Respiratory rate, Fresh gas flows, Ventilator, Alarm, Miscellaneous functionality		May 7, 2001
QR025002S, RevC	Operator interface, Monitor, Ventilator, Alarm functionality		July 6, 2001
QR025101	Ventilation module		April 30, 2001
2009/VV/346	Evaluation of the fresh gas flow measurement accuracy with Fabius GS equipped with fresh gas flow sensors from Honeywell and Sensirion		Sept 14, 2009 (date of test execution)

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DET NORSKE VERITAS MANAGEMENT SYSTEM CERTIFICATE

Certificate No. 70994-2008-AQ-IND-NA

This is to certify that the Management System of

at

has been found to conform to the Quality Management System Standard(s):

ISO 9001:2008 / ISO 13485:2003

This Certificate is valid for the following product or service category:

Design, Manufacture and Supply of Drug Fluting & Bare Metal Vascular Stents, Stent Systems, Inflation Devices and PTCA Balloon Catheters.

Initial Certification date:

This Certificate is valid until:

The audit has been performed under the supervision of:

**NORWEGIAN
ACCREDITATION
MSYS 002**

Place and date:

Hovik, 29 July 2011

For the accredited unit:

DET NORSKE VERITAS
CERTIFICATION AS, NORWAY

Aud Loken Ekklid
Management Representative

Conditions of fulfillment of conditions set out in the Certification Agreement may render this Certificate invalid.
 This Certificate may be applied to any of the following activities:
 www.dnv.com

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■ 업체명/어기번호 (Company name of Applicant / License No.)

[Redacted]

■ 대표자 (Representative)

[Redacted]

■ 업체 소재지 (Company address of Applicant)

[Redacted]

■ 제조소명 (Name of Manufacturer)

[Redacted]

■ 제조소 소재지 (Address of Manufacturer)

[Redacted]

■ 품목군 및 등급 : 불임 심조
(Category and Classification of Product(s) : See attached list)

의료기기 제조 및 품질관리기준에 적합함을 인정합니다.
(We hereby certify that the above manufacturer complies with Korea
Good Manufacturing Practices of Medical Devices for the product(s) listed above.)

발행일자(Date of Issue) : 2013. 01. 17.

유효기한(Date of Expiration) : 2013. 05. 16.

KFDA 식품의약품안전청장
Korea Food & Drug Administration

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Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization

has established and applies a quality management system for medical devices
for the following scope:

**Design, development, manufacture and distribution of medical devices for
cardiovascular interventions, peripheral interventions, minimally invasive
surgical procedures, endoscopy, gastrointestinal, gynaecological,
nephrological, biliary, urological and electrosurgical applications**

Facilities: see attachment

Proof has been furnished that the requirements specified in

EN ISO 13485:2003+AC:2009

are fulfilled. The quality management system is subject to yearly surveillance.

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K O R E A

Stability Test Report of Fresenius Helixone Dialysers


To whom it may concern

We hereby declare that the enclosed stability report is applicable to all below mentioned [redacted] Dialysers [Table I]. The raw materials and the whole manufacturing process, including packing and sterilization are the same as for Fresenius Helixone Dialysers.

[Table I] Product List: Fresenius Helixone Dialysers

Dialyser Type	Article number	Dialyser Type	Article number
[redacted]	5004831	FX CorDiax 40	F00001588
[redacted]	5004731	FX CorDiax 50	F00001589
[redacted]	5004741	FX CorDiax 60	F00001590
[redacted]	5004791	FX CorDiax 80	F00001591
[redacted]	5004801	FX CorDiax 100	F00001592
[redacted]	5008221	FX CorDiax 120	F00002384
[redacted]	5008841	FX CorDiax 600	F00001593
[redacted]	5008851	FX CorDiax 800	F00001594
[redacted]	5008861	FX CorDiax 1000	F00001595
[redacted]	5008881	FX50classix	F00002385
[redacted]	5008901	FX60classix	F00002386
[redacted]	5008881	FX80classix	F00002387
[redacted]	5008901	FX100classix	F00002388
[redacted]	5009721		

Bad Homburg, 08.01.2013
Place and date

i. A. Dorothea Hoffmann

Associate Medical Device Regulatory


Fresenius Medical Care AG & Co. KGaA
D-61346 Bad Homburg - Germany

i. V. Annette Manck

Manager Medical Device Regulatory

Fresenius Medical Care AG & Co. KGaA
D-61346 Bad Homburg - Germany

No. 19 of files for the year 2013

I hereby certify, that the above are the true signatures of

Mrs. **Annette Manck**, born April 18th, 1973

Mrs. **Dorothea Hoffmann**, born March 4th, 1966

both resident Else-Kröner-Straße 1, in 61352 Bad Homburg v.d. Höhe,
Germany

acting on behalf of

acknowledged in my presence on this day at Else-Kröner-Str. 1, 61352 Bad
Homburg v.d. Höhe, after I went there at the request of the appearing.

Mrs. Manck and Mrs. Hoffmann are personally known to me.

The notary asked whether he or any member of his firm or any person sharing
office with him had acted in the matter which is the subject of this instrument,
except in notarial capacity. The question was replied to in the negative.

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Chr. Walther
Notary

PRODUCT CENTER DIALYSERS			
Product verification and validation Development report		01	2011-04-18
Product category Haemodialysers/Haemofilters		Date 11/10/2010	Page 2 / 6

Table 1: Comparison of the performance criteria for type and QC testing of different guidances and standards

	FDA 1998	EN 1283 : 1996	ISO 9637 : 2010	ANSI/AAMI RD16: 2007	Realization with FME dialysers
Form of test	Bench tests 3 batches of each model		Prior to marketing or after changes	Prior to marketing or after changes	Type tests are regulated by internal SOP and must be released by the director of quality assurance and head of filter production after new developments or in the process of design change 3 batches of each filter type must be tested.
Performance Characteristics:					
Ultrafiltration rate (UF)	Test liquid: bovine or human blood - Measure at TMPs 0, 100, 300, 500 and TMP _{max} - K _{UF} should be calculated	Test liquid: bovine or human blood (HCT 52 ± 5%) Sequence of measurement from TMP _{min} to TMP _{max} at each blood flow rate stated by manufacturer	Test solution: bovine or human blood (HCT 52 ± 5%) Test over the manufacturer's specified range of TMP and blood flow rates	Test solution: bovine or human blood (HCT 50 ± 5%) Test over the manufacturer's specified range of TMP and blood flow rates	Only for low flux filter UF is measured with blood under the following condition at 2 different TMPs: A < 1m ² Q _{UF} 200 ml/min A > 1m ² Q _{UF} 300 ml/min Afterwards the UF is tested with H ₂ O to get the relation to the routine measurements FME dialysis machines 2003, 4003 and 5003 are all equipped with balancing chambers allowing precise volumetric flow control. This is a prerequisite for the use of highly permeable membranes in dialysis Also the product leaflet contains the warning: "the dialysers must only be used with dialysis machines which enable precise volumetric fluid control..." K _{UF} cannot be used to set treatment parameters. To provide a general view of the dialysers UF behaviour, the UF is determined according internal SOP at various TMPs.
HighFlux: Requires the use of an ultrafiltration controller	-Measure slope of UF rate versus TMP at UF rates between 600 and 1800 ml/h; (at least 4 data points) - K _{UF} should be measured using either bovine or human blood				
Clearances: Marker molecules	• Urea, • Creatinine • Vitamin B ₁₂ or Inulin	Test liquids with • Urea • Creatinine • Phosphate • Vitamin B ₁₂ • Inulin (for HDF)	Test solutions dialysate with • Urea 15-45 mmol/l • Creatinine 50-100 µmol/l • Phosphate 4-8 mmol/l • Vitamin B ₁₂ 5-10 µmol/l	Test solutions: dialysate with • Urea • Creatinine • Phosphate • Vitamin B ₁₂ • Inulin	Test solutions: RO water with NaCl Creatinine Phosphate Vitamin B ₁₂ Inulin
Measurement at HF-Q _{UF} 200 – 500 (ml/min) and TMP _{max} LF: 0 and TMP _{max}		Measure at each flow rate (blood + dialysing fluid) stated by manufacturer	Blood and dialysate flow rates shall cover the manufacturer's recommended range K _{UA} results can be included	Blood and dialysate flow rates shall cover the manufacturer's recommended range K _{UA} results can be included	In the routine Sodium is measured instead of Urea in water at Q _{UF} 200 ml/min or 300 ml/min dependent on the surface of the dialyser. K _{UA} for urea is included in leaflet and brochure. Clearances at different Q _{UF} can easily be calculated from K _{UA}
K _{UA}					

1. Title:	Determination of clearance performance data F and FX series	Doc. No. FIL 51/AA304E
2. Scope:	FMC-D GmbH, St. Wendel Plant, Production Units Dialysers and Membranes	Rev. No. 01
3. Purpose:	The work instruction specifies the procedure and conditions for the measurement of clearance data on the manual test stand.	Valid from: May 13, 2009

4. Processes and responsibilities:

Measuring principle:

Based on EN 1283, the input and output concentration of a specific key substance is measured during a simulated dialysis under constant conditions and the clearance calculated as follows:

$$C = Q_{B_{in}} \times \left(1 - \frac{C_{B_{out}}}{C_{B_{in}}} \right) + Q_F \times \frac{C_{B_{out}}}{C_{B_{in}}}$$

Measurement conditions:

	Filter types		
	FX paed AV paed	Low flux ≤ 1.5 m ² High flux ≤ 1.1 m ²	Low flux > 1.5 m ² High flux > 1.1 m ²
Temperatures [°C]	37 ± 1	37 ± 1	37 ± 1
Q _B [ml/min]	100	200 ± 4	300 ± 6
Q _D [ml/min]	300	500 ± 10	500 ± 10
Q _F [ml/min]	0	0ml/min ultrafiltration corrected clearance values	

Exception: FX 8, which is also measured at Q_B 300 [ml/min]

Test solutions:

a) Standard test solution:

9 g/l NaCl

50 mg/l Vitamin B₁₂

} dialysate-sided measured against water

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
Confirmation Letter for Performance Test Reports

To Whom It May Concern:


We,  hereby certify that the scope of the performance test report, "Fabius Family Technical Data Verification Specification and Report" (2012. 04.25) is the examination of the technical data of Fabius GS Premium, Fabius GS and Fabius Tiro.

In addition, each test unit used for each performance test listed in below table is basically equivalent to Fabius Tiro in terms of software and performance specification.

Therefore, we can conclude that the results of tests are also applicable to Fabius Tiro and there have been no significant change of product since the date of each test report that affects the safety and essential performance.

Table

Report no	Title Item	Device under Test	Issue date
06-1318-PR07 Rev.00	Climatic Compatibility	Fabius GS (part no 8607000-01)	May 19,2006
QR025003	Environmental condition	Fabius GS (part no 4116397)	September 15,2003
VER2005-MT-1021	Specification Test	Fabius GS and Fabius Tiro gasblock(8605831, 8605578)	November 1,2005
QR025102	Gas delivery requirement	Fabius GS Gas Delivery (part no 4116402)	April 30, 2001
VER/2006/VV/105	Weight	Fabius GS and Fabius Tiro	July 21, 2006
QR025002E	Electrical and power requirement	Fabius GS (part no 4116397)	April 21,2001
QR025002S, RevY	Operator Interface, Oxygen monitor, Pressure monitor, Volume monitor, Respiratory monitor, Gas Delivery	Fabius GS (part no 4116397)	December 2,2003

Report no	Title item	Device under Test	Issue date
	subsystem, Ventilator, Alarm, Miscellaneous functionality		
QR025002S, RevA	Operator interface, Oxygen monitor, Pressure monitor, Volume monitor, Respiratory rate, Fresh gas flows, Ventilator, Alarm, Miscellaneous functionality	Fabius GS (part no 4116397)	May 7, 2001
QR025002S, RevC	Operator interface, Monitor, Ventilator, Alarm functionality	Fabius GS Prototype P3	July 6, 2001
QR025101	Ventilation module	Fabius GS (part no 4116399)	April 30, 2001
2009/VV/346	Evaluation of the fresh gas flow measurement accuracy with Fabius GS equipped with fresh gas flow sensors from Honeywell and Sensirion	Fabius GS	Sept 14, 2009 (date of test execution)

2013-01-07

Ulrich Schröder
Director Quality & Regulatory Affairs

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
Fig. 5 Maximum biting force of 1st premolar of male.


Fig. 6 Maximum biting force of 1st premolar of female.

Table 4. Maximum biting force of female (unit : kg).

Age	No. of Sub.	LM	LP	AT	RP	RM
		Mean	Mean	Mean	Mean	Mean
20	22	33.10 ± 4.70	24.30 ± 2.98	12.10 ± 2.77	21.20 ± 4.21	37.70 ± 3.86
21	21	37.30 ± 4.72	21.00 ± 3.30	13.80 ± 2.39	21.30 ± 3.62	35.50 ± 3.44
22	23	29.50 ± 2.37	20.30 ± 3.80	9.60 ± 3.17	21.90 ± 4.53	29.70 ± 4.37
23	20	32.90 ± 5.26	20.80 ± 4.21	10.70 ± 3.09	20.70 ± 3.37	34.40 ± 3.92
24	22	31.90 ± 3.54	19.50 ± 6.75	10.60 ± 3.34	19.50 ± 3.54	33.00 ± 4.08
25	24	32.50 ± 3.87	19.60 ± 3.92	12.60 ± 2.27	18.90 ± 4.61	34.60 ± 5.25
Average	132	32.87 ± 4.63	20.92 ± 4.46	11.57 ± 3.09	20.58 ± 3.98	34.15 ± 5.71

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
Fig. 7 Maximum biting force of 1st molar of male.


Fig. 8 Maximum biting force of 1st molar of female

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Table 5. Maximum biting force of normal persons.

Region	Anterior	Premolar	Molar
	Mean	Mean	Mean
Male	13.67±2.65	22.02±3.64	32.45±4.12
Female	11.57±3.09	20.75±4.21	33.51±4.70

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Table 6. Difference of male and female.

	Sex	Mean
	AT	Male
Female		11.57±3.09
LP	Male	22.43±3.93
	Female	20.92±4.46
RP	Male	21.60±3.30
	Female	20.58±3.98
LM	Male	32.42±4.06
	Female	32.87±4.63
RM	Male	32.47±4.22
	Female	34.15±4.71

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Table 7. Difference of right and left region.

Sex	Region	Mean
Male	LM	32.42 ± 4.06
	RM	32.47 ± 4.22
	LP	22.43 ± 3.93
	RP	21.60 ± 3.30
Female	LM	32.87 ± 4.63
	RM	34.15 ± 4.71
	LP	20.92 ± 4.46
	RP	20.58 ± 3.98

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Table 8. Maximum biting force of the patient wearing the full denture.

	Mean
AT	5.10 ± 2.02
LP	7.30 ± 3.09
RP	7.50 ± 4.06
LM	13.10 ± 4.38
RM	13.30 ± 4.03

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Table 9. Maximum biting force of the patient wearing the upper & lower partial denture.

	Mean
LP	13.60 ± 4.86
RP	14.10 ± 4.82
LM	23.40 ± 7.79
RM	24.40 ± 8.90

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Table 10. Maximum biting force of the patient wearing the full denture and normal persons.

Region	Wearing the full denture	normal persons
Anterior	5.10 ± 2.02	12.62 ± 2.87
Premolar	7.30 ± 3.51	21.38 ± 3.93
Molar	13.20 ± 4.10	32.98 ± 4.41

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
Fig. 10. Comparison of maximum biting force of the patient wearing the full denture and normal persons.

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Table 11. Maximum biting force of the patient the wearing partial denture and normal persons.

Region	Wearing the Partial denture	normal persons
Premolar	13.85 ± 4.72	21.38 ± 3.93
Molar	23.90 ± 8.16	32.98 ± 4.41

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
Fig. 11. Comparison of maximum biting force of the patient wearing the partial denture and normal persons.

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
Page : 1 of 3

KTR KOREA TESTING & RESEARCH INSTITUTE

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[Redacted]

[Redacted]

has prepared for MFDS multiple consignments of:

Insoluble Type 1 Bovine Collagen (for pharmaceutical use with all tissues entirely originating in New Zealand) from June 2008 through to May 2010. These encompass the following Batches: **001, 002, 003**

Declaration 1: The bovine tendons were all derived from New Zealand born and raised bovines (breed Angus Cross) **aged 24 months or younger**, that had been inspected ante and post mortem and found to be fit for human consumption.

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Certificate of Analysis: [Redacted]

Product: [Redacted]

Type: [Redacted] product has been prepared by a combination of acid and alkali extraction of bovine flexor tendon. Residual tendon after extraction has been freeze dried as flat sheets. – **count 19 sheets.**

Origin: [Redacted] 100% pure New Zealand Origin – farmed bovines slaughtered in licensed premises inspected ante and post mortem and passed as fit for human consumption)

Consignment Net Weight/Volume: 121.4 grams
 Lot number: Batch # GE-025

Parameter:	Specification:	Result:
Description	Frozen bovine flexor tendon has been inspected, trimmed, washed, sliced and extracted in reagent grade acid, then alkali. The product is then neutralised, washed, ground, frozen and lyophilised. The product is a dry flat sheet.	Conforms
Intended use	Sample to confirm production process and suitability for further processing.	Eligible. <i>Note: Not to be fed to ruminant animals.</i>
Colour	Pale and clean	Conforms
Moisture content	Less than 12% (w/w)	Conforms
Additives	None	None
Preservatives	None	None
Bioburden	TVAC	To follow
Storage	Store at room temperature in dry conditions. Best before 28 February 2013 if stored correctly.	

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Process for Preparation of

Process Step	Conditions	Holding Time (Sequence)
Receive tendons from animal slaughter premises	Ensure tendons are frozen on receipt. Store in freezer until required	
Thaw	Thaw overnight at $^{\circ}\text{C}$. Complete thawing by immersing frozen material in running water ¹ if necessary.	hours (Day)
Inspect – confirm freedom from hair and fat	Set out tendons on a white tray and inspect individually.	hours (Day)
Wash in SDS	Wash the tendons by immersing them in a solution of sodium lauryl sulphate (SDS)	hour (Day)
Wash (water)	Wash the tendons in running RO water. Use at least 10 x the volume of SDS – wash until foam bubbles no longer appear.	hours (Day)
Freeze	Freeze the tendons individually (set out on a tray)	
Peel and trim	Remove the tendons from the freezer at a time. Peel by slicing off the outer layer from the two arms and discard the peelings. Cut the two arms from the body of the tendon. Also remove the narrow portion from the bottom of each arm.	minutes (Day)
Freeze	Lay the tendons straight and return to the freezer.	(Day)
Slice	Remove tendons from the freezer at a time. Form these freezers into a bundle in the slicer guide. Create fine slices (mm) and weigh.	hours (Day)
Acid extraction (pH)	Transfer the slices to a container. For each kilogram of sliced wet tendon add litres of water and litres of NH_4Cl / $\text{M Na}_2\text{SO}_4$ at room temperature ² . Mix gently. After hours pour off the mixture and drain the slices over a stainless steel colander.	hours (Day)
Wash (Na_2SO_4)	Return the slices to the container. For each kilogram of sliced wet tendon add litres of $\text{M Na}_2\text{SO}_4$ at room temperature. Mix gently. After minutes pour off the mixture and drain the slices over a stainless steel colander.	minutes (Day)
Alkali extraction (pH)	Return the slices to the container. For each kilogram of sliced wet tendon add 5 litres of M NaOH / $\text{M Na}_2\text{SO}_4$ at room temperature. Mix gently. After hours pour off the mixture and drain the slices over a	hours (Day)

Process Step	Conditions	Holding Time (Sequence)
	stainless steel colander.	
Wash (water) (pH)	Wash the slices in water until the pH of the water, after standing for one hour, is neutral. Pour off the water and drain the slices over a stainless steel colander.	hours (Days)
Wash (acid) (pH)	Wash the slices in N HCl until, after standing for hour, the pH is less than	hours (Day)
Grind/disperse	Grind the slices in a commercial blender with parts of NHCl to part of slices. Sieve out any un-ground material from the viscous fluid. Repeat the grinding process until a smooth dispersion is obtained.	hours (Day)
Freeze	Place the dispersed material in plastic bags – approximately litre per bag. Seal the bags and freeze flat on trays. (These trays are the drying shelves of the freeze dryer.)	hours (Days)
Dry	Transfer frozen packages to the freeze dryer. Cut open each bag so that vapour can move freely from all parts of the frozen sheet. Run the freeze dryer through a normal drying cycle for days. During this period the material warms from frozen (C) to °C.	(Days)
Package	Remove the now dry sheets from the dryer. Weigh each sheet and package with an inner bag which is in turn place in an outer bag which is vacuum sealed. The product is stored at room temperature	Day

Certified as accurate and in current use for the preparation for all

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KUREA TESTING & CERTIFICATION INSTITUTE

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Safety data sheet	08.07.2009 E.3
 Polypropylene HB120FB	
1. Identification of the substance/preparation and the company/undertaking	
Trade name:	<input type="text"/>
Intended use:	<input type="text"/>
Manufacturer:	<input type="text"/>
2. Composition/information on ingredients	
The product is a polypropylene homolyme <input type="text"/>	
Contains no substance classified as hazard <input type="text"/> which should be taken into account according to EC directives.	
3. Hazards identification	
Health: The product is not classified as dangerous preparation. However, inhalation of dust may irritate the respiratory tract. Prolonged inhalation of high doses of decomposition products may give headache or irritation of the respiratory tract.	
Fire: The product burns, but is not classified as flammable.	
Environment: The product is not considered dangerous for the environment.	
4. First-aid measures	
No specific instruction needed.	
5. Fire-fighting measures	
Suitable extinguishing media: Water in spread jet, dry chemicals, foam or carbon dioxide.	
Special exposure hazards: Principal toxicant in the smoke is carbon monoxide.	
6. Accidental release measures	
Suck or sweep up spill. All spill of granules must be removed immediately to prevent slipping accidents.	
7. Handling and storage	
Handling: During processing and thermal treatment of the product, small amounts of volatile hydrocarbons may be released. Provide adequate ventilation. Local exhaust ventilation may be necessary. Avoid inhalation of dust and decomposition fumes. Dust from the product gives a potential risk for dust explosion. All equipment shall be grounded.	
Storage: No other specific storage requirements.	
8. Exposure controls / personal protection	
Provide adequate ventilation. Local exhaust ventilation may be necessary.	
<input type="text"/>	<input type="text"/>
Page 1 of 3	

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Test Facility:

This report is confidential. No part may be used for advertising or public announcement without written permission. Results apply only to the sample(s) tested.

Report Number

March

P.O. #:

CHEMISTRY TEST REPORT

Sample Information:

Tip

Date Received: March 14, 2011

Date in Test: March 14, 2011

Date Completed: March 16, 2011

Test Information: Test Code:

Procedure #:

USP PHYSICOCHEMICAL EXTRACTION PARAMETERS				
Surface Area	Extraction Time	Extraction Media	Volume Used (mL)	Extraction Temp.
602 cm ²	24 Hours	Purified Water	100	70 °C

RESULTS

Parameter	Current USP Specifications	Results	
Non-Volatile Residue	≤ 15 mg	PASS	0.4 mg
Residue On Ignition	≤ 5 mg	Not Required	
		Non-Volatile Residue	< 5 mg
Heavy Metals	≤ 1 ppm (as Pb)	PASS	< 1 ppm
Buffering Capacity	≤ 10 mL 0.01 N Titrant	PASS	0.05 mL

Note: The extraction media was a clear and colorless solution for both the pre and post extracts.

Testing conducted in accordance with current Good Manufacturing Practices.

1285 Kennestone Circle • Marietta, GA 30066 • 888.847.6633 • 770.514.0262 • Fax 770.514.0294

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Page 1/4
Material Safety Data Sheet According to 91/155 EC
Printing date 27.11.2002 Reviewed on 04.10.2002
1 Identification of substance: <ul style="list-style-type: none">. Product details:. Trade name: <u>GRILAMID TR 55 LX NATUR</u>. Article number: 154903. Application of the substance / the preparation <i>Plastics processing industry</i>. Manufacturer/Supplier: <input type="text"/>. Emergency information: <input type="text"/>
2 Composition/Data on components: <ul style="list-style-type: none">. Chemical characterization. Description: <i>Polyamid 12 TR</i> <i>Polyamide 12</i>. Dangerous components: <i>Void</i>. Additional information <i>For the wording of the listed risk phrases refer to section 16.</i>
3 Hazards identification <ul style="list-style-type: none">. Hazard designation: <i>void</i>. Information pertaining to particular dangers for man and environment <i>The product does not have to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.</i>
4 First aid measures <ul style="list-style-type: none">. General information <i>No special measures required.</i>. After inhalation <i>Seek medical treatment in case of complaints.</i>. After skin contact <i>Do not pull solidified product away from the skin.</i> <i>After contact with the molten product, cool rapidly with cold water.</i>. After eye contact <i>Rinse opened eye for several minutes under running water. If symptoms persist, consult doctor.</i>. After swallowing <i>In case of persistent symptoms consult doctor.</i>
5 Fire fighting measures <ul style="list-style-type: none">. Suitable extinguishing agents <i>CO₂, extinguishing powder or water jet. Fight larger fires with water jet or alcohol-resistant foam.</i>. Special hazards caused by the material, its products of combustion or flue gases: <i>Can be released in case of fire</i> <i>Carbon monoxide (CO)</i> <i>Hydrogen cyanide (HCN)</i>
(Contd. on page 2)

Material Safety Data Sheet
According to 91/155 EC

Printing date 27.11.2002

Reviewed on 04.10.2002

Trade name: GRILAMID TR 55 LX NATUR

. Protective equipment: No special measures required.

(Contd. of page 1)

6 Accidental release measures

- . Person-related safety precautions: Not required.
- . Measures for environmental protection: No special measures required.
- . Measures for cleaning/collecting: Collect mechanically.
- . Additional information: No dangerous materials are released.

7 Handling and storage

- . Handling
 - . Information for safe handling: No special measures required.
 - . Information about protection against explosions and fires:
 - Dust can combine with air to form an explosive mixture.
- . Storage
 - . Requirements to be met by storerooms and containers: No special requirements.
 - . Information about storage in one common storage facility: Not required.
 - . Further information about storage conditions: None.

8 Exposure controls and personal protection

- . Additional information about design of technical systems: No further data; see item 7.
- . Components with critical values that require monitoring at the workplace:
 - The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- . Personal protective equipment
 - . General protective and hygienic measures
 - The usual precautionary measures should be adhered to in handling the chemicals.
 - . Breathing equipment: Not required.
 - . Protection of hands: Not required.
 - . Eye protection: Safety glasses recommended during refilling.

9 Physical and chemical properties:

. General Information

Form:	Granulate
Colour:	Different according to colour
Smell:	Nearly odourless

. Change in condition

. Melting point/Melting range:	Glass transition temperature 110°C
. Boiling point/Boiling range:	Not determined

. Flash point: Not applicable

. Self-inflammability: Product is not selfigniting.

. Danger of explosion: Product is not explosive.

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
Material Safety Data Sheet

Date of Preparation: 1/12/00 Revision: Original

Section 1 - Chemical Product and Company Identification		
Product/Chem Manufacturer		
Section 2 – HAZARDOUS INGREDIENTS/IDENTITY INFORMATION		
Hazardous Components/Ingredient Name	CAS Number	% wt or % vol
Dop/Dioctylphthalate (Di(2-ethylhexyl) phthalate)	117817	35%
<p><small>This product contains the following toxic chemical(s) subject to the reporting requirements of Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986, and Subpart C-Supplier Notification Requirement of 40 CFR Part 372.</small></p>		
Section 3 - Physical and Chemical Properties		
Physical State: Solid Appearance and Odor: Pellets with no odor or very mild non-offensive odor. Odor Threshold: N/A Vapor Pressure: N/A Vapor Density (Air=1): N/A Specific Gravity (H₂O=1): 1.20	Water Solubility: Essentially Insoluble Boiling Point: N/A Melting Point: N/A Evaporation Rate: Non-Volatile	
Section 4 - First Aid Measures		
Inhalation: Avoid breathing fine dust. Eye Contact: Wear safety goggles upon handling. Skin Contact: None Ingestion: None <i>After first aid, get appropriate in-plant, paramedic, or community medical support.</i> Note to Physicians: Special Precautions/Procedures: Use self-contained breathing apparatus if degradation occurs due to excessive heat. Stop degradation by water Quench.		
Section 5 - Fire-Fighting Measures		
Flash Point: N/A Extinguishing Media: Water Unusual Fire or Explosion Hazards: N/A Hazardous Combustion Products: HCL liberated in combustion. Fire-Fighting Instructions: Do not release runoff from fire control methods to sewers or waterways. Fire-Fighting Equipment: Because fire may produce toxic thermal decomposition products, wear a self-contained breathing apparatus (SCBA) with a full facepiece operated in pressure-demand or positive-pressure mode.		
Section 6 - Stability and Reactivity		
Stability: Material is stable below temperature of 400°F. Polymerization: Hazardous polymerization cannot occur. Chemical Incompatibilities: none Conditions to Avoid: Excessive heat.		

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SECTION 1 - CHEMICAL PRODUCT AND COMPANY IDENTIFICATION				
Product Name: Stainless Steel Grades: 16 - 6, 21 - 6 - 9, 303, 304 - 304L, 305, 309 - 309S, 310 - 310S, 316 - 316L, 317 - 317L, 321, 347 - 348, 403, 405, 410, 416, 430, 446.				
SECTION 2 - COMPOSITION/INFORMATION ON INGREDIENTS				
Ingredient	CAS #	% Weight	Occupational Exposure Limits (mg / m3) OSHA PEL	ACGIH TLV
Iron (Fe)	7439 - 89 - 8	39-88	10 (oxide, dust & fume)	10 (oxide,dust & fume)
*Chromium (Cr)	7440 - 47 - 3	10-27	1 (Cr)	0.5 (Cr)
*Nickel (Ni)	7440 - 02 - 0	0-22	1.0 (Ni)	1.0 (Ni)
*Manganese (Mn)	7439 - 96 - 5	0.3 - 10	5.0 (fume)	0.2 (Mn)
*Copper (Cu)	7441 - 50 - 8	0.04 - 4.4	0.1 (fume)	0.2 (fume)
Molybdenum (Mo)	7439 - 98 - 7	0-4	15 (total dust)	15 (insoluble)
**Aluminum (Al)	7429 - 90 - 5	0-2	15 (total dust)	5 (welding fumes)
Calcium (Ca)	7440 - 70 - 2	0-2		
Silicon (Si)	7440 - 21 - 3	0-2	15 (dust)	10 (inhaled dust)
Columbium (Cb)	7440 - 03 - 1	0-1		
*Cobalt (Co)	7440 - 48 - 4	0-0.75	0.1 (Co)	0.02 (Co)
Titanium (Ti)	7440 - 32 - 8	0-0.7	15 (TiO2)	10 (TiO2)
Carbon (C)	7440 - 44 - 0	0-0.5		
Sulphur (S)	7704 - 34 - 8	0.001 - 0.4	13 (SO2)	5 (SO2)
*Selenium (Se)	7482 - 49 - 2	0-0.4	0.2 (Se)	0.2 (Se)
*Phosphorus (P)	7723 - 14 - 0	0.001 - 0.2	0.1 (P)	0.1 (P)
Tantalum (Ta)	7440 - 25 - 7	0.01 - 1.1	5 (Ta)	5 (Ta)
Note: Those elements identified by an * (including red or white phosphorous) and those elements capable of generating highly toxic fumes or dusts (identified by a **) are classified as toxic by EPA in 40 CFR 372.65 and subject to reporting requirements of SARA Title III Section 313 and 40 CFR 372				
SECTION 3 - HAZARDS IDENTIFICATION				
EMERGENCY OVERVIEW: Stainless steel products in their solid state present no inhalation, ingestion or contact health hazard. However, inhaling dusts, fumes or mists which may be generated during certain manufacturing procedures (burning, melting, welding, sawing, brazing, grinding and machining) may be hazardous to your health. Dusts may also be irritating to the unprotected skin or eyes.				
ACUTE EFFECTS: Excessive exposure to dusts / fumes may cause irritation of eyes, nose or throat. Inhalation of dusts / fumes may result in metal fume fever (metallic taste in mouth, dryness and irritation of throat, chills and fever).				
CHRONIC EFFECTS: Prolonged inhalation of fumes or dusts may cause a variety of adverse health effects to the respiratory system, including (but not necessarily limited to) lesions of the mucous membrane, bronchitis, pneumonia and cancers to the nasal cavity and respiratory tract.				
POTENTIAL HEALTH EFFECTS/MEDICAL CONDITIONS AGGRAVATED BY EXPOSURE: Any pre-existing chronic respiratory condition (asthma, chronic bronchitis, emphysema).				
ROUTES OF ENTRY: Inhalation (dusts / fumes / mists), Contact with Skin and Eyes (dusts / mists), Ingestion (dusts).				
SECTION 4 - FIRST AID MEASURES				
INHALATION: Immediately remove victim to fresh air. If condition persists, consult physician.				
EYE CONTACT: Immediately flush with running water to remove particulates, consult physician.				
SKIN CONTACT: If irritation develops, remove clothing and wash with soap and water. If condition persists, consult physician.				
INGESTION: Consult physician.				
NOTE TO PHYSICIAN: None.				
SECTION 5 - FIRE FIGHTING MEASURES				
FLASH POINT: Nonflammable.				
SPECIAL FIRE FIGHTING INSTRUCTIONS AND EQUIPMENT: None required.				
AUTOIGNITION TEMPERATURE: NA				
FLAMMABLE LIMITS: Nonflammable.				
EXTINGUISHING MEDIA: Use what is appropriate for surrounding fire.				
HAZARDOUS COMBUSTION PRODUCTS: None.				
UNUSUAL FIRE AND EXPLOSION HAZARDS: None known.				
SECTION 6 - ACCIDENTAL RELEASE MEASURES				
CLEAN UP PROCEDURES: No special procedures needed.				
SPECIALIZED EQUIPMENT: None.				
SECTION 7 - HANDLING AND STORAGE				
PRECAUTIONS TO BE TAKEN IN HANDLING: Minimize activities which may generate dusts, mists or fumes. Keep areas well ventilated. Use suitable equipment to move materials.				
PRECAUTIONS TO BE TAKEN IN STORAGE: None required.				
SECTION 8 - EXPOSURE CONTROLS, PERSONAL PROTECTION				
RESPIRATORY PROTECTION: Wear NIOSH approved dust / mist / fume respirator when welding or burning this metal.				
EYE/FACE PROTECTION: Face shields (welding or burning), Safety glasses (cutting or grinding).				
OTHER PROTECTIVE EQUIPMENT: Use appropriate protective clothing such as welding aprons and gloves when welding or burning.				

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<div style="border: 1px solid black; width: 200px; height: 40px; margin: 0 auto;"></div>																													
1.	<p>Chemical product and company identification</p> <p>Article number: </p> <p>Product name: </p> <p>Company: </p>																												
2.	<p>Components, formulae</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Glycerol (C₃H₈O₃) 84.4% (w/v)</td> <td style="text-align: right;">CAS-Nr.: [56-81-5]</td> </tr> <tr> <td>Potassium chloride (KCl) 11.0%</td> <td style="text-align: right;">CAS-Nr.: [7447-40-7]</td> </tr> </table>	Glycerol (C ₃ H ₈ O ₃) 84.4% (w/v)	CAS-Nr.: [56-81-5]	Potassium chloride (KCl) 11.0%	CAS-Nr.: [7447-40-7]																								
Glycerol (C ₃ H ₈ O ₃) 84.4% (w/v)	CAS-Nr.: [56-81-5]																												
Potassium chloride (KCl) 11.0%	CAS-Nr.: [7447-40-7]																												
3.	<p>Possible hazards</p> <p>No hazardous material.</p>																												
4.	<p>First aid</p> <p>Eye contact: Flush with water.</p> <p>Skin contact: Flush with water.</p> <p>Ingestion: Seek medical advice.</p>																												
5.	<p>Procedures for fire</p> <p>Fire extinguishing media: appropriate: water, powder, foam, CO₂.</p> <p style="padding-left: 40px;">Inappropriate: not known.</p>																												
6.	<p>Procedures for unintended release</p> <p>Wipe up small amounts with absorbent material and dispose with regular waste.</p>																												
7.	<p>Storage and handling</p> <p>No special measures are necessary.</p>																												
8.	<p>Exposition limits and personal protective measures</p> <p>Respiratory protection: not required.</p> <p>Eye protection: wear face shield.</p> <p>Skin protection: wear rubber gloves.</p>																												
9.	<p>Physical and chemical properties</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">State of matter:</td> <td>liquid</td> </tr> <tr> <td>Colour:</td> <td>colourless</td> </tr> <tr> <td>Odor:</td> <td>odorless</td> </tr> <tr> <td>pH-value (25°C):</td> <td>7.0</td> </tr> <tr> <td>Melting point:</td> <td>ca. -30° C</td> </tr> <tr> <td>Boiling point:</td> <td>ca. 120° C</td> </tr> <tr> <td>Ignition point:</td> <td>not applicable</td> </tr> <tr> <td>Flash point:</td> <td>not flammable</td> </tr> <tr> <td>Explosion limits: lower</td> <td>not applicable</td> </tr> <tr> <td style="padding-left: 40px;">upper</td> <td>not applicable</td> </tr> <tr> <td>Density:</td> <td>1.19 g/cm³</td> </tr> <tr> <td>Solubility in ...:</td> <td>not determined</td> </tr> <tr> <td>Vapour pressure:</td> <td>not determined</td> </tr> <tr> <td>Viscosity:</td> <td>ca. 30 mPa s</td> </tr> </table>	State of matter:	liquid	Colour:	colourless	Odor:	odorless	pH-value (25°C):	7.0	Melting point:	ca. -30° C	Boiling point:	ca. 120° C	Ignition point:	not applicable	Flash point:	not flammable	Explosion limits: lower	not applicable	upper	not applicable	Density:	1.19 g/cm ³	Solubility in ...:	not determined	Vapour pressure:	not determined	Viscosity:	ca. 30 mPa s
State of matter:	liquid																												
Colour:	colourless																												
Odor:	odorless																												
pH-value (25°C):	7.0																												
Melting point:	ca. -30° C																												
Boiling point:	ca. 120° C																												
Ignition point:	not applicable																												
Flash point:	not flammable																												
Explosion limits: lower	not applicable																												
upper	not applicable																												
Density:	1.19 g/cm ³																												
Solubility in ...:	not determined																												
Vapour pressure:	not determined																												
Viscosity:	ca. 30 mPa s																												

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10. Conclusion/結論

Soxhlet extraction test was conducted to NC IOL in accordance with ISO11979-5:2006 Annex A. By comparing with the measured value of NY-60 before storage, the results were identified as follows:

- 1) The extract analysis^{*1)} of NC IOL identified MMA only, which is identical to that of NY-60 as comparison sample. MMA elution volume per lens was 0.37µg for NC IOL and 0.93µg for NY-60. MMA elution weight from NC IOL was less than half weight in comparison with NY-60. However, considering the elution weight identified in very small amount, it seems the same level of elution weight. The following table 1 shows the safety factor calculated by no-observed-effect concentration (NOEC) of the cytotoxicity test data^{*2)} and that identified the elution weight to be harmless.
- 2) For the constituents except for MMA, all indicated the elution weight below detection limit (N.D.). When compared the constituents with the solvent control under HPLC, no peak other than MMA peak was identified increasing or decreasing.
- 3) Extration rate calculated from sample weight difference before and after extraction was 0.28% and was equivalent^{*3)} to 0.25% of NY-60.


- * 1): For calibration curve used for quatitative analysis of constituent, refer to "Re-creation of the calibration curve with the change of HPLC detector lamp (R-DCI-20080919-1)". Refer to the attachment 1 for the analysis condition.
- * 2): The cytotoxicity test of methyl methacrylate using V79 cells (refer to HRI Report No. 09-I-444, Test Number: F-09-086, issued Nov. 18, 2009)
- * 3): The constant mass of the measured weight is defined as within ±0.1mg. Therefore, the difference of the above 0.03% was determined to be within the margin of measurement error.

As the results of Soxhlet extraction test above, the elution characteristic of NC IOL was identified equivalent to that of NY-60.

Table 1. MMA Safety factor

Test Sample	Detection concentration [ppm]	Elution weight in extract [µg]	Elution weight per lens [µg]	Elution weight/ lens/1mL [µg/mL]	no-observ ed-effect concentration (NOEC) [µg/mL]	Safety factor (NOEC/Elution weight per lens per 1mL)
NC IOL	1.47	14.70	0.368	0.0368	5000	Approx. 13600-fold
NY-60	3.70	37.00	0.925	0.0925	5000	Approx. 54000-fold

* No-observed-effect concentration (NOEC) in Table 1 in the cytotoxicity test was set to the maximal concentration of test system, because it is reported that the concentration of 5000µg/mL as maximal concentration of test system do not inhibit the colony forming.

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10. Conclusion

As a result of the test for leachable on NC IOL in accordance with ISO 11979-5 : 2006 Biocompatibility Annex B, the following was indicated. Furthermore, the measurement of spectrum transmittance of the lens was performed with the test method defined in ISO11979-2 : 1999/Cor 1:2003 and to be according to the Work Instruction (Control No.: WIDS/AC).

- 1) From the extraction with distilled water, no extraction of constituents was verified. (Below detection limit: N.D.)
- 2) From the extraction with acetone, MMA, PBMA, and AVBAA were verified to be extracted. However, the extraction volume was assured to be safety with the result of cytotoxicity test (NOEC: no-observed-effect concentration). (see Table 1) Furthermore, in comparison with the previous same model of yellow lens (NY-60) is given in the attachment 2.
- 3) Extraction rate calculated from sample weight difference before and after extraction was 0.01% for distilled water and 0.25% for acetone. It was equivalent to this of yellow lens (NY-60) having the same model. (0.00% for distilled water, 0.23% for acetone)
- 4) As per the following Fig. 1, the transmittance curve corresponded to the curves before and after extraction. Also, the transmittance value around transition point (400nm to 500nm) was determined to be within the margin of measurement error.

*: Extractable report from NH IOL (lens sample): R-SYT-20080710-1

Table 1. Safety factor calculated from the results of cytotoxicity tests for each substance.

Test Sample	Detection concentration [ppm]	Elution weight in extract [µg]	Elution weight per lens [µg]	Elution weight/ lens/mL [µg/mL]	no-observed-effect concentration (NOEC) [µg/mL]	Safety factor (NOEC/Elution weight per lens per 1mL)
MMA	21.07	299.19	4.214	0.2968	—	No-cytotoxicity
EDMA	<i>0.079</i>	<i>1.12</i>	<i>0.016</i>	<i>0.0011</i>	30	Approx. 27000-fold
n-BA	<i>0.073</i>	<i>7.04</i>	<i>0.015</i>	<i>0.0010</i>	25	Approx. 24000-fold
PBMA	0.81	11.50	0.162	0.0114	50	Approx. 4400-fold
BRM	<i>0.748</i>	<i>2.07</i>	<i>0.029</i>	<i>0.0021</i>	—	No-cytotoxicity
T-150	<i>0.707</i>	<i>7.43</i>	<i>0.020</i>	<i>0.0014</i>	—	No-cytotoxicity
AVBAA	0.137	1.95	0.027	0.0019	310	Approx. 16000-fold

Note: The italic types in Table 1 indicated the concentration to be below the detection limit and as such the detection limit was applied. (calculated from 14.2mL for volume and 71 pieces for lenses)


Fig. 1 Spectral transmittance before and after extraction (average of n=5)

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Supersedes Date: 07/11/2000

1. PRODUCT IDENTIFICATION

Product Name: PVC Sheet
Product Code: InteDur, InteFoam, InteClear, InteCel, InteCel PW, Tuf Board
Chemical Family: Polymer of Chlorinated Hydrocarbon
Chemical Name: Polyvinyl Chloride
CAS No: 9002-86-2
Synonyms: Rigid PVC Sheet (Type I and II), Expanded Foam PVC Sheet, Integral Skin Expanded Foam PVC Sheet, Clear PVC Sheet, Wood PVC sheet, Celuka PVC Sheet, Wood/PVC Composite Sheet
Formula: Proprietary
Technical Information:

2. PRODUCT INGREDIENTS

<u>No.</u>	<u>Components</u>	<u>CAS No.</u>	<u>Percent (%)</u>	<u>OSHA PEL</u>
1	PVC	9002-86-2	50 - 100%	5 mg/M ³ (respirable dust)
2	Proprietary	Mixtures	0 - 50%	Not established

3. PHYSICAL/CHEMICAL PROPERTIES

Physical Form: Solid Sheet
Color: Finished sheet with colors specified
Odor: Insignificant
Boiling Point: Not applicable
Melting Point: Not established
Freezing Point: Not applicable
Solubility in Water: None
Specific Gravity: 0.4 - 2.0 (water = 1)
Vapor Density: Not applicable (air = 1)
Evaporation Rate: None (Butyl Acetate = 1)
Vapor Pressure: Not applicable
% Volatile: None
pH: Not applicable

The physical data presented above are typical values and should not be construed as a specification.

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Page : 1 of 2

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8.0 Materials of Construction

Table 2A: Materials Matrix of Device Assessed

Device Name: WallFlex Esophageal Stent Delivery System – Fully and Partially covered

Component/ Material Name	Part Number	Patient Contact? (Yes/No)	Chemical Name/ Formulation/Grade	Processing aids*		% Weight of device
				Supplier	BSC	
Stent- Fully and Partially covered						
Stent	90114613-XX	Yes	See Below	UN	UN	4.2%
Consists of: 0.110" Light Oxide Nitinol Wire	90114559-01		Nitinol Wire: Nickel: 55.0-56.0 Wt% Titanium: Balance Carbon: ≤500ppm or 0.05Wt% Oxygen: ≤500ppm or 0.05Wt% Hydrogen: ≤50ppm or 0.005 Wt% Others: ≤0.85 Wt%			
Coating	MS0063-010- 013	Yes	Mix Silicone Solution	UN	UN	1.73%
Consists of: Inhibitor Sil/DisS 2MET-3BUT 2OL	MS0039-001		Inhibitor Sil/DisS 2MET-3BUT 2OL			
MED4820 Liquid Silicone Elastomer	MS0043-07		Silicone, Part A 4820 NG PERM			
Xylene	MS0012-011		Xylene			
Primer	MS0042-001	Yes	Primer, Nusil, CF8-135	UN	UN	Trace
Suture	90070568-02	Yes	Polyester: 95.6% PTFE: 4% D&C Green # 6: 0.4%	UN	UN	Trace
UV Adhesive	90124956-01	Yes	Dymax 1128A-M Adhesive	UN	UN	Trace
Delivery Device- Fully covered only						
9-fla reconstrainment band/Turlux 2812TR	90451399-01	Yes	Methyl Methacrylate Acrylonitrile/ Butadiene/ Styrene (MABS) (Molded at Protek)	UN	UN	0.02%
Delivery Device- Partially covered only						
Compression sleeve	90384025-01	Yes	Pebax 2533 SA01 (Molded at Protek)	UN	UN	0.02%

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Molecular Mass: 310.18

CAS Number: 7758-87-4

TEST	SPECIFICATION	LOT BCBD8691V RESULTS
APPEARANCE (COLOR)	WHITE/COLORLESS	
APPEARANCE (FORM)	POWDER	
PHASE-COMPOSITION (BY XRD ANALYSIS)	CORRESPONDS TO REQUIREMENTS	
TRACE ANALYSIS (ADDITIONAL TESTS)	TOTAL HEAVY METALS AS LEAD ≤ 50MG/KG	
METAL TRACE ANALYSIS (ICP)	CORRESPONDS TO REQUIREMENTS	
BARIUM (ICP)	≤ 20 MG/KG	
CHROMIUM (ICP)	≤ 15 MG/KG	
NICKEL (ICP)	≤ 5 MG/KG	
ZINC (ICP)	≤ 20 MG/KG	
ARSENIC TRACES (MHS-AAS)	≤ 0.5 MG/KG	
MERCURY TRACES (AFS)	≤ 500 UG/KG	
TOTAL SULFUR AS SO ₄ (ICP)	≤ 100 MG/KG	
CHLORIDE (CL)	≤ 50 MG/KG	
QC RELEASE DATE	14/FEB/11	

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Shelf-life Report
(accelerated, 1 and 2 year real time)
for

Date: September 2010

	Signature	Date
Prepared by 		
Approved by 		
Approved by 		
Approved by 		

Description of change
- Addition of 2 year real time shelf life data

CONTENTS

1. SUMMARY	3
2. OBJECTIVE	4
3. INTRODUCTION	4
4. DEFINITIONS	4
5. PROCEDURES	5
5.1 <i>Execution of the study</i>	5
5.2 <i>Test methods and acceptance criteria</i>	5
6. AMENDMENTS TO THE PROTOCOL	5
6.1 <i>Syringe performance</i>	5
6.2 <i>Extrusion force test method improvement</i>	6
6.3 <i>Packaging change</i>	6
7. RESULTS	6
8. CONCLUSIONS	7
9. RECOMMENDATIONS	7
10. REFERENCES	7
11. APPENDICES	7

Confidential

1. **SUMMARY**

Product lots of [REDACTED] have been tested according to the Shelf life stability protocol 1601-029.00.

Lots have been tested for chemical characteristics, label and packaging integrity, appearance, cytotoxicity and sterility.

Test results simulating a 1 year and a 2 year shelf life of [REDACTED] demonstrate that all tested lots meet the test specifications.

One and two year real time aging test results of [REDACTED] demonstrate that all tested lots meet the specifications.

Based on the accelerated and real time aging results it is concluded that the [REDACTED] [REDACTED] meets a shelf life claim of 2 years.

2. **OBJECTIVE**

The objective of this study is to demonstrate that [REDACTED] meets its expected shelf life claim of 2 years. To verify the expected shelf life claim of 2 years an accelerated and a real time shelf life study have been performed. The test results are compared to the product and packaging specifications defined in table 1.

3. **INTRODUCTION**

This report describes the results of the shelf life study of final product [REDACTED] as directed by the shelf life study protocol 1601-029.00 [REDACTED] lots have been used for testing chemical characteristics, label and packaging integrity, appearance and sterility. Tests have been performed internally and at external QC testing laboratories [REDACTED]. The accelerated shelf life study results corresponding with a shelf life of 1 and 2 years, including the 1 and 2 years real time shelf life study results are presented in this report.

4. **DEFINITIONS**

[REDACTED] Is a sterile, latex-free, non-pyrogenic, totally bioresorbable, non-permanent implant, whose principle component is synthetic Poly-ε-Caprolacton (PCL) microspheres suspended in a gel carrier of phosphate buffered saline (Ph.Eur.), glycerin (Ph.Eur.) and carboxymethylcellulose (USP). [REDACTED] has a particle size range of 25-50 microns and should be injected with a 27 gauge needle. [REDACTED] is available in different initial average molecular weights of the PCL microspheres. Currently the following product family is available; [REDACTED]

Initial packaging system:

Packaging: Consists of primary and secondary packaging

Primary packaging: Polypropylene syringe with closure cap. The primary packaging serves as a sterile barrier.

Secondary packaging: Validated single packaging combination of medical grade [REDACTED] blister sealed on an automatic thermoforming packaging machine, using validated seal time, seal temperature and seal pressure parameters.

Current packaging system (see also §6.3):

Packaging: Consists of primary, secondary and tertiary packaging

Primary packaging: Polypropylene syringe with closure cap. The primary packaging serves as a sterile barrier.

Secondary packaging: Protective laminated (aluminum) pouch packaging

Tertiary packaging: Protective cardboard box including the IFU

In all cases, the shelf life study has been performed without using the secondary and tertiary packaging.

Worst case product: [REDACTED] in a 1 ml syringe is considered as representative for the entire molecular weight range of [REDACTED] because the degradation mechanism and therefore degradation behavior is the same (predictable, controlled and total) for the entire product family. The only difference between the various molecular weights of the PCL microspheres within the product range is the initial molecular weight, meaning that the time needed for hydrolysis increases linearly with increased initial molecular weight. All other process steps and product characteristics are identical, thereby supporting the representative characteristics of [REDACTED]. Therefore, results are extrapolated to the product family.

Sterilization: Product components, microspheres and carrier, are routinely gamma irradiated at a dose of 25 kGy - 40 kGy, and steam sterilized (autoclave), respectively. The final product is aseptically processed. In order to avoid any sterility effects, batches are chosen randomly.

5. PROCEDURES

5.1 Execution of the study

The shelf life study is performed on [REDACTED] final product according to the shelf life study protocol 1601-029.00.

Accelerated shelf life study was performed at 40°C ± 1°C for minimally 71 and 140,5 days corresponding with a shelf life of 1 and 2 years at 15-25°C, respectively.

Real time shelf life study was performed at 15-25°C for 1 and 2 years.

5.2 Test methods and acceptance criteria

Table 1 shows a listing of performed test methods and the acceptance criteria.

Table 1.

Parameter	Test method	Specification
<i>Final Product</i>		
Inherent viscosity	Ubbelohde method ASTM D2857, ISO 1628-1	Conform
Syringe performance*	[REDACTED]	[REDACTED]
Sterility	EP	no growth
Cytotoxicity	ISO 10993-5	Acceptable
<i>Packaging</i>		
Legibility labels	Readable primary packaging	Label clearly legible from a distance of 30 ± 5 cm No spots Non peelable No air bubbles and loose ends
Seal integrity	Burst test, ISO 11067-1, ASTM F1140	≥ 17,0 kPa*** Visually OK
Appearance	Visual inspection	Product specification

* Syringe performance by means of extrusion tests will be performed on accelerated aged syringes "for information only" since the temperature raise during aging could influence the force needed for extrusion of the product.

** Specification has changed due to test method improvement, for explanation see §6.2.

*** Seal integrity testing is omitted at 1 and 2 year real time aging, seal integrity specifications for the blister are not critical as the blister does not serve as microbial barrier and only functions as a protective packaging for the syringe. Seal integrity was verified by opening the blister by hand and visual inspection of the seal after opening.

6. AMENDMENTS TO THE PROTOCOL

6.1 Syringe performance

The acceptance criteria for syringe performance testing has been adjusted compared to the shelf life protocol 1601-029.00. Initially, it was believed that the extrusion force using a 27G ¾ needle should not exceed the 50 N force to be applied on the syringe using the test set up at [REDACTED]. However, practical extrusion force assessments using [REDACTED]

performed at the cosmetic clinic who is executing the clinical study revealed that higher extrusion forces were acceptable (< 70 N). Even when smaller needles (e.g. 30G) were used still acceptable extrusion forces could be applied on the syringes using the . Furthermore it became clear that there is a direct relation between the time needed for pushing the content out of the syringe and the force needed to do so. It was therefore decided to recalculate the maximum allowable extrusion force when using a 27G ¼ needle to ≤ 7 N per second. For 26G ½ needles an extrusion force of < 5 N per second is acceptable.

6.2 Extrusion force test method improvement

Extrusion force testing of 1 year real time aged syringes were performed using a different apparatus when compared to the 1 and 2 year accelerated aged syringes. For the accelerated aged syringes an was used whereas for the 1 year real time aged syringes a Sintech 200/M tensile compression tester was used. Extrusion force measured with the and the extrusion force measured with the tensile/compression tester 200/M have found to be comparable (see studyplan R08017). It has been concluded in document R08017 that because the experiments with the tensile/compression tester Sintech 200/M are more accurate, standardized and give a consistent pattern of the force needed to start and maintain the flow of dermal filler, the extrusion force of products will be tested with the tensile/compression tester Sintech 200/M with a extrusion rate of 1ml/min. The specification is adjusted into <40 N for a 27G ¼ needle with a flow rate of 1 ml/min. (These forces are far below the maximum lateral finger/thumb forces for men and women which are 103N and 69N, respectively (www.fda.gov/ohrms/dockets/ac/03/briefing/03034b1_02_precl-fda.doc)).

6.3 Packaging change

The packaging concept of has changed during the shelf life study. The initial packaging concept comprised a polypropylene syringe serving as microbial barrier (primary packaging) packed in a packaging (secondary packaging) serving as a protective barrier during transport and storage. The new packaging concept comprises the identical primary packaging (i.e. polypropylene syringe) and function, packed in a protective laminated (aluminum) pouch packaging (thickness 119 micron), which is sealed and packed in a protective cardboard box (thickness 1,2 mm). Since the primary packaging and function does not change and the new secondary and tertiary packaging are believed to have equal or even better protective properties when compared to the blister, it is concluded that the results of this shelf life study are also valid for the packed product when packed according to the new packaging concept.

7. RESULTS

The results of the accelerated shelf life study on are described in Appendix 1-3.

Results are analyzed at T=71 and T=140,5 days corresponding with a shelf life simulating of 1 and 2 years respectively. Results are compared to the specifications initially set and analyzed at T=0 days.

Syringe performance testing by means of extrusion force testing was performed "for information only" since the influence of higher temperatures during accelerated aging (40°C) could not be foreseen.

Cytotoxicity testing of 71 days aged devices was omitted since the cytotoxicity test performed after 140,5 days aging did not demonstrate any cytotoxic effects.

The results of the 1 year real time aging study on are described in Appendix 4. Results demonstrate that all tested lots meet the requirements.

The results of the 2 year real time aging study on are described in Appendix 5. Results demonstrate that all tested lots meet the requirements.

Appendix 6 provides an overview of all tested lots (accelerated and real time aging) in time. It is shown that there is a good correlation between accelerated aging and real time aging results for Ellansé™.

8. CONCLUSIONS

Based on the results demonstrated in Appendix 1-6 it can be concluded that all tested lots meet the specifications and pass the tests performed. Based on data obtained from accelerated and real time aging it is claimed that [REDACTED] meets a shelf life of 2 years. Furthermore, it is concluded that the results of this shelf life study are also valid for the [REDACTED] packed product when packed according to the new packaging concept.

9. RECOMMENDATIONS

Changes in the production and/or packaging process or concept will be evaluated for impact on the shelf life results and could lead to a revalidation of the product shelf life.

10. REFERENCES

1. 1601-029.00 Shelf-life testing protocol (Accelerated & real-time) for [REDACTED]
2. Studyplan 08005 Extrusion force testing
3. Studyplan 08017 Extrusion force testing of [REDACTED] with Sintech 200/M
4. EP
5. EN 868-5:1999 Packaging materials and systems for medical devices which are to be sterilized - Part 5: Heat and self-sealable pouches and reels of paper and plastic film construction - Requirements and test methods
6. NEN EN ISO 11607-1&2:2006 Packaging for terminally sterilized medical devices part 1 and part 2
7. ASTM F1140-00 (Reapproved 2005) Standard Test Methods for Internal Pressurization Failure Resistance of Unrestrained Packages for Medical Applications
8. ASTM D2857-95 Standard Practice for Dilute Solution Viscosity of Polymers
9. EN-ISO 1628-1:1998 Plastics - Determination of the viscosity of polymers in dilute solution using capillary viscometers - Part 1: General principles
10. ISO 10993-5: 1999 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
11. ASTM F1980-02 Standard Guide for Accelerated Aging of Sterile Medical Device Packages
12. Updated 510(k) Sterility Review Guidance K90-1; Guidance for Industry and FDA, August 30, 2002
13. Guideline for Submitting Documentation for the Stability of Drugs and Biologics, FDA Center for Drug Evaluation and Research, February 1987
14. Other applicable standards and guidelines

11. APPENDICES

- Appendix 1: Acc. aging shelf-life test results of [REDACTED] T=0 days
Appendix 2: Acc. aging shelf-life test results of [REDACTED] T=71 days
Appendix 3: Acc. aging shelf-life test results of [REDACTED] T=140,5 days
Appendix 4: Real time aging shelf-life test results of [REDACTED] T=1 year
Appendix 5: Real time aging shelf-life test results of [REDACTED] T =2 years
Appendix 6: Overview accelerated and real time aging results [REDACTED]

Appendix 1 ACCELERATED AGING SHELF LIFE TEST RESULTS T = 0 days

Analytical tests			Evaluation time												
			0 days												
Parameter	Test method	Specification	LOT 08004-S3	LOT 08004-S3	LOT 08004-S3	AVG	LOT 08004-S4	LOT 08004-S4	LOT 08004-S4	AVG	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	AVG	AVG Test
Inherent viscosity	Ubbelohde method ASTM D2457, ISO 1129-1	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Pass
Syringe performance	Extrusion force test (EMCM BV)	<5 Ns (25G 1/4 needle) <7 Ns (27G 1/4 needle)	3.6 Ns (26G 1/4)	3.6 Ns (26G 1/4)	3.2 Ns (25G 1/4)	3.7 Ns (26G 1/4)	5.5 Ns (27G 1/4)	5.5 Ns (27G 1/4)	5.7 Ns (27G 1/4)	5.6 Ns (27G 1/4)	5.2 Ns (27G 1/4)	5.5 Ns (27G 1/4)	5.7 Ns (27G 1/4)	5.5 Ns (27G 1/4)	Pass
Sterility	EP	No growth	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Pass
Cytotoxicity	ISO10993-5	Acceptable	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass
Legibility labels	Readable primary packaging	Label clearly legible from a distance of 30 ± 5 cm No spots Non peelable No air bubbles and loose ends	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass
Seal integrity	Burst test, ISO 11067-1, ASTM F1140	≥ 17.0 kPa	36.6 kPa (avg. start and end of seal process)			36.6 kPa	35.3 kPa (avg. start and end of seal process)			36.6 kPa	36.6 kPa (avg. start and end of seal process)			36.6 kPa	Pass
Appearance	Visual inspection	Product specification	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass

Appendix 2 ACCELERATED AGING SHELF LIFE TEST RESULTS T = 71 days

Analytical tests			Evaluation time												
			71 days												
Parameter	Test method	Specification	LOT 08004-S3	LOT 08004-S3	LOT 08004-S3	AVG	LOT 08004-S4	LOT 08004-S4	LOT 08004-S4	AVG	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	AVG	AVG Test
Inherent viscosity	Ubbelohde method ASTM D2457, ISO 1625-1	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Pass
Syringe performance	Extrusion force test (EMCM BV)	<5 Ns (26G 1/4 needle) <7 Ns (27G 1/4 needle)	6.0 Ns (27G 1/4)	6.0 Ns (27G 1/4)	5.9 Ns (27G 1/4)	6.0 Ns (27G 1/4)	6.5 Ns (27G 1/4)	7.0 Ns (27G 1/4)	6.7 Ns (27G 1/4)	6.7 Ns (27G 1/4)	7.0 Ns (27G 1/4)	6.3 Ns (27G 1/4)	6.3 Ns (27G 1/4)	6.7 Ns (27G 1/4)	Pass
Sterility	EP	No growth	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Pass
Cytotoxicity	ISO10993-5	Acceptable	Omitted, see results 2 year			N/A	Omitted, see results 2 year			N/A	Omitted, see results 2 year			N/A	N/A
Legibility labels	Readable primary packaging	Label clearly legible from a distance of 30 ± 5 cm No spots Non peelable No air bubbles and loose ends	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass
Seal integrity	Burst test, ISO 11067-1, ASTM F1140	≥ 17.0 kPa	38.4 kPa	38.1 kPa	36.1 kPa	38.2 kPa	40.3 kPa	38.6 kPa	37.7 kPa	38.8 kPa	35.0 kPa	38.3 kPa	35.2 kPa	37.6 kPa	Pass
Appearance	Visual inspection	Product specification	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass

Appendix 3 ACCELERATED AGING SHELF LIFE TEST RESULTS T = 140.5 days

Analytical tests			Evaluation time: 140.5 days												
Parameter	Test method	Specification	LOT 08004-S3	LOT 08004-S3	LOT 08004-S3	AVG	LOT 08004-S4	LOT 08004-S4	LOT 08004-S4	AVG	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	AVG	AVG T=140.5
Inherent viscosity	Ubbelohde method ASTM D2897, ISO 1628-1	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Pass
Syringe performance	Extrusion force test (EMCM BV)	<5 Ns (26% needle) <7 Ns (27% needle)	6.0 Ns (27% needle)	7.0 Ns (27% needle)	6.4 Ns (27% needle)	5.6 Ns (27% needle)	6.4 Ns (27% needle)	6.5 Ns (27% needle)	6.9 Ns (27% needle)	6.6 Ns (27% needle)	Blockage*	6.7 Ns (27% needle)	6.4 Ns (27% needle)	5.6 Ns (27% needle)	Pass
			LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5
Sterility	EP	No growth	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Pass
			LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5
Cytotoxicity	ISO10993-5	Acceptable	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass
			LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5
Legibility labels	Readable primary packaging	Label clearly legible from a distance of 30 ± 5 cm No spots Non peelable No air bubbles and loose ends	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass
Seal integrity	Burst test, ISO 11067-1, ASTM F1140	> 17.0 kPa	40.9 kPa	36.3 kPa	36.2 kPa	39.1 kPa	37.3 kPa	36.7 kPa	36.8 kPa	37.6 kPa	38.7 kPa	34.6 kPa	39.2 kPa	37.5 kPa	Pass
Appearance	Visual inspection	Product specification	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass

Appendix 4 REAL TIME AGING SHELF LIFE TEST RESULTS T = 1 YEAR

Analytical tests			Evaluation time: 1 year												
Parameter	Test method	Specification	LOT 08004-S3	LOT 08004-S3	LOT 08004-S3	AVG	LOT 08004-S4	LOT 08004-S4	LOT 08004-S4	AVG	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	AVG	AVG T=1 YR
Inherent viscosity	Ubbelohde method ASTM D2897, ISO 1628-1	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Pass
Syringe performance	Extrusion force test (EMCM BV, in house)	<5 Ns (26% needle) <7 Ns (27% needle) <40 N, flow rate 1ml/min (27% needle)	31.7 N (27% needle)	31.3 N (27% needle)	30.4 N (27% needle)	31.1 N (27% needle)	31.1 N (27% needle)	31.0 N (27% needle)	30.1 N (27% needle)	30.7 N (27% needle)	29.9 N (27% needle)	32.0 N (27% needle)	32.2 N (27% needle)	31.4 N (27% needle)	Pass
			LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5
Sterility	EP	No growth	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Sterile	Pass
			LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5
Cytotoxicity	ISO10993-5	Acceptable	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass
			LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5	LOT 08004-S5
Legibility labels	Readable primary packaging	Label clearly legible from a distance of 30 ± 5 cm No spots Non peelable No air bubbles and loose ends	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass
Seal integrity	Burst test, ISO 11067-1, ASTM F1140	> 17.0 kPa Visually OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	Pass
Appearance	Visual inspection	Product specification	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass

* Needle blockage

Appendix 5 REAL TIME AGING SHELF LIFE TEST RESULTS [REDACTED] T = 2 YEARS

Analytical tests			Evaluation time												
			2 years												
Parameter	Test method	Specification	LOI 0904-S3	LOI 0904-S3	LOI 0904-S3	AVG	LOI 8804-S4	LOI 8804-S4	LOI 8804-S4	AVG	LOI 0904-S6	LOI 0904-S6	LOI 0904-S6	AVG	AVG 142 YR
Inherent viscosity	Ubbelohde method ASTM D2857, ISO 1628-1	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Conform	Pass
Syringe performance	Extrusion force test (EMCM BV, in house)	<5 N/s (26G ½ needle) >7 N/s (27G ½ needle) <40 N, flow rate 1ml/min (27G ½ needle)	LOI 0904-S3	LOI 0904-S3	LOI 0904-S3	21,0 N (27G ½)	LOI 8804-S4	LOI 8804-S4	LOI 8804-S4	21,4 N (27G ½)	LOI 0904-S6	LOI 0904-S6	LOI 0904-S6	21,5 N (27G ½)	Pass
			31,1 N (27G ½)	30,4 N (27G ½)	31,4 N (27G ½)		31,6 N (27G ½)	31,4 N (27G ½)	31,2 N (27G ½)		31,6 N (27G ½)	32,0 N (27G ½)	31,7 N (27G ½)		
Sterility	EP	No growth	LOT 0904-S5	LOT 0904-S5	LOT 0904-S5	Sterile	LOT 8804-S4	LOT 8804-S4	LOT 8804-S4	Sterile	LOT 0904-S6	LOT 0904-S6	LOT 0904-S6	Sterile	Pass
			Sterile	Sterile	Sterile		Sterile	Sterile	Sterile		Sterile				
Cytotoxicity	ISO10993-5	Acceptable	LOI 0904-S5	LOI 0904-S5	LOI 0904-S5	Pass	LOI 8804-S4	LOI 8804-S4	LOI 8804-S4	Pass	LOI 0904-S6	LOI 0904-S6	LOI 0904-S6	Pass	Pass
			Pass	Pass	Pass		Pass	Pass	Pass						
Legibility labels	Readable primary packaging	Label clearly legible from a distance of 30 ± 5 cm No spots Non peelable No air bubbles and loose ends	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass
Seal integrity	Visual inspection	Visually OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	Pass
Appearance	Visual inspection	Product specification	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass

Appendix 6 OVERVIEW ACCELERATED AND REAL TIME AGING SHELF LIFE TEST RESULTS [REDACTED]

Analytical tests			Evaluation time				
Parameter	Test method	Specification	AVG T=0 DAYS	AVG T=71 DAYS	AVG T=140,5 DAYS	AVG T=1 YEAR	AVG T=2 YEARS
Inherent viscosity	Ubbelohde method ASTM D2857, ISO 1628-1	Conform	Conform	Conform	Conform	Conform	Conform
Syringe performance	Extrusion force test (EMCM BV)	<5 N/s (26G ½ needle) ≤7 N/s (27G ½ needle)	5,5 N/s (27G ½)	6,5 N/s (27G ½)	6,6 N/s (27G ½)	N/A	N/A
Syringe performance	Extrusion force test (in house)	<40 N, flow rate 1ml/min (27G ½ needle)	N/A	N/A	N/A	31,1 N (27G ½)	31,4 N (27G ½)
Sterility	EP	No growth	Sterile	Sterile	Sterile	Sterile	Sterile
Cytotoxicity	ISO10993-5	Acceptable	Pass	Pass	Pass	Pass	Pass
Legibility labels	Readable primary packaging	Label clearly legible from a distance of 30 ± 5 cm No spots Non peelable No air bubbles and loose ends	Pass	Pass	Pass	Pass	Pass
Seal integrity	Burst test, ISO 11067-1, ASTM F1140	≥ 17,0 kPa	36,6 kPa	38,2 kPa	38,1 kPa	N/A	N/A
Seal integrity	Visual inspection	Visually OK	N/A	N/A	N/A	Pass	Pass
Appearance	Visual inspection	Product specification	Pass	Pass	Pass	Pass	Pass

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[Redacted]

Test Report

Title: 18-Month Real Time Stability Report for 0.8cc [Redacted]

Written by: [Redacted] Date: 07/07/09

Approvals:

<u>Function</u>	<u>Signature</u>	<u>Date</u>
Eng/Development	[Redacted]	07/07/09
Regulatory Affairs	[Redacted]	07/07/09
Quality Control	[Redacted]	07/07/09
Manufacturing	[Redacted]	07/07/09
Quality Assurance	[Redacted]	07/07/09

1.0 PURPOSE

This report documents the results of the initial eighteen months stability of [REDACTED] product stored at 2-8°C conducted under [REDACTED] Stability Protocol for [REDACTED].

2.0 BACKGROUND

- 2.1 The [REDACTED] product is filled into a sterile polymer syringe with closure. A tamper-evident seal is placed on the syringe-closure tip. The [REDACTED] syringes are placed into a thermoformed tray and the tray heat-sealed using a foil lid. The trayed product is then placed into a carton and labeled. The [REDACTED] product is stored and transported at 2-8°C.
- 2.2 The Stability Protocol [REDACTED] applies to [REDACTED] product filled in [REDACTED] [REDACTED] challenged at three temperature ranges for up to three years in duration. A single real-time and two accelerated temperature challenges are being executed for the study.
- 2.3 This report only addresses the stability results at 2-8°C since the 2-8°C stability results are the only stability data to support actual [REDACTED] product stability. The accelerated stability studies at 25°C and 38°C are used for information purpose only and the results are documented in a separate report.

3.0 SCOPE

This report is limited to the analysis of the first eighteen calendar months stability for [REDACTED] product stored at 2-8°C conducted under Stability Protocol [REDACTED]. Twenty-four and thirty six month stability time points will be addressed in future reports.

4.0 REFERENCES

- 4.1 [REDACTED] *Release Testing - Final Product*
- 4.2 [REDACTED] *Stability Protocol for [REDACTED]*

5.0 DEVIATIONS/JUSTIFICATIONS

- 5.1 The syringe extrusion force of [REDACTED] products, as reported in the "Viscometry Measurement Report" of [REDACTED] was measured at a constant crosshead speed up to a maximum of 12.7 mm/min. However, the current syringe extrusion force measurement is performed at a constant speed of 100 mm/min as specified in [REDACTED].

The extrusion force limit of 10 lbs (equivalent to 44.5 N) was determined in LRE 05-2 "Syringe Extrusion Force Feasibility Report". The extrusion force is a function of crosshead speed as shown in [REDACTED]. With the crosshead speed increase from 12.7 to 100 mm/min as to standardize the measurement condition, the extrusion force limit increases from 10 lbs (or 44.5 N) to 16.9 (or 75.3 N) using the ratio of 34.92/20.63 in Table 2 of LRE 05-2 (Instron Calc at 12.7 mm/min vs. 100 mm/min).

The specification of 44.5 N in [REDACTED] was an error and the revised extrusion force limit of 75.3 N as discussed above was used instead in this report.

- 5.2 A void at the cap end of the syringe was incorrectly classified as an "air bubble". As a result, it caused the false sense not meeting an acceptance criterion. Therefore the definition is revised for "air bubble" to be applied only to "air space that is encapsulated within the product" and will be reflected in [REDACTED] and/or related document. In this stability study, no interstitial air bubbles were noted at any time point.

A void at the cap end of the syringe is cosmetic in nature and not related to chemical stability of the product. In fact, this type of voids at the cap end is common in the injectables and the air in the voids (including the void in the needle) is extruded before product injection as general medical practices.

A void at the stopper tip end of the syringe is classified as a "fill volume failure." This void at the stopper tip end was overcome with overfill as defined in [REDACTED]

6.0 MATERIALS

- 6.1 [REDACTED] Packaging [REDACTED] Product Syringes into Product Cartons Lot F061000
- 6.1.1 [REDACTED] Lot S061003, fill date 02/13/06
- 6.2 [REDACTED] into Product Cartons Lot F061001
- 6.2.1 [REDACTED] fill date 02/10/06
- 6.3 [REDACTED] into Product Cartons Lot F061002
- 6.3.1 [REDACTED] Lot S061004, fill date 02/14/06

7.0 REQUIREMENTS AND INITIAL RESULTS

- 7.1 All the stability samples specified above were manufactured in accordance with [REDACTED] and [REDACTED]
- 7.2 Stability Test Requirements and Initial Release Testing Results: The stability test requirements and acceptance criteria are outlined in [REDACTED] Appendix II and listed in Table 1 below. The initial release results (time point = 0 months) for the three product lots in the stability study are shown in Table 2. The initial release test data are stored with the batch records. All results at the initial release passed the acceptance criteria.

Table 1: Product Release Tests and Acceptance Criteria

Test No.	Test Name	Specification
[REDACTED]	Visual Appearance	Pass
[REDACTED]	PMMA Homogeneity	Pass
[REDACTED]	Syringe Extrusion Force	< 75.3N
[REDACTED]	Syringe Mass Determination	For Information Only
[REDACTED]	Hydroxyproline	0.36 - 0.48%
[REDACTED]	Collagen Content	3.0 - 4.0%
[REDACTED]	SDS-PAGE	Banding pattern conforms to standard
[REDACTED]	Lidocaine HCl	0.27 - 0.33%
[REDACTED]	pH	6.8 - 7.8
[REDACTED]	Heavy Metals	< 20 ppm
[REDACTED]	Endotoxins	< 10 EU/ml
[REDACTED]	Sterility	Pass
[REDACTED]	Packaging Inspection	Pass
[REDACTED]	ArteFill Packaging Leak/Burst Test Method	≤0.02, ≥0.69 PSI

* These assays are not product release attributes

Table 2: Results (time point = 0 month) for Stability Lots

Lot Number	Visual Inspection / PMMA Homogeneity	Hydroxyproline Content	Collagen Content	Lidocaine Content	pH at 40C	Heavy Metals	Endotoxin	Sterility	SDS-PAGE	Syringe Extrusion Force
	Pass / Pass	0.36 - 0.48%	3.0 - 4.0%	0.27 - 0.33%	6.8 - 7.8	< 20 ppm	< 10 EU/ml	Pass	Conforms to Reference	< 75.3 N
S061002	Pass / Pass	0.42	3.5	0.32	6.9	ND	< 5	Pass	Conforms	37.2
S061003	Pass / Pass	0.42	3.5	0.32	6.9	ND	< 5	Pass	Conforms	33.8
S061004	Pass / Pass	0.42	3.5	0.32	6.9	0.8	< 5	Pass	Conforms	36.6

8.0 Real Time Stability Data at Actual Product Storage Condition – 2-8°C

8.1 Stability sample testing was performed according to the procedures listed on Appendix III. Test results for Day 90, Day 182, Day 281, Day 366, and Day 553 at 2-8°C are summarized in Table 3 to Table 7, respectively. Additionally, Day 366 and Day 744 samples from all three lots were tested and passed for sterility.

Table 3: Day 90 (3 months) at 2-8°C Testing Result Data for Stability Lots

Lot Number	Visual Inspection / PMMA Homogeneity	Hydroxyproline Content	Collagen Content	Lidocaine Content	pH at 40C	Heavy Metals	Endotoxin	SDS-PAGE	Syringe Extrusion Force	Leak Test	Burst Test	Syringe Mass Determination (mg)	Packaging Inspection
	Pass / Pass	0.36 - 0.48%	3.0 - 4.0%	0.27 - 0.33%	6.8 - 7.8	< 20 ppm	< 10 EU/ml	Conforms to Reference	< 75.3 N	≤ 0.02 PSI	≥ 0.69 PSI	-	Pass
S061002	Pass / Pass	0.41	3.4	0.30	6.9	ND	<5	Conforms	41.6	0.01	1.36	17.1	Pass
S061003	Pass / Pass	0.42	3.5	0.31	6.9	ND	<5	Conforms	39.1	0.01	1.38	17.7	Pass
S061004	Pass / Pass	0.42	3.5	0.30	6.9	ND	<5	Conforms	39.8	0.00	1.33	20.2	Pass

ND = Not Detected

Table 4: Day 182 (6 months) at 2-8°C Testing Result Data for Stability Lots

Lot Number	Visual Inspection / PMMA Homogeneity	Hydroxyproline Content	Collagen Content	Lidocaine Content	pH at 40C	Heavy Metals	Endotoxin	SDS-PAGE	Syringe Extrusion Force	Leak Test	Burst Test	Syringe Mass Determination (mg)	Packaging Inspection
	Pass / Pass	0.36 - 0.48%	3.0 - 4.0%	0.27 - 0.33%	6.8 - 7.8	< 20 ppm	< 10 EU/ml	Conforms to Reference	< 75.3 N	≤ 0.02 PSI	≥ 0.69 PSI	-	Pass
S061002	Pass / Pass	0.41	3.4	0.31	6.9	ND	<5	Conforms	44.1	0.01	1.41	24.8	Pass
S061003	Pass / Pass	0.41	3.4	0.30	6.9	2.1	<5	Conforms	41.6	0.00	1.42	27.7	Pass
S061004	Pass / Pass	0.42	3.5	0.30	6.8	ND	<5	Conforms	42.3	0.00	1.39	34.9	Pass

ND = Not Detected

Table 5: Day 281 (9 months) at 2-8°C Testing Result Data for Stability Lots

Lot Number	Visual Inspection / PMMA Homogeneity	Hydroxyproline Content	Collagen Content	Lidocaine Content	pH at 40C	Heavy Metals	Endotoxin	SDS-PAGE	Syringe Extrusion Force	Leak Test	Burst Test	Syringe Mass Determination (mg)	Packaging Inspection
	Pass / Pass	0.36 - 0.48%	3.0 - 4.0%	0.27 - 0.33%	6.8 - 7.8	< 20 ppm	< 10 EU/ml	Conforms to Reference	< 75.3 N	≤ 0.02 PSI	≥ 0.69 PSI	-	Pass
S061002	Pass / Pass	0.43	3.6	0.27	6.9	ND	-	Conforms	43.3	0.01	1.38	44.3	Pass
S061003	Pass / Pass	0.42	3.5	0.28	6.9	ND	-	Conforms	38.0	0.01	1.40	37.3	Pass
S061004	Pass / Pass	0.44	3.7	0.29	6.9	ND	-	Conforms	47.5	0.01	1.30	59.9	Pass

ND = Not Detected

Table 6: Day 366 (12 months) at 2-8°C Testing Result Data for Stability Lots

Lot Number	Visual Inspection / PMMA Homogeneity	Hydroxyproline Content	Collagen Content	Lidocaine Content	pH at 40C	Heavy Metals	Endotoxin	SDS-PAGE	Syringe Extrusion Force	Leak Test	Burst Test	Syringe Mass Determination (mg)	Packaging Inspection
	Pass / Pass	0.36 - 0.48%	3.0 - 4.0%	0.27 - 0.33%	6.8 - 7.8	< 20 ppm	< 10 EU/ml	Conforms to Reference	< 75.3 N	≤ 0.02 PSI	≥ 0.69 PSI	-	Pass
S061002	Pass / Pass	0.45	3.7	0.32	6.9	ND	<5	Conforms	40.8	0.00	1.40	38.6	Pass
S061003	Pass / Fail	0.45	3.7	0.33	6.9	ND	<8.8	Conforms	44.9	0.01	1.40	40.4	Pass
S061004	Pass / Pass	0.47	3.9	0.30	6.9	ND	<5	Conforms	42.4	0.00	1.28	50.1	Pass

ND = Not Detected

Table 7: Day 553 (18 months) at 2-8°C Testing Result Data for Stability Lots

Lot Number	Visual Inspection / PMMA Homogeneity	Hydroxyproline Content	Collagen Content	Lidocaine Content	pH at 40C	Heavy Metals	Endotoxin	SDS-PAGE	Syringe Extrusion Force	Leak Test	Burst Test	Syringe Mass Determination (mg)	Packaging Inspection
	Pass / Pass	0.36 - 0.48%	3.0 - 4.0%	0.27 - 0.33%	6.8 - 7.8	< 20 ppm	< 10 EU/ml	Conforms to Reference	< 75.3 N	≤ 0.02 PSI	≥ 0.69 PSI	-	Pass
S061002	Pass / Pass	0.40	3.3	0.29	6.9	0.2	<5	Conforms	42.7	0.00	1.33	58.9	Pass
S061003	Pass / Pass	0.40	3.4	0.29	6.9	ND	<5	Conforms	42.0	0.00	1.38	51.2	Pass
S061004	Pass / Pass	0.41	3.4	0.29	6.9	ND	<5	Conforms	40.3	0.00	1.28	72.3	Pass

ND = Not Detected

8.2 All lots passed for PMMA homogeneity at all time points with the exception of one syringe of Lot S061003 at Day 366. This failure represents one out of two hundred twenty five syringes inspected during the course of the stability study to date, and one out of forty five syringes at Day 366 time point at 2-8°C. This single failure result at 2-8°C is a suspected aberration as there is no apparent trend.

8.3 All three lots passed Lidocaine Hydrochloride content and Hydroxyproline content up to 18 months.

8.4 All three lots passed the revised extrusion force acceptance criterion specified in section 5.1 at all time points tested.

9.0 DISCUSSIONS

9.1 The twelve and twenty-four month sterility test results demonstrate that the container closure system is adequate to assure product sterility and integrity at the product storage condition.

9.2 At 2-8°C, chemical stability indicating parameters pH and protein purity by [REDACTED] have passed the acceptance criteria up to 18 months.

9.3 All three 0.8cc fill volume lots passed the acceptance criteria for lidocaine hydrochloride and hydroxyproline contents up to the eighteen month time point.

- 9.4 There was one PMMA homogeneity visual inspection failure at the twelve month time point. This result is suspected an aberration as there is no apparent trend.
- 9.5 All three lots passed [REDACTED] up to 18 months.
- 9.6 All other parameters other than stated above in this section (Heavy Metals, Endotoxin, Leak Test, Burst Test and Packaging Inspection) passed the acceptance criteria at 2-8°C up to Day 553 (18 months).

10.0 CONCLUSION

- 10.1 The data generated in this study support up to 18 months expiry dating for [REDACTED] fill volume [REDACTED] product.

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1.0 PURPOSE / OBJECTIVE

The purpose of this paper is to document a rationale for the aging study of [redacted] to justify that aging of [redacted] does not require interim (additional) time points to support labelling of the [redacted] for a 2 year shelf life.

2.0 SCOPE:

This position paper applies to all sizes of [redacted] s which would be processed through the standard [redacted] process for bioprosthesis product lines.

Table 1.0 – [redacted] in the scope of this paper

[redacted] Model Numbers Affected	
[redacted]	[redacted]
[redacted]	[redacted]
[redacted]	[redacted]
[redacted]	[redacted]

3.0 REFERENCE DOCUMENTS

- [redacted] - Two Year Real Time Shelf Life for [redacted] and [redacted]
- [redacted] - Two Year Real Time Shelf Life Report of [redacted], Model [redacted]
- [redacted] - Three Year Real Time Shelf Life Protocol (US Study) For [redacted] Revalving System [redacted] of the [redacted]
- [redacted] - Two Year Real Time Shelf Life Qualification for [redacted]
- [redacted] - Interim Report of 6 month shelf life study for [redacted] mm Model [redacted]
- [redacted] - Interim Report of 1 year shelf life study for [redacted] mm Model [redacted]

4.0 DEVICE DESCRIPTION

4.1 The [redacted] is a three component medical device (Figure 1.0) which is designed to replace a native aortic heart valve concomitant surgical removal of the failed native valve. It is composed of a [redacted], an [redacted] pericardium, [redacted], and [redacted] fixed [redacted].

Figure 1.0

5.0 RATIONALE

5.1 [redacted] Configurations

The [redacted] are approved in two configurations. The originally approved configuration of the [redacted] was without [redacted] anti-mineralization treatment. The more recently approved [redacted] configuration includes the [redacted] anti-mineralization treatment.

The [redacted] with [redacted] Treatment is supported by shelf life studies performed for other approved [redacted] tissue valves [redacted] since the [redacted] tissue and [redacted] treatment processing used for [redacted] including the packaging solution and chemical sterilization parameters, are identical to the approved [redacted] tissue heart valves. These approved [redacted] tissue valves have demonstrated that the [redacted] treatment did not have an effect on the product performance of the [redacted] tissue valves aged up to 5 years. It follows therefore that because the material and its interaction with the [redacted] treatment are comparable, that the data generated for [redacted] and [redacted] Tissue valves is applicable to [redacted].

5.2 [redacted] conducted a 2-year aging study for the [redacted] under protocol [redacted] "Two Year Real Time Shelf Life for [redacted] and [redacted]; and corresponding report [redacted] "Two Year Real Time Shelf Life Report of [redacted] Model [redacted]". This report presented the data for tissue properties and functional performance of the valves real-time aged to 2 years. The results showed that the valves met all the acceptance criteria as set forth in protocol [redacted]. The test was designed to compare the aged samples with the control samples (freshly made).

If there was negligible difference in tissue properties and if the valve functional performance met the hydrodynamic requirements, the products passed the acceptance criteria required to be labeled for the 2 year shelf life.

5.3 The results of aged [redacted] tested in protocol [redacted] were compared to results of freshly made [redacted] products at the 2-year time point passed all the acceptance criteria. These results indicate that the valves were not impacted by the aging process and would pass at any interim time point from time zero up to 2 years. It should be noted that the collection of data at interim time points (i.e., 3 months, 6, months, 1 year, 1-1/2 years etc.)

5.4 [redacted] Standard Operating Procedure, “Cardiovascular Shelf Life Procedure”, states that: “Real time aging involves the storing of test units at product storage conditions (typically in a monitored storage environment) for the desired shelf life period prior to testing”. The aged product is required to pass the performance specifications at the specified time frame. Successful demonstration of aged product testing indicates that aging has not impacted the product performance and the device can be labelled to the desired shelf life.

5.5 In addition to the 2 year real-time protocol [redacted] additional studies were conducted for the [redacted]. These protocols are: [redacted] - *Two Year Real Time Shelf Life*
[redacted] - *Three Year Real Time Shelf Life*
Qualification for Protocol (US Study) For [redacted] These studies have not completed the aging process but the interim data points will be collected and reported as the aging is completed.

5.6 A table showing a comparison of data collected at the various time points is presented in Table 2.0. As shown by the results in the table, all performance specifications were met for the aged non-[redacted] product. Additionally, the [redacted] s that were aged at 6 months and 1 year met the established acceptance criteria. This indicates that the aging process does not negatively impact the non-[redacted] device up to 2 year shelf life.

6.0 CONCLUSIONS

Based on the historical data performed on [redacted] s and the recent study of protocol [redacted] it can be concluded that [redacted] met the performance requirements at the stated time point up to 2 year time frame and additional interim time points are not deemed necessary.

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Data Sheet for Accelerated Aging

Insert Assembly

Requested Test Temperature 85°C
 Requested by: Benny Henderson Date: 10/2/08

Part numbers	45-0187-01	40-0183-01	45-0188-01	40-0184-01	44-0172-01
Description					
Quantity	8	8	8	8	8
Lot Number(s)	P00453	P00453	P00453	P00453	P00453

Action	Date	Time	Temperature (Degrees C)	Accumulated Test Time (years)	Functional Test Result Pass/Fail	Comments
Pre-test	10/2/08	1100-1500	21.5°C	0	PASS	B. Tucker 10/2/08
Start	10/2/08	1508	85°C	N/A	N/A	B. Tucker 10/2/08
Stop	10/8/08	0808	85°C	1 YEAR	N/A	B. Tucker 10/8/08
Test at 1 year	10/8/08	1100-1200	23.5°C	1 YEAR	PASS	T. Brinifits 10/8/08
Start	10/8/08	1306	85°C	N/A	N/A	B. Tucker 10/8/08
Stop	10/14/08	0808	85°C	2 years	N/A	B. Tucker 10/14/08
Test at 2 year	10/14/08	9:30 - 10:00 AM	24.2°C	2 years	PASS	T. Brinifits 10/14/08

The accumulated shelf life successfully demonstrated was 2 years.

Start date of real time aging study 10/13/08 T. Brinifits

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Report #:

CONFIDENTIAL -

- 5.2.4 Functional and Structural Testing –Sample size is based on ISO 5840 and FDA Heart Valve guidance for performing durability testing. Three (3) samples for mm and three (3) samples for mm were tested.
- 5.2.5 Control Samples - to be processed per the current manufacturing processes and sterilized once. These control samples were used for comparative data.
- 5.2.6 Histology testing was conducted for characterization purposes only. A single leaflet from nine separate valves were sufficient to provide information on intrinsic tissue characteristics and characterize any changes that may have occurred during the real-time aging.
- 5.3 Worst-Case
 - 5.3.1 The shelf life samples were sterilized two times prior to aging to represent worst case manufacturing conditions.

6.0 ACCEPTANCE CRITERIA

- 6.1 Visual Inspection

Each of the samples must be inspected visually in order to ensure the visual acceptance criteria of the valve and frame will be met. The support will be used as reference.
- 6.2 Storage solution tests

The pH and concentration data is for informational purposes only.
- 6.3 Shrink Temperature Test

Each of the samples tested shall have a shrink temperature greater than 80 degrees Celsius. A lower tolerance bound will be calculated based on the test results. The lower bound must be greater than 80 degrees Celsius.
- 6.4 Pronase

Perform t and F tests at 95% confidence level between the aged samples and the controls. Acceptance will require the test group to show either no significant difference or be significantly more resistant to protease digestion when compared to the control group.
- 6.5 Hydroxyproline Assay


Perform t and F tests at 95% confidence level between the aged test samples and the controls. Acceptance will require the test group to show either no significant

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
SIGNATURE OF STUDY PERSONNEL
Package Shelf Life Validation:
Straight Grasping Forceps

Performed by:

Tim / Study Executor
Ltd. 2011/07/19
Date Completed

Reviewed by:

Howard / Study Director
Ltd. 2011/07/19
Date Completed

Approved by:

Yuanmin / Manager 2011/07/19
Date Completed

Microbial Ranking Test - Subcontract Reviewed by:

Laboratory Number

Procedure Number

/ Study Director

LAB

TABLE OF CONTENTS

Statement of GLP Compliance 6

Quality Assurance Statement..... 7

Archiving 8

Introduction 9

Test method and procedure..... 10

Test Results..... 13

References 15

Sample photo..... 16

Attachment 1. Seal Peel Strength Test (Before / After Aging)..... 17

Attachment 2. Burst / Creep Test (Before / After Aging)..... 27

Attachment 3. Dye Penetration Test (Before / After Aging)..... 37

Attachment 4. Microbial Ranking Test (After Aging)..... 47

TESTING METHOD AND PROCEDURE

Test Article

Straight Grasping Forceps (LOT:677100909)

Aging procedure

5 years Accelerated Aging Incubation:

The shelf life study was conducted on the original packaging of sterile finished products, Forceps". After sterilization, the whole sterile packages were placed into the aging Incubator, and incubated for 57 days under controlled conditions (75°C) simulating the real time aging of 5 years. After accelerated aging study, the integrity of product packaging will be demonstrated by investigating the seal peel strength, burst & creep, dye penetration, the microbiological ranking properties of packages.

Package test:

I. Seal Peel Test:

Prepare the aging testing packages (containing intrinsic product inside) which was made of plastic film. Then the packages were conditioned at $22\pm 2^{\circ}\text{C}$ and humidity of $55\pm 5\%$ prior to the test at least 24 hours.

- 1.1 According to ASTM F88/F 88M-09, we use the universal test machine which connects to a computer (INSTRON 5581) to perform the seal peel strength on 10 sets of each testing packages.
- 1.2 Technique of holding the tail: Technique C. (The width of specimen:25mm)
- 1.3 Clamp the plastic sides of the test strip in the upper grip of universal test machine, and there was no support for package tail while testing.
- 1.4 Set the separation rate of the machine at 200 mm/min, then measure and record the maximum, minimum and average force required to peel each seal apart.

2. Burst and Creep Test:

Prepare the aging testing packages (containing intrinsic product inside) which is made of plastic film. Then the packages are conditioned at $23\pm 2^{\circ}\text{C}$ and humidity of $55\pm 5\%$ prior to the test at least 24 hours.

2.1 The burst test was performed on 10 packages based on ASTM F1140-07 using the LSSD-01, Labthink.

The flow rate and sensitivity were both set as 2 and 1.

2.2 The packages were inflated until rupture, it was the point that burst pressure was captured and recorded.

Burst volumes are expressed in inches of water.

2.3 Creep test was performed on a second batch of 10 samples at 80% of the above average burst pressure obtained from the burst testing (Creep Pressure).

2.4 The flow rate was set as the same as burst test (set as "2"). The packages were inflated until the creep pressure achieved and held for 30 seconds (Creep duration).

3. Dye penetration Test:

3.1 Conditions:

- 3.1.1 Prepare the aging testing packages (containing intrinsic product inside) which was made of plastic film.
- 3.1.2 Ensure that there is no any liquid water inside packaging. If there is any indication that the package has been exposed to any liquid water, it must be dried at applicable temperature before testing.
- 3.1.3 The test condition shall be controlled at $23 \pm 2^{\circ}\text{C}$ ($73.4 \pm 3.6^{\circ}\text{F}$) and $55 \pm 2\%$ relative humidity except the specific conditioning requirements.

3.2 Procedure

- 3.2.1 Inject sufficient dye penetrant solution into the package to cover the longest edge to a depth of approximately 5 mm by using hypodermic needle.
- 3.2.2 Rotate the package to ensure each seal edge has been exposed to the dye penetrant solution.
- 3.2.3 Each seal edge shall be exposed to dye penetrant solution for 20 seconds.
- 3.2.4 Observe the package whether any stain appeared on each seal edge or not.
- 3.2.5 The package which has no any stain appeared is considered to pass the test.

4. Microbial Ranking Test:

The samples were placed into filter units and attached to an exposure chamber. The chamber was subjected to a 15 minute aerosol of *Bacillus atrophaeus* spores at a target concentration of 1×10^6 CFU/mL ± 0.5 log. A flow rate of 2.8 L/min was drawn through each sample for the duration of the test. Spores which passed through the test samples were collected on membrane filters and enumerated. The logarithm reduction values (LRV) were calculated as the logarithm of the ratio of the challenge concentration to the permeate concentration.

TEST RESULTS

Package test:

1. Seal Peel Test:

Package site	Test Site	Average Strength (N/m)	Standard Deviation	Maximum (N/m)	Minimum (N/m)
Upper	Before aging		-----		
	After aging				
Down	Before aging		-----		
	After aging		-----		-----
Left	Before aging				-----
	After aging		-----		
Right	Before aging		-----		-----
	After aging		-----		

2. Burst and Creep Test:

Before aging:

	Average Strength (kPa)	Standard Deviation	Maximum (kPa)	Minimum (kPa)	80% of Burst pressure (kPa)
Burst pressure					
Creep pressure (kPa)	Pass		Fail		
	10		0		

Creep test: using Creep Pressure (80% of Burst pressure) at 3.1 kPa for 30 seconds.

After aging:

	Average Strength (kPa)	Standard Deviation	Maximum (kPa)	Minimum (kPa)	80% of Burst pressure (kPa)
Burst pressure					
Creep pressure (kPa)	Pass		Fail		
	10		0		

Creep test: using Creep Pressure (80% of Burst pressure) at 4.5 kPa for 30 seconds.

3. Dye Penetration:

	Result	Leaking
Before aging	Pass	No
After aging	Pass	No

4. Microbial Ranking Test:

After aging:

Test Article Number	Colony Forming Unit (CFU)	Log Reduction Value (LRV)
1	<1	>5.568
2	<1	>5.568
3	<1	>5.568
4	<1	>5.568
Reference Material	<1	>5.568

Negative control (CFU/membrane): <1

Challenge control (CFU): 3.7×10^5

REFERENCES

1. Good Laboratory Practice for Nonclinical Laboratory Studies. Part 58. Title of the US Code of Federal Regulations (1997) United States Food and Drug Administration.
2. ASTM F1980-07 Standard Guide for Accelerated Aging of Sterile Medical Device Packages.
3. ASTM F88 -85 Standard Test Method for Seal Strength of Flexible Barrier Materials. American Society of Testing Materials.
4. ASTM F1140:2000 Standard Test Methods for Internal Pressurization Failure Resistance of Unrestrained Packages for Medical Applications.
5. ASTM F1929-98 Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration.
6. ASTM D4332: 2001 Standard practice conditioning containers, packages or packaging components for testing.
7. ASTM 1608 Standard Test Method for Microbial Ranking of Porous Packaging Materials (Exposure Chamber Method).

Confirmation Letter For Shelf life validation

27 November 2012



To Whom it May Concern:

It is recommended that at least three replicate samples are tested at each sample interval, however due to its custom built manufacturing process, we selected at least two batches for shelf life

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Topic: needle Aging report	Document code	Approved by	Checked by	Checked by
Sub-topic:	Report date		1	
	Department and author	Quality Assurance Department		

1. Purpose

This trial was conducted to validate the shelf life of the adhesive specifications of the puncture needle, taking into account the characteristics of quality items and their stability which can be considered relevant to aging deterioration.

2. Conclusion

These results show that the various characteristics of the puncture needle() with the adhesive specifications showed no conspicuous aging deterioration, demonstrating that product quality is stable throughout the shelf life indicated.

The quality characteristics of the product after two sterilization cycles and a lapse of five years showed no conspicuous deterioration, and had equivalent capability compared with the product immediately after sterilization.

In addition, since performance complies with various standards, the safety of the product was also reconfirmed.

3. Products validated

Puncture needle

Table 1 List of products validated

	O.D.	I.D.	Grade	Lot No.	Sterilization cycles	Number of years after sterilization
①	0.72	0.49	ETW		1	Immediately after
②	↑	↑	↑		1	5 years
③	0.42	0.28	TW		1	Immediately after
④	↑	↑	↑		1	5 years
⑤	↑	↑	↑		2	5 years

Recital: For spinal needle and epidural needle with Hub Type B, quality of the various constituent materials is equivalent. However, some of the epidural needle outer cannulae(SUS304) are graduated by electrolytic marking (oxide layer).

4. Items examined: Quality evaluation items

4-1. Visual inspection

Complies with ISO 7864: 1993 4. Cleanliness
Check for discoloring, cracking and the like

4-2. Bonding strength inspection

Complies with ISO 7864: 1995.
13.1 Bonding between hub and needle tube
Spec. : 40N ≤ (22G 0.72mm)
: 22N ≤ (27G 0.42mm)
Equipment:

Fig.1 Overview of test machine

4-3. Limitation of acids and alkalis

Complies with ISO 7864: 1995 Section 5 Limits for acidity or alkalinity.
Products validated : ①, ②, ③, ⑤

4-4. Limitation of eluted substances

Complies with ISO 7864: 1995 Section 6 Limits for extractable metals.
Products validated : ①, ②, ③, ⑤

Recital: Sterility test
 Applicable standard :
 Assessment : No growth of bacteria or fungus
 Sample size : 1 (Bacteria and fungus respectively)
 Conditions : Details shown in Table 2.

Table 2 Sterility test conditions

	Culture medium			Culture conditions	
	Culture medium	Maker name	Lot No.	Temperature	Period
Bacteria	TGC medium	Nissui Pharmaceutical Co., Ltd.		32±2°C	14 days
Fungi	Tryptic Soy Broth	Beckon, Dickinson and Company		22±2°C	14 days

Reference: Strains used in performance tests of the above culture medium
 Indicator bacterium : Bacillus atrophaeus 10⁶ spores/strip Spore strip type
 Manufacturer : NAMSA
 Lot No. : N28200
 Test certificate : Yes (ISO 11138-1 compliant)
 (Using similar strains of bacteria and fungus)

5. results

The validation results for each item are shown in Table 3.
 The results for bonding strength for each item are shown in Tables 4 and 5 (Recital: The sterility test results are shown in Table 6).
 The various characteristics of each item meet the standards with no problem.

Table 3 Validation results

	①	②	③	④	⑤
External check	Passed	Passed	Passed	Passed	Passed
Adhesive Strength	Passed	Passed	Passed	Passed	Passed
Limitation of acids and alkalis	Passed	Passed	-	Passed	Passed
Limitation of eluted substances	Passed	Passed	-	Passed	Passed
Recital: Sterility test	Passed	Passed	Passed	Passed	Passed

Table 4 Outer needle bonding strength test results

(N)					
	①	②	③	④	⑤
AV					
MAX					
MIN					
SPEC.					

Table 5 Stylet bonding strength test results

(N)					
	①	②	③	④	⑤
AV					
MAX					
MIN					
SPEC.					

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Sterility test report

* Sterility test report.

* Media performance test documentation.

Sterility Test Report

Sample Information	
Location	_____
Product	_____
Lot No	_____
Period	_____

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**Introducer Set
Design Qualification Report
Report #**

1. Purpose

1.1 The purpose of this report is to document the findings of the
Introducer Set Design Qualification Protocol #
1) for the time = 5 year accelerated aging samples.

2. Summary

- 2.1. A total of 70 packages were assessed according to transportation testing standards. Post transportation simulation samples 31-40 were aged for 190 days at 55C / 70-80% RH, which is equivalent to 5 years of storage. After the aging cycle, the samples were visually inspected for physical damage on each component (sterile barrier pouch and c packages were then tested for seal strength and five were tested for package integrity (bubble test).
- 2.2. Post Transportation Results (Time = 0)
- 2.2.1. The sterile barrier pouch passed both the package integrity testing (bubble test) and the seal strength test. As a result the package design is considered to be adequate for use in a typical shipping and distribution environment (Refer to r complete report).
- 2.3. 1 year Accelerated Aging Results
- 2.3.1. The sterile barrier pouch passed the seal strength test, but failed integrity testing (bubble test). One pouch failed the bubble leak test as a result of a puncture on the Polyethylene/polyester side of the pouch, near the stopcock.
- 2.3.2. This damage was not noted in the visual inspection post transportation simulation, and therefore must have been the result of mishandling post transportation simulation, the effects of heat from the accelerated aging process, or excessive pressure during the bubble leak test causing the weakened spot in the plastic to fail (Refer to for complete report).
- 2.3.3. Therefore, it was inconclusive at that time as to whether or not the packaging was adequate for use after 1 year of storage. To ensure that the transportation/aging did not contribute to the failures described above, package approval was pending the results of the 3 year accelerated aging report. In addition, data from this 5 year accelerated aging samples and 1, 3, and 5 year samples shall be analyzed to support these findings.
- 2.4. 3 year Accelerated Aging Results
- 2.4.1. The sterile barrier pouch passed both the package integrity testing (bubble test & visual inspection) and the seal strength test. As a result the package design is considered to be adequate for use in a typical shipping and distribution environment and after 3 years of storage (Refer to DV.23.R.2 for complete report).
- 2.5. 5 year Accelerated Aging Results
- 2.5.1. The five year accelerated aging results are summarized in section 4.

Company Confidential

Introducer Set
Design Qualification Report
Report # -----

2.6. Packaging Configuration:

- 2.6.1.** Several components were packaged into inner 1059B Tyvek pouches (CPN's 1194, 1206, & 2939), which are packaged into the outer 1073B Tyvek pouch (CPN 4438). Fifteen (15) pouches were then packaged into a 58" x 9" x 5" standard size to fill the headspace of the shipper to prevent uncontrolled internal movement of product.

3. Deviations

- 3.1.1.** None

4. Results

4.1. Sterile Barrier Pouch

- 4.1.1.** All five pouches passed the visual inspection and bubble leak test.
4.1.2. The average seal strength post transportation was 2.75 lbf, which is higher than the 0.5lbf acceptance criteria stated in the attached protocol.

4.2. Device

- 4.2.1.** No damage was observed that would interfere with the proper use of the product.

5. Conclusion

- 5.1.** The sterile barrier pouch passed both the package integrity testing (bubble test & visual inspection) and the seal strength test. As a result the package design is considered to be adequate for use in a typical shipping and distribution environment and after 5 years of storage.
5.2. This report shall be documented in the project binder for project # "MPC 010" "Master Packaging Validation."

6. Attachments

- Attachment 1 – Executed Protocol
Attachment 2 – Bubble Test Results
Attachment 3 – DDL Accelerated Aging Report

Company Confidential

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1, 3 & 5 Year Accelerated Aging Package Evaluation Testing
for the

Check Flow (

Report Number 0608043A

PO Number: 232604

14 May 2007

Customer:

Testing Performed by

Tested and proven.

DDL, Inc.
10200 Valley View Road, Suite 101
Eden Prairie, MN 55344
952-941-9226
Fax 952-941-9318
1-800-229-4235

DDL West, Inc.
3303 Harbor Boulevard, Suite B9
Costa Mesa, CA 92626
714-979-1712
Fax 714-979-1721
1-800-229-4235

www.testedandproven.com

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Package Integrity Testing

Bubble Leak Test

This test was performed based on ASTM F 2096-04; *Detecting Gross Leaks in Medical Packaging by Internal Pressurization (Bubble Test)* and DDL, Inc. S.O.P. #6047, Rel. 03.

The test method is designed to detect open pathways, channels, or voids across the seal area intended as a primary sterile barrier, and to detect pinholes in non-porous materials. The test may be used for flexible or rigid packages with porous or non-porous materials.

A set-up specimen was used to determine the set pressure by creating a 0.008 inch hole in the Tyvek side of the package. The test apparatus described in the 'materials and equipment' section was used for air pressure. The package was submerged in a water bath for at least five (5) seconds with no air pressure being applied to the pouch. The pressure was then increased slowly until a steady flow of bubbles came from the 0.008 inch hole created in the Tyvek. The packages were observed for any evidence of leakage in the form of a steady stream of bubbles emanating from the package.

All raw data can be found in Appendix D of this test report. A tabular summary of the results is presented in the results section of this test report.


Figure 7: This is a photograph of the bubble leak test set-up.

**Aerosol Spore Challenge (Exposure Chamber Method)
ASTM F1608 GLP Report**

Test Article: 1)
2)
3)

Purchase Order: 379783
Laboratory Number: 646180
Study Received Date: 26 Jul 2012
Test Procedure(s): Standard Test Protocol (STP) Number: STP0112 Rev 05
Protocol Detail Sheet (PDS) Number: 201202957 Rev 01

Summary: This test employs a low flow rate, extended exposure time, and increased challenge level, which allows testing of materials with high differential pressure values. This test method is not intended to predict the performance of a given material in a specific sterile packaging application. All test method acceptance criteria were met.

Results: Any results with an approximate sign (~) or a greater than sign (>) are considered an approximation based on counts outside of the statistically accurate range [25-250 colony forming units (CFU) per plate].

Test Article	CFU	LRV
1	1.3 x 10 ³	2.8
2	4.5 x 10 ²	3.3
3	3.8 x 10 ²	3.4
4	~2.1 x 10 ²	~3.6

LRV = Log Reduction Value
Negative Control (CFU/membrane): <1
Challenge Control (CFU): 8.5 x 10⁵

Test Article	CFU	LRV
1	~1.6 x 10 ²	~3.6
2	~9.3 x 10 ¹	~3.8
3	~1.0 x 10 ¹	~4.8
4	~4.0 x 10 ¹	~4.2

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Monopolar is the most commonly used electrosurgical modality. This is due to its versatility and clinical effectiveness. In monopolar electro-surgery, the active electrode is in the surgical site. The patient return electrode is somewhere else on the patient's body. The current passes through the patient as it completes the circuit from the active electrode to the patient return electrode.

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Monopolar Circuit

This picture represents a common monopolar circuit. There are four components to the monopolar circuit: Generator, Active Electrode, Patient, and Patient Return Electrode.

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
[redacted] Coagulators ([redacted])

The proprietary design of [redacted] the shaft temperature as compared to existing [redacted] coagulators and may reduce the potential for oral burns.

Inadvertent thermal injuries during routine otolaryngology procedures can result from inadequate device insulation.¹

The shaft temperature of [redacted] is lower than that of the competitive devices evaluated (see chart below).^{2,3}

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1. Richter [redacted], 2008;72(7):1013-1021.
2. Internal testing of [redacted] coagulator, [redacted] devices. Data on file at [redacted]
3. [redacted] and the [redacted] are used in procedures that require a smaller diameter and do not have the additional insulation.

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	<input style="width: 100%;" type="text"/>	Measurements & Monitors
ANSI/AAMI SP10 Clinical Evaluation Report		

Performed by: The performer of the test report is listed below:

Function, Title	Name	Signature	Date
Clinical Liaison	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>

Document Approvers: The document approvers are listed below:

Function, Title	Name	Signature	Date
NBP Program Manager	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>	<input style="width: 80%;" type="text"/>

ANSI/AAMI SP10 Clinical Evaluation Report

1 Summary

The [] in combination with Mobile [] cuffs ([]), [], Mobile CL Reusable cuffs ([]), [] and the cradle [] was evaluated according ANSI/AAMI SP10:2002/(R)2008+A1:2003/(R)2008 using the auscultatory reference method.

485 valid paired measurements were obtained from 125 subjects including 19 pediatric subjects and compared to indirect cuff/stethoscope auscultation. Two trained observers determined systolic and diastolic blood pressure using the Korotkoff method. The 5th Korotkoff sound was used to determine diastolic pressure. The sequential single-arm method was used because the stepwise cuff deflation of the instrument didn't allow simultaneous instrument / observer measurements. Differences between the two observers were verified to meet the criteria of ANSI/AAMI SP10:2002/(R)2008+A1:2003/(R)2008. The measurement results of the device under test had to yield a mean difference versus the averaged observer measurements of no greater than ± 5 mmHg, with a maximum standard deviation of 8 mmHg.

Summary data for data evaluation according to ANSI/AAMI SP10:2002/(R)2008 +A1:2003/(R)2008:

The limb circumferences ranged from [] cm to 41 cm, with a distribution of 20.8% below 25 cm and 2.8% above 35 cm. The auscultatory systolic pressure ranged from 71.5 mmHg to 153.5 mmHg, with an average of 114,6 mmHg and with a distribution of 19,2% below 100 mmHg and 0,0% above 160 mmHg. The auscultatory diastolic pressure ranged from 40.5 mmHg to 102.0 mmHg, with an average of 72.0 mmHg and with a distribution of 17,6% below 60 mmHg and 0,8% above 100 mmHg. 16,0% of the subjects were measured supine, 54,4% were measure seated and 29,6% were measured standing. The mean errors and standard deviations in the table below demonstrate that the mean error and the standard deviation of all device measurements compared to the respective observer reference measurements are within the limits of the ANSI/AAMI SP10:2002/(R)2008 +A1:2003/(R)2008:

	Method 1 (individual measurements)		Method 2 (subject averages)	
	Systolic [mmHg]	Diastolic [mmHg]	Systolic [mmHg]	Diastolic [mmHg]
Mean Difference	4.10 (± 5)	-1.63 (± 5)	4.11 (± 5)	-1.65 (± 5)
Standard Deviation	6.81 (<8)	6.87 (<8)	4.49 (<5.24)	5.03 (<6.74)
# of measurements	485	485	125	125

The values in () are the limits according to ANSI/AAMI SP10:2002/(R)2008+A1:2003/(R)2008. For method 2 the maximum standard deviation depends on the actual mean difference.

The obtained clinical data together with the laboratory validation data demonstrate that the [] in combination with Mobile [] cuffs [], [] cuffs ([]), [] and the cradle [] meets the requirements (mean error of less than ± 5

ANSI/AAMI SP10 Clinical Evaluation Report

2 Data Collection and Clinical Evaluation Plan

2.1 Non-Clinical Data Collection

Site	
Department	

2.2 Locations for Clinical Data Collection

Hospital	
Department	

2.3 Objectives

The objective of this clinical evaluation was to determine that the overall accuracy of the [redacted] [redacted] meets the requirements as expressed in the ANSI/AAMI SP10:2002/(R)2008 +A1:2003/(R)2008 for adult and pediatric subjects. The auscultatory method is used as the reference method.

For systolic and diastolic pressures, treated separately, the mean difference of paired measurements of the device under test and the auscultatory method shall be ± 5 mmHg or less, with a standard deviation of 8 mmHg or less.

2.4 Legal and Ethical Considerations

This study conforms to the terms of the Declaration of Helsinki (latest revision), [redacted] [redacted] Device Directive, FDA, applicable IIR requirements, and each site's country, state, and local IRB [redacted] The patient safety of equipment [redacted] [redacted]

Laboratory performance evaluation of the device under test has been conducted using an advanced pressure simulator as a source. The overall accuracy remains to be evaluated and compared under clinical conditions.

There were no additional risks for the volunteer or patient because:

- No invasive procedures, or drawing blood or tissue samples were carried out due to this study
- No additional, or different medical treatment to the subjects, of any kind, was required in conjunction with this study

ANSI/AAMI SP10 Clinical Evaluation Report

3 Results

Subjects

Overall 125 subjects including 19 pediatric subjects entered the clinical study. All collected data were evaluated.

Number of subjects:	125
Number of measurements:	485

Table 1 - Number of subjects and paired measurements

Range and distribution of limb circumferences

Figure 1 summarizes the distribution of limb circumferences as seen in this clinical evaluation. The limb circumferences range from cm to 41 cm, with a distribution of 20.8% below 25 cm and 2.8% above 35 cm.


Figure 1 - Distribution of limb circumference

Range and distribution of systolic blood pressure

Figures 2 summarizes the distribution of systolic blood pressures as seen in this clinical evaluation. The auscultatory systolic pressure ranges from 71.5 mmHg to 153.5 mmHg, with an average of 114,6 mmHg and with a distribution of 19,2% below 100 mmHg and 0.0% above 160 mmHg.

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ANSI/AAMI SP10 Clinical Evaluation Report

Bland-Altman graphical plots of systolic differences against the systolic reference values


Figure 4 - graphical display of systolic differences against systolic values determined with the Korotkoff method

Bland-Altman graphical plots of diastolic differences against the diastolic reference values


Figure 5 - graphical display of diastolic differences against diastolic values determined with the Korotkoff method

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ANSI/AAMI SP10 Clinical Evaluation Report

Percentages of readings with differences within 5, 10 and 15 mmHg

	Systolic	Diastolic
NBP readings within ± 5 mmHg of reference	54.6%	64.1%
NBP readings within ± 10 mmHg of reference	81.2%	87.0%
NBP readings within ± 15 mmHg of reference	94.0%	95.3%

Table 2 - Distribution of paired measurement differences in percentages

Summary discussion:

The obtained clinical data together with the laboratory validation data demonstrate that the [redacted] [redacted] in combination with [redacted] cuffs [redacted], [redacted], [redacted] cuffs [redacted], [redacted] and the cradle [redacted] meets the specified accuracy (mean error of less than ± 5 mmHg and a standard deviation of less than 8 mmHg) over the specified blood pressure and pulse rate range.

4 Adverse Events

As expected no adverse events caused by the study itself have been reported.

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Artificial Organs
20(12):1278-1281, Blackwell Science, Inc., Boston
© 1996 International Society for Artificial Organs

Clinical Use of Omniflow Vascular Graft as Arteriovenous Bridging Graft for Hemodialysis

Shoei-Shen Wang and Shu-Hsun Chu

Department of Surgery, National Taiwan University, Taipei, Taiwan

Abstract: From April 1992 to November 1995, 109 bridging arteriovenous fistulas were performed in 89 patients for hemodialysis at National Taiwan University Hospital. These included 61 Omniflow vascular grafts (OVG) in 47 patients and 48 GoreTex polytetrafluoroethylene (PTFE) grafts in 42 patients. There were 36 male and 53 female patients, whose ages ranged from 24 to 84 years with a mean of 59 ± 2 years. The operation was performed under lidocaine local anesthesia with antibiotic prophylaxis. The anatomic configurations included 50 (45.9%) forearm loop grafts, 19 (17.4%) forearm straight grafts, and 40 (36.7%) thigh loop grafts. Hemodialysis via the bridging graft was started 2 weeks after the operation. The patients received hemodialysis 3 times a week or more when indicated. No patients suffered from limb loss. Pseudoan-

eurysm and graft infection requiring operation occurred 0.96% and 1.92% per patient year, respectively, in the OVG group, and pseudoaneurysm and graft infection both occurred 2.38% per patient in the PTFE group. In the OVG group, the thigh loop graft had better patency than the forearm straight graft. The forearm loop graft had a better patency rate than the forearm straight graft. In the OVG group, the 1 to 4 year patency rate was $77 \pm 6\%$, $58 \pm 8\%$, $48 \pm 9\%$, and $34 \pm 13\%$, respectively, whereas in the PTFE group, the 1 to 4 year patency rate was $69 \pm 7\%$, $50 \pm 9\%$, $43 \pm 10\%$, and $26 \pm 15\%$, respectively. Therefore, the Omniflow is a good alternative bridging graft for hemodialysis. **Key Words:** Blood vessel prosthesis—Vascular patency—Omniflow—Arteriovenous fistula—GoreTex—Polytetrafluoroethylene.

In chronic renal failure patients, the autogenous Brescia-Cimino fistula is usually the most popular and satisfactory vascular access for hemodialysis (1,2). When autogenous arteriovenous fistula (AVF) is no longer feasible, the use of a prosthetic conduit to form a bridging AVF is the best alternative. The expanded polytetrafluoroethylene graft was the most commonly used vascular prosthesis (3-5), yet it had a high rate of thrombosis, infection, or puncture complications (6,7). Because the patency rate of biosynthetic grafts was reported to be higher than the patency rate of expanded polytetrafluoroethylene graft (8,9), we introduced the biosynthetic graft as a bridging graft for hemodialysis.

The Omniflow II vascular graft (OVG) is a biosynthetic vascular prosthesis composed of glutaraldehyde-tanned ovine collagen which is grown around a polyester mesh template. The polyester

mesh is an integral part of the graft wall and is encapsulated by Type 3 collagen whereas the outer layers are formed by longitudinally oriented Type 1 collagen (10-12). The collagen may provide the biocompatibility and hemocompatibility whereas the polyester mesh may enhance the durability. Therefore, we started to use the OVG as an arteriovenous bridging graft for hemodialysis in 1992. In this article, we report our 3 year experiences with this biosynthetic composite vascular graft.

MATERIALS AND METHODS

From April 1992 to November 1995, at National Taiwan University Hospital, 591 operations in 482 patients were performed to create hemodialysis vascular access, including 109 (18.4%) bridging AVF and 482 (81.6%) autogenous AVF. This included all patients with new-onset renal failure and many other patients with complicated access histories who were referred for final care after many operations had been tried elsewhere and had failed. When our operations failed, many patients shifted from hemodialysis to continuous ambulatory peri-

Received August 1996.

Address correspondence and reprint requests to Dr. Shoei-Shen Wang, National Taiwan University Hospital, No. 7 Chung-Shan South Road, Taipei, Taiwan.

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■ Articles

➤ Percutaneous closure of the left atrial appendage versus warfarin therapy for prevention of stroke in patients with atrial fibrillation: a randomised non-inferiority trial

David R Holmes, Vivek Y Reddy, Zoltan G Turi, Shephal K Doshi, Horst Sievert, Maurice Buchbinder, Christopher M Mullin, Peter Sick, for the PROTECT AF Investigators*

Lancet 2009; 374: 534-42

See Editorial page 501

See Comment page 504

*Members listed at end of paper

Mayo Clinic College of Medicine, Rochester, MN, USA (Prof D R Holmes MD); Mount Sinai School of Medicine, New York, NY, USA (V Y Reddy MD); Cooper Hospital, Camden, NJ, USA (Z G Turi MD); Pacific Heart Institute/St Johns Hospital, Santa Monica, CA, USA (S K Doshi MD); Cardiovascular Center Frankfurt, Sankt Katharinen, Frankfurt, Germany (Prof H Sievert MD); Foundation for Cardiovascular Medicine, La Jolla, CA, USA (M Buchbinder MD); Integra Group, Brooklyn Park, MN, USA (C M Mullin MS); and Krankenhaus der Barmherzigen Brüder, Regensburg, Germany (P Sick MD)

Summary

Background In patients with non-valvular atrial fibrillation, embolic stroke is thought to be associated with left atrial appendage (LAA) thrombi. We assessed the efficacy and safety of percutaneous closure of the LAA for prevention of stroke compared with warfarin treatment in patients with atrial fibrillation.

Methods Adult patients with non-valvular atrial fibrillation were eligible for inclusion in this multicentre, randomised non-inferiority trial if they had at least one of the following: previous stroke or transient ischaemic attack, congestive heart failure, diabetes, hypertension, or were 75 years or older. 707 eligible patients were randomly assigned in a 2:1 ratio by computer-generated randomisation sequence to percutaneous closure of the LAA and subsequent discontinuation of warfarin (intervention; n=463) or to warfarin treatment with a target international normalised ratio between 2.0 and 3.0 (control; n=244). Efficacy was assessed by a primary composite endpoint of stroke, cardiovascular death, and systemic embolism. We selected a one-sided probability criterion of non-inferiority for the intervention of at least 97.5%, by use of a two-fold non-inferiority margin. Serious adverse events that constituted the primary endpoint for safety included major bleeding, pericardial effusion, and device embolisation. Analysis was by intention to treat. This study is registered with Clinicaltrials.gov, number NCT00129545.

Findings At 1065 patient-years of follow-up, the primary efficacy event rate was 3.0 per 100 patient-years (95% credible interval [CrI] 1.9-4.5) in the intervention group and 4.9 per 100 patient-years (2.8-7.1) in the control group (rate ratio [RR] 0.62, 95% CrI 0.35-1.25). The probability of non-inferiority of the intervention was more than 99.9%. Primary safety events were more frequent in the intervention group than in the control group (7.4 per 100 patient-years, 95% CrI 5.5-9.7, vs 4.4 per 100 patient-years, 95% CrI 2.5-6.7; RR 1.69, 1.01-3.19).

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Intracorneal inlay to correct presbyopia: Long-term results

Ömer F. Yılmaz, MD, Neşe Alagöz, MD, Gökhan Pekel, MD, Engin Azman, MD,
Ebru F. Aksoy, MD, Hanefi Çakır, MD, Ercüment Bozkurt, MD, Ahmet Demirok, MD

PURPOSE: To evaluate the long-term visual results of Acufocus ACI-7000 (now Kamra) intracorneal inlay implantation in presbyopic phakic patients.

SETTING: Beyoğlu Eye Training and Research Hospital, Istanbul, Turkey.

DESIGN: Clinical trial.

METHODS: This study comprised patients with emmetropic or post-laser in situ keratomileusis (LASIK) presbyopia. Patients had an uncorrected near visual acuity (UNVA) of 20/40 or worse, correctable to 20/25 or better at distance. The inlay was implanted on the stromal bed after the LASIK flap was relifted or a flap created. The inlay was centered on the presurgical position of the first Purkinje reflex. The main outcome measures were distance and near vision and the complication rate.

RESULTS: The study enrolled 39 patients aged 45 to 60 years. At the 4-year follow-up, all patients (N = 22) had 2 or more lines of improvement in UNVA with no significant loss in distance vision. The mean final UNVA was 20/20 (Jaeger [J1]); 96% of patients could read J3 or better. The uncorrected distance acuity was 20/40 or better in all eyes. Five patients had cataract progression, and 2 had a change in refractive status. No eye with an intracorneal inlay had intraoperative complications during cataract extraction. Four inlays were explanted during the study. There were no severe corneal complications that affected final vision.

CONCLUSION: Intracorneal inlay implantation was an effective, safe, and reversible procedure for the long-term surgical treatment of presbyopia.

Financial Disclosure: No author has a financial or proprietary interest in any material or method mentioned. Additional disclosure is found in the footnotes.

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DEPARTMENT OF HEALTH & HUMAN SERVICES Public Health Service

Food and Drug Administration
10903 New Hampshire Ave.
Silver Spring, MD 20993

Certificate No. 4317-8-2010


CERTIFICATE TO FOREIGN GOVERNMENT

In order to allow the importation of United States products into foreign countries, the U.S. Food and Drug Administration (FDA) certifies the following information concerning the product(s) to be exported listed below:

<u>Name of Product(s)</u>	<u>Name of Manufacturer/Distributor Address</u>
See Attached List (1 page)	<div style="border: 1px solid black; width: 200px; height: 30px;"></div>

The product(s) described above (and the manufacturing/distribution site(s) which produces/distributes it) is subject to the jurisdiction of the FDA under the Federal Food, Drug, and Cosmetic Act.

It is certified that the above product(s) may be marketed in, and legally exported from, the United States of America at this time. The manufacturing plant(s) in which the product(s) is produced is subject to periodic inspections. The last such inspection showed that the plant(s), at that time, appeared to be in substantial compliance with current good manufacturing practice requirements for the products(s) listed above.


Chester T. Reynolds
Acting Deputy Director,
Division of Risk Management Operations
Office of Compliance
Center for Devices and Radiological Health

This certificate expires 24 months from the date notarized.

COUNTY OF MONTGOMERY
STATE OF MARYLAND

Subscribed and sworn to before me this 21 day of sept month 2010 year. ·

Regulatory Status

The has been approved for commercial distribution in the following countries or regulatory jurisdictions:

Country or Regulatory Jurisdiction	Commercial Approval
European Union (CE Mark)	November 2005
Bahamas	<input type="checkbox"/> 2006
Israel	October 2006
South Africa	March 2007
Saudi Arabia	October 2007
Bulgaria	February 2008 ¹
United States	April 2008 <input type="checkbox"/> January 2010 <input type="checkbox"/>
Canada	May 2009 <input type="checkbox"/> October 2010 <input type="checkbox"/>
Thailand	July 2009
Taiwan	May 2010
Australia	<input type="checkbox"/>
Turkey	March 2011
Russia	March 2011

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순번	인정 번호	기관명	공통규격 (전기안전)	보조규격	전자파규격	개별규격
8	KT123	한국에스지에스 (주) 군포시험소	CAN/CSA-C22.2No.6 0601.1:2008 UL60601-1:2003(R2006) EN60601-1:1990+A1:1993+A2:1995 EN60601-1:2006 IEC60601-1:1998+A1:1991+A2:1995 IEC60601-1:2005(ThirdEdition) IEC60601-1:2005(ThirdEdition) KSPIEC60601-1:2002 KSCIEC60601-1:2008	EN60601-1-1:2001 EN60601-1-3:1994 EN60601-1-8:2004+A1:2006 IEC60601-1-1:2000 IEC60601-1-8:2003A1:2006 (Firstedition) KSCIEC60601-1-1:2002 KSCIEC60601-1-3:2002 KSCIEC60601-1-8:2005	EN60601-1-2:2007 IEC60601-1-2:2007	EN60601-2-10:2000+Am1:2001 EN60601-2-10:2000+Am1:2001 EN60601-2-2:2001 EN60601-2-2:2007 EN60601-2-22:1996 EN60601-2-27:2006 EN60601-2-28:1998 EN60601-2-32:1994 EN60601-2-35:1996 EN60601-2-38:1996+A1:2000 EN60601-2-42:2003 EN60601-2-44:2001+A1:2003 EN60601-2-49:2001 EN60601-2-7:1998 EN60601-2-9:1996 IEC60601-2-10:1987(FirstEd)+Am1:2001 IEC60601-2-2:1998(Thirdedition) IEC60601-2-2:2003(Fourthedition) IEC60601-2-22:1996(Secondedition) IEC60601-2-2-28:1998(Firstedition) IEC60601-2-27:2003(Secondedition) IEC60601-2-32:1994(Firstedition) IEC60601-2-35:1996(Firstedition) IEC60601-2-38:1996(Firstedition)+A1:1999 IEC60601-2-42:2002(Secondedition) IEC60601-2-44:2001+A1:2002(Ediior2) IEC60601-2-49:2001(Firstedition) IEC60601-2-7:1998(Secondedition) IEC60601-2-9:1996(Secondedition) KSCIEC60601-2-10:2002 KSCIEC60601-2-22:2002 KSCIEC60601-2-27:2005 KSCIEC60601-2-28:2002 KSCIEC60601-2-32:2002 KSCIEC60601-2-35:2002 KSCIEC60601-2-38:2002 KSCIEC60601-2-42:2003 KSCIEC60601-2-44:2005 KSCIEC60601-2-49:2005 KSCIEC60601-2-7:2002 KSCIEC60601-2-9:2002 KSPIEC60601-2-2:2002

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